

NASA-CR-189,543

NASA Contract Report 189543

NASA-CR-189543
19920011041

MEASUREMENT OF VORTEX FLOW FIELDS

T. Kevin McDevitt, Todd A. Ambur, Gary M. Orngard, F. Kevin Owen


COMPLERE INC.
Palo Alto, CA

Contract NAS1-18667
January 1992



National Aeronautics and
Space Administration

Langley Research Center
Hampton, Virginia 23665-5225



FOR REFERENCE

NOT TO BE TAKEN FROM THIS ROOM

Abstract

The objective of Phase II was to design, build and demonstrate a three dimensional laser fluorescence anemometer for use in the Langley 16- by 24- Inch Water Tunnel. Innovative optical design flexibility combined with compact and portable data acquisition and control systems have been incorporated into the instrument. This will allow its use by NASA in other test facilities. The final instrument and support systems differ in several significant aspects from the design envisaged during our Phase II proposal preparations. Our original mirror traverse alignment concept has been replaced by a more versatile fiber optic system. This facilitates normal and off-axis beam alignment, removes mirror losses and improves laser safety. This added optical flexibility will also enable simple adaptation for use in the adjacent jet facility. New proprietary concepts in transmitting color separation, light collection and novel prism separation of the scattered light have also been designed and built into the system. Off-axis beam traverse and alignment proved much more complex than initially conceived. This led to the requirement for a specialized, programmable traverse controller and the inclusion of an additional traverse for the off-axis arm. To meet this challenge, an "in-house" prototype unit was designed and built and traverse control software developed specifically for the water tunnel traverse applications. A specialized data acquisition interface was also required. This was designed and built for the Laser Fluorescence Anemometer system.

Introduction

At present, significant efforts are being made to effect design changes which will improve aircraft agility, maneuverability and performance. But, although significant progress has been made in computational aerodynamics, reliable design changes still cannot be made without recourse to experiment. Attempts to extend tactical flight envelopes still require extensive preflight ground based model testing. Unfortunately, conventional wind tunnel testing is expensive and time consuming and most facilities were built before present-day optical methods for quantitative flow field measurements were envisaged. Consequently, the non-intrusive detailed documentations of lee-side vortex flow-fields which are often required to support design evaluation and optimization are few.

However, in the past, qualitative water tunnel simulations have guided many practical designs and, since most of these facilities have been built with excellent optical access, they are ideally suited for use in advanced flow field diagnostics. Since the performance of most lifting and maneuvering bodies is governed by extensive transitional and turbulent viscous wakes and vortical lee-side flows, non-intrusive optical measurement techniques are required. Consequently, water tunnels offer the opportunity to obtain inexpensive, detailed flow field measurements to support "cut and try" designs and basic research.

To realize this capability, a two dimensional laser fluorescence anemometer was built and tested in the Ames-Dryden Water Tunnel during Phase I. The instrument was used in an experimental study of vortex flow fields designed to determine the mechanisms and feasibility of controlling vortex breakdown by introducing relatively low rates of jet blowing along the vortex core. The objective of Phase II was to build a three dimensional instrument for studies in the

Background

When a slender delta wing is at an angle of attack to an oncoming stream, the upper and lower surface boundary layers flow outward and separate from the leading edges to form two free shear layers that roll up into a pair of vortices above the wing. Increasing angle of attack strengthens the vortices until the induced wing pressure field and associated adverse streamwise pressure gradients cause vortex breakdown. The flow is further complicated as the leading edge vortices mix with the wake from the trailing edge downstream of the wing. The phenomenon of vortex breakdown (or vortex bursting) can have a significant influence on control surface performance and unsteady loading. The inherent unsteadiness of the breakdown process compounds the problem as it continually moves the breakdown region back and forth along the vortex axis. This creates serious time dependent flow problems and asymmetrically disposed breakdown positions above the wing that are aggravated with side-slip.

Wide variations of breakdown patterns have been observed, and with increasing swirl the patterns change from spiral to near axisymmetric (Ref. 1). Spiral breakdown most commonly occurs over delta wings. In this breakdown process, the filament of fluid along the axis does not spread out symmetrically from a fixed stagnation point but, instead, takes on a spiral form around an unsteady "stagnation point" which varies in both space and time. Axisymmetric breakdown over delta wings, although rare, can also occur (Ref. 2). In this case, the vortex has a roughly axisymmetric breakdown pattern with a characteristic bubble which can have single or multiple tails (Ref. 3).

Unfortunately, the parameters and conditions that result in vortex breakdown are poorly understood, because reliable quantitative experimental data are difficult to obtain. With limited experimental information to guide flow field modeling, numerical studies of vortex breakdown and control have met with only limited success (Ref. 4). There have been two principal reasons for this. In the first place, flow field unsteadiness associated with breakdown produces directional intermittency. This leads to large uncertainties in mean and unsteady flow measurements obtained with conventional pitot and hot wire probes. Secondly, and perhaps more important, is the fact that vortex breakdown is known to be extremely sensitive to any form of introduced disturbance. Probes, due to their blockage, may drastically alter the breakdown position. For these reasons, almost complete reliance has been placed on flow visualization techniques to determine flow field characteristics. In the past, air or hydrogen bubbles have been used as tracers to visualize flow patterns in water tunnels. For steady flows, streak lines can be identified with streamline patterns. However, in more complex flows of practical interest, the use of bubbles for flow visualization has distinct drawbacks. First of all, their introduction acts as a fluid lubricant which alters the apparent fluid viscosity and so its turbulent structure. Secondly, light refraction at the gas/water interfaces will destroy laser beam coherence and make it impossible to obtain laser velocimeter measurements in the regions where the tracer is present. But, with the advent of the laser fluorescence anemometer, there are now opportunities to determine accurate quantitative flow field velocity measurements of the vortex bursting process.

The Laser Fluorescence Anemometer

The principle of operation of the laser fluorescence anemometer is shown schematically in Fig. 1. The mean velocity and turbulence measurements are made with a dual-beam velocimeter utilizing a Bragg cell that enables moving interference fringes to be generated in the focal volume so that instantaneous velocity magnitude and direction measurements can be achieved from the frequency shift (f_D) around the incident and modulated laser beam interference frequency (f_0). i.e. $U = \lambda(f_D - f_0) / 2 \sin(\theta/2)$ where λ is the wavelength of the incident laser light. Mean and fluctuating concentration measurements are achieved by observing the intensity of fluorescent light emitted from the focal volume at a different wavelength (λ_C). At correct levels, tracer fluorescence is linearly proportional to the trace material concentration and, therefore, fluorescent intensity is directly proportional to the concentration of fluid from the seeded flow in the focal volume. The cathode current from a second photomultiplier tube is coupled to a high-gain current to voltage converter to produce a continuous voltage proportional to the instantaneous concentration.

Since fluorescence is such a complex phenomenon dependent upon many parameters, only those of particular importance to the present application have been considered. In this context, the principal requirements were linearity combined with adequate sensitivity (i.e., signal/noise ratio), frequency response, and spatial resolution. Since the fluorescent intensity of a particular organic dye can be a strong function of the solvent, which in this case was room temperature water, other dye-solvent-temperature combinations would produce different (and possibly increased) fluorescence. However, the fluorescent output was more than adequate for the present tests. In addition, the relationship between fluorescent intensity and concentration is of course exponential but, at the extremely low dye concentrations used in these experiments, a linear approximation could be made without introducing significant errors.

The data in Fig. 2 show this linearity of the present technique, which employed dysodium fluorescein dye (a sodium salt of fluorescein which has been used for flow visualization for many years) as the trace material in water. Since, for a given dye concentration, the measured fluorescence is a function of laser beam intensity and collection optics, the ordinate of Fig. 2 has been plotted in arbitrary (voltage) units. The high signal-to-noise ratio that can be obtained for dye concentrations as low as 0.04 ppm (by weight) is illustrated in Fig. 3. The rise time ($\approx 50 \mu\text{sec}$), which corresponds to the time taken to chop the laser beam, shows the adequate frequency response of the system.

Since there is usually an overlap of the absorption and emission curves for most dye-solvent combinations it is possible that fluorescent photons emitted from the probe volume could be reabsorbed by the dye molecules that are between the probe volume and the collection optics. This effect could produce a range dependent signal. However, fluorescent intensity measurements previously obtained across an entire test section when filled with stagnant water at maximum dye concentration showed that these effects and those of possible beam attenuation were negligible. Thus, unlike absorption techniques which measure integrated ("line of sight") properties, the receiving optical arrangement primarily governs the spatial resolution of the fluorescence technique. In the present experiments, off-axis light collection and multimode fiber aperture size resulted in a maximum focal volume length dimension of approximately 0.5 mm although smaller

spatial resolution can be achieved by appropriate choice of collection optics without affecting the LDV, since the velocimeter interference fringes are moving.

The instrument delivered to NASA LaRC comprises three primary elements namely: the optical, traverse control, and data acquisition systems. The Laser Optical system (Fig. 4) uses a 6 watt Argon-ion laser. The transmitting optical arrangement (Fig. 5) is straightforward with a few unique features addressing the common problem of beam distortion or thermal blooming at higher laser powers. Frequency shifting is done before the color separation prisms using a single Acousto-Optic modulator made of a selected flint glass which can handle the full laser power with minimal distortion. (For significantly greater laser powers we also found a water cooled Bragg cell which was more than required for this application.) This is followed by the color separation prisms, the first of which are fused silica for power handling capacity. A final prism of dense flint provides maximum dispersion once the laser beam has been split into at least eight beams. Final color selection is made using right angle prisms. The lines used for this experiment were 514.5nm, 488nm, 476.5nm. Although the emission spectrum of the fluorescence is centered about 515nm, there was sufficient higher wavelength emission for the edge filter to select wavelengths above 525nm. Other laser lines could have been selected if needed for fluorescence excitation or for separation from emission lines.

Pure fused-silica core single-mode polarization-preserving fibers are used for light transmission; two fibers per color. The use of optical fibers not only avoids the tedium of mirror-traverse alignment, but also greatly simplifies the transmission of light to the third axis. The pure fused silica core fibers seem immune from the progressive transmission losses which are found in other fibers. Polarization preserving fibers provide greater modal stability when the fibers are flexed or manipulated. For mechanical protection the fibers are armored and contained within a conduit. Upon exiting the fibers the beams are collimated at 4.4 mm dia. with a separation of 60mm.

Forward-scattered light is collected with a single 80mm diameter lens and focused into a 200 μm multi-mode optical fiber which conducts it to the color separation and signal detection box (Fig. 6). The fluorescence signal is split off with an edge filter. For maximum efficiency, a prism separation scheme is used rather than di-chroic filter and interference filters for the LDV signals.

Experience has shown that accurate positioning is vital to a successful test program. Accurate positioning is complicated in an air-window-water environment due to the difference in refractive indices on either side of the water tunnel window. Refractive index problems are particularly acute for the third component beams which are transmitted at 45 degrees to the normal incidence of the four beam axis. Fig. 7 illustrates the problem, which may not be intuitively obvious. In order to traverse the measurement volume horizontally some distance in water, the two orthogonal optical components, axial(x) and model vertical(z), must be traversed some lesser amount in air. In order for the third or off-axis pair of beams to intersect the measurement volume at the required angle of thirty degrees, the beams must strike the window at an angle of approximately forty-five degrees. In addition, linear horizontal movement of the focal point of the third pair of beams requires lens movement on a sloping line. As the lens gets closer to the window, it will have to rise. The length of movement on this sloping path will be only about half

of the movement of the focal point in water. Thus, two lens systems must move on different paths at different speeds in order to maintain a coincident focal point. To sort this all out, two three-axis traverse tables were installed for computer controlled, algorithm driven positioning of the LDV probe volume and forward scatter collecting optics. A fourth traversing axis supports the off-axis transmitted beams. Position is maintained by a custom designed eight axis traverse controller with micro-stepping drives, optical encoder feedback, and limit switch safety stops. Details of the Traverse Control System are given in Appendix A.

A Laser Velocimeter Data Acquisition System (LVDAS) has been designed. This instrument processes one to three channels of LDV data and digitizes up to four channels of analog data, one of which represents the concentration of dye. The instrument ensures coincidence and multiplexes the data to the computer. Velocities and analog channel values are displayed as well as data rates. Details are given in Appendix B.

Fig. 8 shows the modified data handling system for the 3D laser fluorescence anemometer. The continuous though not necessarily non-zero output from the high-gain current to voltage converter is fed directly into an analog to digital converter to provide 12 bits of digital information at 50 kHz. In water flows, this was more than sufficient to provide essentially real time point concentration data in digital form. But data from the three component LV system were not continuous wave since particle arrival times in the focal volume were random. However, whenever valid and essentially coincident data were received on all LV channels, a necessary requirement for shear stress measurement, these velocities, along with the instantaneous concentration voltage, was recorded. From a series of such readings, mean and turbulent velocity and concentration profiles were determined along with the turbulent shear stress and velocity/concentration cross correlations. These latter cross products provide new information on turbulent mixing rates in complex flows. In addition, we are able to determine details of the concentration/turbulent intermittencies from ensemble averages generated for selected instantaneous concentration levels. This will shed quantitative light on the turbulent structure and entrainment of fluid originating from different points in the flow field.

Experimental Details

Test Facility

The NASA Langley 16- by 24- Inch Water Tunnel is shown in Fig. 9. The tunnel has a vertical test section with an effective working length of about 4.5 ft. The velocity in the test section can be varied from 0 to 0.75 ft/sec., resulting in unit Reynolds numbers from 0 to $7.73 \times 10^4 \text{ ft}^{-1}$ based on a water temperature of 78°F. The normal test velocity yielding smooth flow is 0.25 ft/sec, resulting in unit Reynolds numbers of $2.58 \times 10^4 \text{ ft}^{-1}$ at 68°F. The model support system has deflection ranges of $\pm 33^\circ$ and $\pm 15^\circ$ in two planes of rotation. Rotation is accomplished via electronic remote control, and visual indicators allow the user to set angles within about $\pm 0.25^\circ$.

The fluorescence seeding method for this investigation used fluorescein dye injected into the jet flow field from inside the model. Naturally-occurring particles in the water were used for seeding for the LDV part of the instrument. A representative size distribution provided by NASA is shown in Fig. 10.

Model

The test model selected by NASA was a non-axisymmetric afterbody propulsion model and is shown in Fig. 11. This model is a scaled-down version of a model to be tested in the NASA Langley 16-foot Transonic Wind Tunnel as part of a computational fluid dynamics code validation study. It can be used to simulate nozzle exit velocity ratios typical of those in the wind tunnel study. It consists of a generic forebody with a non-axisymmetric boattail and nozzle. Water is injected into the interior of the model and exhausted through flow-conditioning foam ahead of the throat and exits through the nozzle. Fluorescein dye is introduced upstream of the model in the water supply tube for the jet.

Sample Test Results

Unfortunately, optical beam refraction problems caused by complex test section wall deformations ($\approx 0.10''$) under hydrostatic loading impeded laser beam alignment although limited data were obtained at zero angle of attack and a nominal jet exit to free stream velocity ratio of 1.5. For these test conditions, the cross flow velocity components (v, w) were negligible as expected. However, the model did provide a flowfield in which the basic instrument concepts could be verified. Fig. 12 shows the measured axial velocity distributions. It shows the extent of the jet and the location of the velocity defects in the model wall wakes. Clearly, the jet is highly skewed, but this has since been corrected by subsequent modification of the internal model flow treatment devices. The mean concentration profile (Fig. 13), which is much more symmetric, defines the extent of the jet at this axial station. Some insight into the unsteady features of the flow can be obtained from the fluctuation measurements. Figs. 14 and 15 show there is significant mixing in the outer shear layers between the jet and freestream flows. Peak fluctuation levels in the regions of maximum mean gradients indicate that small scale mixing is dominant. A measure of the level of streamwise mixing in the jet can be determined from Fig. 16. As expected the $\overline{u'c'}$ cross correlation function is positive as faster moving jet fluid is associated with higher concentration whereas fluid originating from the slower moving freestream has lower or zero concentration.

Concluding Remarks

A Laser Fluorescence Anemometer which comprises a three component laser doppler velocimeter system with a fourth channel to measure fluorescent dye concentration has been installed in the NASA Langley Water Tunnel. The system includes custom designed optics, data acquisition, and traverse control instruments and a custom software package.

Feasibility studies clearly demonstrate how water tunnels can be used in conjunction with advanced optical techniques to provide non-intrusive detailed flow field measurements of complex fluid flows with a minimum of expense. The measurements show that the laser fluorescence anemometer will provide new insight into the structure, entrainment and mixing of vortical and shear layer flows.

References

1. Sarpkaya, T., *Vortex Breakdown in Swirling Conical Flows* AIAA Journal, Vol. 9, No. 9, p. 1792, 1971.
2. Lambourne, N. C. and Bryer, D. W., *The Bursting of Leading Edge Vortices--some Observations and Discussion of the Phenomenon*. Brit. A.R.C.R. & M., No. 3282, 1962.
3. Faler, J. H. and Leibovich, S., *Disrupted States of Vortex Flow and Vortex Breakdown*. The Physics of Fluids, Vol. 20, No.9, p. 1385, 1977
4. Karashima, K. and Kitama, S., *The Effect of a Small Blowing on Vortex Breakdown of a Swirling Flow*. Computational Techniques & Applications: CTAC 83, p. 553, 1984.




Fig. 1 Laser Fluorescence Anemometer




Fig. 2 Relationship Between Dye Concentration and Fluorescent Output




Fig. 3 Fluorescence Sensitivity and Frequency Response




Fig. 4a Schematic of Plan View LFA System




Fig. 4b Schematic of Side View LFA System




Fig. 5 Sending Optics




Fig. 6 Receiving Optics




Fig. 7 The Refraction Problem




Fig. 8 Data Acquisition and Traverse Control Systems




Fig. 9 Langley 16- by 24- Inch Water Tunnel




Fig. 10 Size Distribution Particle Concentration




Fig. 11 Propulsion Model




Fig. 12 Axial Velocity Profile




Fig. 13 Concentration Profile




Fig. 14 Velocity Fluctuations




Fig. 15 Concentration Fluctuations




Fig. 16 Velocity-Concentration
Cross-Correlation

Appendix A. The Traverse Control System

The traverse control system is made up of four subsystems, see Fig. A1. The first subsystem is the main data taking computer, the HP 9000-330. The second subsystem, the TCS8 (Traverse Control System 8 Axis), receives high level traverse commands from the HP 9000. The full duplex serial communications that links these two subsystems allows the HP 9000 to monitor the position and status of each axis in the system, see Appendix A.3 TCS8's Serial Interface Command Descriptions. The TCS8 can also function as a stand alone traverse controller. Through the use of the TCS8's front panel, an operator can execute all of the commands that the HP 9000 can, plus the operator can control all axes in jog mode, see Appendix A.1 TCS8's Front Panel Descriptions and Appendix A.2 TCS8's Local Command Descriptions. The third subsystem, the MDS (Motor Drive System), is controlled solely by the TCS8. The TCS8 translates the high level commands from the HP 9000 and its front panel into low level indexer commands, see The Compumotor AX Drive User Manual previously delivered. The TCS8 also receives encoder pulses from the traverses via the MDS. This allows the TCS8 to display realtime position information on its front panel. The fourth and final subsystem of the traverse control system is the slide, motor, encoder, and limit switches that make up each axis. A drawing of each cable that is used to connect the traverse control system is included in Appendix A.4 Traverse Control System Cables.

The TCS8

The TCS8 is a microprocessor controlled system designed to interface an operator with a traverse system. The operator can utilize the TCS8 through the front panel, see Appendix A.1 TCS8's Front Panel Descriptions and Appendix A.2 TCS8's Local Command Descriptions, and/or with one or two host computers over serial interfaces, see Appendix A.3 TCS8's Serial Interface Command Descriptions. The TCS8 stores all the critical parameters of motion, for each of the eight axes that it controls, in non-volatile memory. The critical parameters of motion being: position, encoder counts per unit travel, encoder counts per motor revolution, velocity, and acceleration. All of these parameters may be viewed, set, and saved. The TCS8 has three modes of motion. They are absolute, relative, and jogged. With absolute movements, the operator specifies the final location; with relative movements, a distance is specified; and with jogged movements the operator presses a jog key on the front panel of the TCS8 until the desired location is obtained.

The Motor Drive System

Each of the two MDS's have 4 indexer/drivers contained within them. The TCS8 communicates with the indexers in the MDS's over a closed loop serial daisy chain. When two MDS's are used, as in this setup, the first MDS in the chain must be set to 8 and the second set to 4. By setting the first MDS to 8, the operator is opening the closed loop serial daisy chain allowing the second MDS to be included in the chain. The 4/8 switch is located on the back panel of each MDS, see Fig. A2. This figure also shows the location of all the motor, limit, and encoder

connections. Channels X1, X2, Y1, and Y2 of the TCS8 control axis 1 through 4 on the first MDS and channels Z1, Z2, A1, and A2 control axis 1 through 4 on the second MDS. The TCS8 Encoders connector on the back of each MDS has a corresponding connector of the back of the TCS8, see Fig. A3 Schematic of TCS8 Back Panel. The interconnecting cable is detailed in Appendix A.4 Traverse Control System Cables.

Positioning Resolution

The indexer/drivers that are used in the MDS can drive the motors at 12,800 steps/revolution. The encoder used on each axis are 1000 pulses/revolution with quadrature encoding. Quadrature encoding adds a factor of 4 to the number of pulses/revolution to make this number 4000 pulses/revolution. This number, 4000 pulses/revolution, is well within the limit of 12,800 steps/revolution set by the indexer. The final factor in the product of the resolution of an axis is the number of threads/inch of the lead screw. All of the axes of the traverse system, except the auxiliary axis, have lead screws of 5 threads/inch, the auxiliary axis has a lead screw of 10 threads/inch. So the positioning resolution of the axes with a 5 threads/inch lead screw is 0.00005 inches and the auxiliary axis has a resolution of 0.000025 inches.




Fig. A1 Langley Traverse Control System




Fig. A2 Schematic of Motor Drive System Back Panel





Fig. A3 Schematic of TCS8 Back Panel

Appendix A.1 TCS8's Front Panel Description



1. POSITION DISPLAY WINDOWS

There are eight windows corresponding to the eight axes that the TCS8 is capable of controlling. The position of each axis is continuously updated, by monitoring its encoder, and displayed in a fixed format of a sign, two digits, a decimal point, and four digits.

2. POWER KEY

The power key is used to store the current configuration to nonvolatile memory before turning off power to the TCS8. Pressing the power key turns the displays off and saves the current configuration. Pressing it again turns the displays back on. This key can be used to implement a screen saver function.

3. JOG CONTROL KEYS

These keys are used to control up to eight axes in a jog mode. The mode (slaved, one's only, or two's only) can be set through the jog menu. When the operator presses a jog key, the respective axis will begin to move. The direction that the axis moves is determined by the operator pressing either a plus or minus jog key. A plus jog key will turn the lead screw in a clockwise direction (away from the motor), a minus jog key will turn it in the counter-clockwise direction (towards the motor). By releasing the jog key the operator stops motion on that axis. Motion will also stop, if the axis reaches the limit for the direction it is moving, or if the indexer determines that the axis has stalled.

4. SCROLL KEYS

These keys are used to scroll items through the MENU, COMMAND, and CHANNEL windows. All of the menus, their commands, and channel variations will be detailed in another appendix.

5. COMMAND WINDOWS

These three windows (MENU, COMMAND, and CHANNEL) are used, in tandem with their respective scroll keys, to formulate a command to be executed by the TCS8.

6. EXECUTE KEY

This key is used to execute the command currently formulated in the MENU, COMMAND, and CHANNEL windows.

7. DATA WINDOW

Many of the TCS8's commands require some added data, e.g. the distance to move or an axis' encoder counts per unit. Data for these commands is entered from the numeric key pad on the lower right of the TCS8 into the DATA window. Only a valid real number can be entered into the DATA window. If the operator enters an invalid real number the character that is invalid will flash until the operator presses backspace or a valid character.

8. STOP KEY

The stop key, when pressed, will stop motion on all axes. The TCS8 will not loose track of the position of any axis. A move command started by the host computer and stopped by the stop key will finish normally with the position being reported. The position reported is the instantaneous position when the stop key was pressed. The final position of the axis being moved could be different than what was reported thus the host computer should read the position again after a panic stop.

9. STATUS WINDOW

The STATUS window reflects the result of all commands. For commands that are not instantaneous, this window displays a busy status and then when the command completes it displays a ready status. The results of all view commands are displayed in the STATUS window. The STATUS window also displays the activity over the COM interfaces. For example, when the command for viewing position is sent over the COM1 interface, the STATUS window will display "COM1 VP" and when the command is completed the window will display "COM1 vp".

10. NUMERIC KEY PAD

The numeric key pad is used to enter a number into the data window. The user may backspace in the window or clear (shift-backspace) the window. Only a valid real number can be entered into the data window. If the operator enters an invalid real number the character that is invalid will flash until the operator presses backspace or a valid character.

Appendix A.2 TCS8's Local Command Descriptions

This appendix describes the command set that can be executed from the front panel of the TCS8. Using the up and down keys under the MENU, COMMAND, and CHANNEL windows, the operator can formulate a command and then execute it by pressing the EXECUTE key. Some commands require extra information to be entered into the DATA window through the use of the numeric key pad. Each description includes a list of related commands that should be refer to to enhance the operator's understanding of the command. Also where applicable, the default setting is given.

MENU: MOVE

COMMAND: TO ZERO

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The MOVE TO ZERO command is an easy way to move some or all of the axes to the zero position. This command can also be accomplished with the MOVE ABSOLUTE command and a zero in the DATA window. Before using this command, selected axes must be initialized with the INIT DRIVE ON command. This command can be canceled by pressing the STOP key. When the STOP key is pressed, all axes will stop immediately. If an axis encounters a limit before reaching zero, the rest of its movement is aborted.

RELATED COMMANDS: MOVE ABSOLUTE, MOVE RELATIVE, INIT Drive ON

MENU: MOVE

COMMAND: ABSOLUTE

CHANNELS: X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The MOVE ABSOLUTE command requires a position to be entered in the DATA window. This position and the current position of the axis is used to calculate the relative distance the axis must move. Before using this command, selected axes must be initialized with the INIT DRIVE ON command. This command can be canceled by pressing the STOP key. When the STOP key is pressed, all axes will stop immediately. If an axis encounters a limit before reaching the position entered in the DATA window, the rest of its movement is aborted.

RELATED COMMANDS: MOVE TO ZERO, MOVE RELATIVE, INIT Drive ON

MENU: MOVE

COMMAND: RELATIVE

CHANNELS: X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The MOVE RELATIVE command requires a distance to be entered in the DATA window. This position is used to calculate the relative distance the axis must move. Before using this command, selected axes must be initialized with the INIT DRIVE ON command. This command can be canceled by pressing the STOP key. When the STOP key is pressed, all axes will stop immediately. If an axis encounters a limit before moving the distance entered in the DATA window, the rest of its movement is aborted.

RELATED COMMANDS: MOVE TO ZERO, MOVE ABSOLUTE, INIT Drive ON

MENU: JOG

COMMAND: MODE

CHANNELS: SLAVED, ONE'S, TWO'S

DESCRIPTION: The JOG MODE command sets the way the JOG keys operate. When SLAVED is the setting, both the one and two axis of the X, Y, Z, or A coordinate will move the same amount. When ONE'S is the setting, only the one axes of the X, Y, Z, or A coordinate will move. And finally, when TWO'S is the setting, only the two axes of the X, Y, Z, or A coordinate will move. The current mode is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: SLAVED

MENU: SET

COMMAND: CPU

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CPU command allows the user to change the counts per unit travel. The CPU for an axis is determined by multiplying the encoder resolution (counts/revolution) by the lead screws resolution (revolutions/unit of travel). A units conversion can be added here to change for example from inches to centimeters. When the CPU for an axis is changed, the position is automatically converted. This command requires a value to be entered in the DATA window.

RELATED COMMANDS: SET CPR, SET POSITION

DEFAULT:	X1	20000
	X2	20000
	Y1	20000
	Y2	20000
	Z1	20000
	Z2	20000
	A1	40000
	A2	40000

MENU: SET

COMMAND: CPR

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CPR command allows the user to change the encoder counts per motor revolution. The CPR for an axis is determined by dividing the encoder resolution (counts/revolution) by the lead screws resolution (revolutions/unit of travel). The encoder counts per motor revolution, that is entered in the DATA window, must be a positive integer.

RELATED COMMANDS: SET CPU

DEFAULT:	X1	4000
	X2	4000
	Y1	4000
	Y2	4000
	Z1	4000
	Z2	4000
	A1	4000
	A2	4000

MENU: SET

COMMAND: POSITION

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET POSITION command allows the user to change the current position of an axis. The new position must be entered in the DATA window before executing the command.

RELATED COMMANDS: SET CPU

MENU: SET

COMMAND: VELOCITY

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET VELOCITY command allows the user to change the maximum speed at which an axis will travel. The range of valid velocities is 0.002 to 50.000 revolutions per second. The default is 5 revs/sec. An axis may stall at velocities higher than the default. The new velocity must be entered in the DATA window before executing the command.

RELATED COMMANDS: SET ACCEL.

DEFAULT:	X1	5.000
	X2	5.000
	Y1	5.000
	Y2	5.000
	Z1	5.000
	Z2	5.000
	A1	5.000
	A2	5.000

MENU: SET

COMMAND: ACCEL.

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET ACCEL. command allows the user to change the maximum acceleration for an axis. The range of valid accelerations is 0.01 to 999.99 revolutions per second per second. The default is 5 revs/sec/sec. The new acceleration must be entered in the DATA window before executing the command.

RELATED COMMANDS: SET VELOCITY

DEFAULT:	X1	5.00
	X2	5.00
	Y1	5.00
	Y2	5.00
	Z1	5.00
	Z2	5.00
	A1	5.00
	A2	5.00

MENU: SET

COMMAND: CrntsOn

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CrntsOn command allows the user to turn the motor currents on. The motor current must be on for an axis to be moved. The information in the DATA window is ignored.

RELATED COMMANDS: SET CrntsOff

MENU: SET

COMMAND: CrntsOff

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CrntsOff command allows the user to power down motors when they will not be used for long periods of time. The information in the DATA window is ignored.

RELATED COMMANDS: SET CrntsOn

MENU: SET

COMMAND: INITS ON

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET INITS ON command allows the user to initialize the indexers without turning on the power to the motors. The information in the DATA window is ignored.

RELATED COMMANDS: INIT Drive ON

MENU: VIEW

COMMAND: Cnt/Unit

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW Cnt/Unit command displays the current setting of the encoder counts per unit travel parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPU

MENU: VIEW

COMMAND: Cnt/MRev

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW Cnt/MRev command displays the current setting of the encoder counts per motor revolution parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPR

MENU: VIEW

COMMAND: VELOCITY

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW VELOCITY command displays the current setting of the velocity parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET VELOCITY

MENU: VIEW

COMMAND: ACCEL.

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW ACCEL. command displays the current setting of the acceleration parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET ACCEL.

MENU: VIEW

COMMAND: INIT

CHANNELS: none

DESCRIPTION: The VIEW INIT command uses the STATUS window to display a one(initialized) or a zero(uninitialized) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of: X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: SET INITS, INIT Drive ON

MENU: VIEW

COMMAND: CURRENTS

CHANNELS: none

DESCRIPTION: The VIEW CURRENTS command uses the STATUS window to display a one(current on) or a zero(current off) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of: X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: SET CrntsOn, SET CrntsOff, INIT Drive ON, INIT Drive OFF

MENU: VIEW

COMMAND: Plus LMT

CHANNELS: none

DESCRIPTION: The VIEW Plus LMT command uses the STATUS window to display a one(on limit) or a zero(not on limit) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of: X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: none

MENU: VIEW

COMMAND: Minus LMT

CHANNELS: none

DESCRIPTION: The VIEW Minus LMT command uses the STATUS window to display a one(on limit) or a zero(not on limit) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of: X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: none

MENU: VIEW

COMMAND: HOME

CHANNELS: none

DESCRIPTION: The VIEW HOME command uses the STATUS window to display a one(on limit) or a zero(not on limit) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of: X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: none

MENU: VIEW

COMMAND: STALL

CHANNELS: none

DESCRIPTION: The VIEW STALL command uses the STATUS window to display a one(stalled) or a zero(not stalled) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of: X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: none

MENU: INIT

COMMAND: DEFAULT

CHANNELS: none

DESCRIPTION: The INIT DEFAULT command restores the initial factory defaults (CPU, CPR, VELOCITY, ACCELERATION, BAUD RATE, BITS/CHAR, PARITY, STOP BITS, HANDSHAKE) of the TCS8. After executing this command, execute the command INIT Drive ON to initialize the indexers. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPU, SET CPR, SET VELOCITY, SET ACCEL.

MENU: INIT

COMMAND: Drive ON

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The INIT Drive ON command initializes the selected axes for movement. After executing this command the currents are on to the motors. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPU, SET CPR, SET VELOCITY, SET ACCEL., SET CntrsOn, SET CntrsOff, INIT DEFAULT

MENU: INIT

COMMAND: Drive OFF

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The INIT Drive OFF command is an alias for SET CntrsOff.

RELATED COMMANDS: SET CntrsOff

MENU: COM1/COM2

COMMAND: BaudRate

CHANNELS: 19.2K, 9600, 4800, 2400, 1200, 300, 110

DESCRIPTION: The COM1/COM2 BaudRate command set the baud rate for the selected communication channel. The information in the DATA window is ignored. The current baud rate is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: 9600

MENU: COM1/COM2

COMMAND: Bit/Char

CHANNELS: SEVEN, EIGHT

DESCRIPTION: The COM1/COM2 Bit/Char command set the bits per character for the selected communication channel. The information in the DATA window is ignored. The current number of bits per character is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: EIGHT

MENU: COM1/COM2

COMMAND: Parity

CHANNELS: NONE, EVEN, ODD

DESCRIPTION: The COM1/COM2 Parity command set the parity for the selected communication channel. The information in the DATA window is ignored. The current parity is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: EVEN

MENU: COM1/COM2

COMMAND: StopBits

CHANNELS: 1, 1.5, 2

DESCRIPTION: The COM1/COM2 StopBits command set the stop bits for the selected communication channel. The information in the DATA window is ignored. The current number of stop bits is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: 1

MENU: COM1/COM2

COMMAND: HandShak

CHANNELS: NO, YES

DESCRIPTION: The COM1/COM2 HandShak command set the handshake for the selected communication channel. The information in the DATA window is ignored. An asterisk marks whether there is handshaking or not.

RELATED COMMANDS: none

DEFAULT: YES

Appendix A.3 TCS8's Serial Interface Command Descriptions

This appendix describes the command set that can be executed through the serial interfaces of the TCS8. Each description includes a code section that outlines the characters that must be sent to execute the command. The vertical bar in this section is used as a separator and is not sent as part of the command code. The symbol "CRLF" stands for the two characters carriage return and line feed. Also where applicable, the default setting is given.

COMMAND: CHANGE SERIAL CONFIGURATION

CODE: CS COM;CATEGORY;ATTRIBUTE;

PARAMETERS:	COM:	1/COM1 2/COM2
	CATEGORY:	0/BAUDRATE
	ATTRIBUTE:	0/19.2K 1/9600 2/4800 3/2400 4/1200 5/300 6/110
	CATEGORY:	1(BITS PER CHARACTER)
	ATTRIBUTE:	0/SEVEN 1/EIGHT
	CATEGORY:	2(PARITY)
	ATTRIBUTE:	0/NONE 1/EVEN 2/ODD
	CATEGORY:	3(STOP BITS)
	ATTRIBUTE:	0/ONE 1/ONE AND A HALF 2/TWO
	CATEGORY:	4(HANDSHAKE)
	ATTRIBUTE:	0/NO 1/YES

DESCRIPTION: This command must be executed with extreme caution and thought. If the user changes an attribute of the same COM port that he is sending the command, he must change to that attribute on the host computer before sending the next command. The best way to change the serial configuration of a COM port is to utilize the front panel commands.

DEFAULT: 9600 baud, EIGHT bits/char, EVEN parity, ONE stop bit, handshaking YES

EXAMPLE: To change the baudrate of COM1 to 2400 the user must send CS1;0;3;

COMMAND: MOVE TO ABSOLUTE POSITION AND REPORT FINAL POSITION

CODE: MA CHANNEL:POSITION,CHANNEL:POSITION,...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 POSITION: Real number free format

DESCRIPTION: This command moves selected channels to absolute positions.

EXAMPLES:

To move all channels to zero the user may send MA0:0,CRLF or MA12345678:0,CRLF

To move channel X1 to zero the user must send MA1:0,CRLF

To move channels X1 and X2 to zero the user may send MA12:0,CRLF or MA1:0,2:0,CRLF or MA1:0,CRLF and MA2:0,CRLF

COMMAND: MOVE TO RELATIVE DISTANCE AND REPORT FINAL POSITION

CODE: MR CHANNEL:DISTANCE,CHANNEL:DISTANCE,...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 POSITION: Real number free format

DESCRIPTION: This command moves selected channels relative distances.

EXAMPLES:

To move all channels one unit the user may send MR0:1,CRLF or MR12345678:1,CRLF

To move channel X1 one unit the user must send MR1:1,CRLF

To move channels X1 and X2 one unit the user may send MR12:1,CRLF or MR1:1,2:1,CRLF or MR1:1,CRLF and MR2:1,CRLF

COMMAND: SET ACCELERATION

CODE: SA CHANNEL:ACCELERATION,CHANNEL:ACCELERATION,...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 ACCELERATION: Real number free format between
 0.01 and 99.99 inclusive.

DESCRIPTION: This command sets the acceleration for selected channels.

DEFAULT: All channels 5.00 revolutions/second/second

EXAMPLES:

To set the acceleration for all channels to 4.00 revolutions/second/second the user may send SA0:4.00,CRLF or SA12345678:4.00,CRLF

To set the acceleration for channel X1 to 4.00 revolutions/second/second the user must send SA1:4.00,CRLF

To set the acceleration for channels X1 and X2 to 4.00 revolutions/second/second the user may send SA12:4.00,CRLF or SA1:4.00 ,2:4.00,CRLF or SA1:4.00,CRLF and SA2:4.00,CRLF

COMMAND: VIEW ACCELERATION

CODE: VA CHANNEL|CHANNEL...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command views the acceleration for selected channels. The TCS8 transmits each of the accelerations requested back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the acceleration for all channels the user may send VA0CRLF or VA12345678CRLF

To view the acceleration for channel X1 the user must send VA1CRLF

To view the acceleration for channels X1 and X2 the user may send VA12CRLF or VA1CRLF and VA2CRLF

COMMAND: SET VELOCITY

CODE: SV CHANNEL:VELOCITY,CHANNEL:VELOCITY,...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 VELOCITY: Real number free format between
 0.001 and 50.000 inclusive.

DESCRIPTION: This command sets the velocity for selected channels.

DEFAULT: All channels 5.000 revolutions/second

EXAMPLES:

To set the velocity for all channels to 4.00 revolutions/second the user may send SV0:4.00,CRLF or SV12345678:4.00,CRLF

To set the velocity for channel X1 to 4.00 revolutions/second the user must send SV1:4.00,CRLF

To set the velocity for channels X1 and X2 to 4.00 revolutions/second the user may send SV12:4.00,CRLF or SV1:4.00 ,2:4.00,CRLF or SV1:4.00,CRLF and SV2:4.00,CRLF

COMMAND: VIEW VELOCITY

CODE: VV CHANNEL|CHANNEL...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command views the velocity for selected channels. The TCS8 transmits each of the velocities requested back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the velocity for all channels the user may send VV0CRLF or VV12345678CRLF

To view the velocity for channel X1 the user must send VV1CRLF

To view the velocity for channels X1 and X2 the user may send VV12CRLF or VV1CRLF and VV2CRLF

COMMAND: SET ENCODER COUNTS PER UNIT TRAVEL

CODE: SU CHANNEL:CPU,CHANNEL:CPU,...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 CPU: Non-zero real number free format.

DESCRIPTION: This command sets the encoder counts per unit travel for selected channels.

DEFAULT: X1,2,Y1,Y2,Z1,Z2 20000 counts/inch and A1,A2 40000 counts/inch

EXAMPLES:

To set the encoder counts per unit travel for all channels to 5000 the user may send

SU0:5000,CRLF or SU12345678:5000,CRLF

To set the encoder counts per unit travel for channel X1 to 5000 the user must send

SU1:5000,CRLF

To set the encoder counts per unit travel for channels X1 and X2 to 5000 the user may send

SU12:5000,CRLF or SU1:5000 ,2:5000,CRLF or SU1:5000,CRLF and SU2:5000,CRLF

COMMAND: VIEW ENCODER COUNTS PER UNIT TRAVEL

CODE: VU CHANNEL|CHANNEL...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command views the encoder counts per unit travel for selected channels. The TCS8 transmits each of the encoder counts per unit travel requested back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the encoder counts per unit travel for all channels the user may send VU0CRLF or

VU12345678CRLF

To view the encoder counts per unit travel for channel X1 the user must send VU1CRLF

To view the encoder counts per unit travel for channels X1 and X2 the user may send VU12CRLF or VU1CRLF and VU2CRLF

COMMAND: SET ENCODER COUNTS PER MOTOR REVOLUTION

CODE: SR CHANNEL:CPR,CHANNEL:CPR,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
CPU: Non-zero integer free format.

DESCRIPTION: This command sets the encoder counts per motor revolution for selected channels.

DEFAULT: X1,2,Y1,Y2,Z1,Z2 and A1,A2 4000 counts/inch

EXAMPLES:

To set the encoder counts per motor revolution for all channels to 500 the user may send SR0:500,CRLF or SR12345678:500,CRLF

To set the encoder counts per motor revolution for channel X1 to 500 the user must send SR1:500,CRLF

To set the encoder counts per motor revolution for channels X1 and X2 to 500 the user may send SR12:500,CRLF or SR1:500,2:500,CRLF or SR1:500,CRLF and SR2:500,CRLF

COMMAND: VIEW ENCODER COUNTS PER MOTOR REVOLUTION

CODE: VR CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command views the encoder counts per motor revolution for selected channels. The TCS8 transmits each of the encoder counts per motor revolution requested back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the encoder counts per motor revolution for all channels the user may send VR0CRLF or VR12345678CRLF

To view the encoder counts per motor revolution for channel X1 the user must send VR1CRLF

To view the encoder counts per motor revolution for channels X1 and X2 the user may send VR12CRLF or VR1CRLF and VR2CRLF

COMMAND: SET POSITION

CODE: SP CHANNEL:POSITION,CHANNEL:POSITION,...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 POSITION: real number.

DESCRIPTION: This command sets the position for selected channels.

EXAMPLES:

To set the position for all channels to 1.5 the user may send SP0:1.5,CRLF or SP12345678:1.5,CRLF

To set the position for channel X1 to 1.5 the user must send SP1:1.5,CRLF

To set the position for channels X1 and X2 to 1.5 the user may send SP12:1.5,CRLF or SP1:1.5,2:1.5,CRLF or SP1:1.5,CRLF and SP2:1.5,CRLF

COMMAND: VIEW POSITION

CODE: VP CHANNEL|CHANNEL...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command views the position for selected channels. The TCS8 transmits each of the positions requested back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the position for all channels the user may send VP0CRLF or VP12345678CRLF

To view the position for channel X1 the user must send VP1CRLF

To view the position for channels X1 and X2 the user may send VP12CRLF or VP1CRLF and VP2CRLF

COMMAND: SET CURRENT TO MOTOR WINDINGS

CODE: SC CHANNEL:ON/OFF,CHANNEL:ON/OFF,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2
 ON/OFF: 1/ON
 0/OFF

DESCRIPTION: This command sets the current to the motor windings for selected channels on or off.

EXAMPLES:

To set the current to the motor windings for all channels on the user may send SC0:1,CRLF or SC12345678:1,CRLF to set them off the user may send SC0:0,CRLF or SC12345678:0,CRLF

COMMAND: VIEW CURRENT TO MOTOR WINDINGS

CODE: VC CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command views the current to the motor windings for selected channels. The TCS8 transmits each response of on/off (1/0) back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the current to the motor windings for all channels the user may send VC0CRLF or VC12345678CRLF

To view the current to the motor windings for channel X1 the user must send VC1CRLF

To view the current to the motor windings for channels X1 and X2 the user may send VC12CRLF or VC1CRLF and VC2CRLF

COMMAND: SET INITIALIZATION OF INDEXER/DRIVERS

CODE: SI CHANNEL|CHANNEL...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command sends the current value of the acceleration, velocity, and the encoder counts per motor revolution to the indexer/driver for the selected channels. This command must be sent before any move commands may be sent.

EXAMPLES:

To initialize all channels the user may send SI0CRLF or SI12345678CRLF

To initialize channel X1 the user must send SI1CRLF

To initialize channels X1 and X2 the user may send SI12CRLF or SI1CRLF and SI2CRLF

COMMAND: VIEW INITIALIZATION OF INDEXER/DRIVERS

CODE: VI CHANNEL|CHANNEL...|CRLF

PARAMETERS: **CHANNEL:** 0/ALL CHANNELS
 1/X1
 2/X2
 3/Y1
 4/Y2
 5/Z1
 6/Z2
 7/A1
 8/A2

DESCRIPTION: This command returns "1" if the indexer/driver has been initialized since the TCS8 was turned on and "0" if it has not. The TCS8 transmits each of the responses back to the host computer separated by carriage return line feeds.

EXAMPLES:


To check the initialization of all channels the user may send VI0CRLF or VI12345678CRLF

To check the initialization of channel X1 the user must send VI1CRLF

To check the initialization of channels X1 and X2 the user may send VI12CRLF or VI1CRLF and VI2CRLF

Appendix A.4

Traverse Control System Cables



COMPLERE INC.

NASA LaRC LASER FLUORESCENCE ANEMOMETER LFA TRAVERSE SYSTEM TRAVERSE SYSTEM CABLE ROUTING DIAGRAM

DRAWING DATE

JULY 8, 1991

DESIGN ENGINEER

TODD A. AMBUR

FILE NAME


LANGLEY TRAVERSE DIAGRAM

DRAWING NUMBER

LFA-TAA-001

**COMPLERE INC.
MOTOR DRIVE SYSTEM**

**COMPUMOTOR
STEPPER MOTOR**



SIGNAL	DESCRIPTION
A+	Motor Winding
A-	Motor Winding
B+	Motor Winding
B-	Motor Winding
SHIELD	Motor Case Ground
NC	No Connection


COMPLERE INC.

**NASA LaRC LASER FLUORESCENCE ANEMOMETER
LFA TRAVERSE SYSTEM
MOTOR DRIVE SYSTEM TO COMPUMOTOR STEPPER MOTOR**

DRAWING DATE JULY 3, 1991	FILE NAME LANGLEY MOTOR
DESIGN ENGINEER TODD A. AMBUR	DRAWING NUMBER LFA-TAA-002

**COMPLERE INC.
MOTOR DRIVE SYSTEM**

**DYNAMICS RESEARCH
ENCODER**



SIGNAL	DESCRIPTION
A+	Quadrature Encoder Signal
A-	Logical Complement of A+
B+	Quadrature Encoder Signal
B-	Logical Complement of B+
Z+	Once Per Revolution
Z-	Logical Complement of Z+
SHIELD	Case Ground
GROUND	Logical Ground
NC	No Connection


COMPLERE INC.

**NASA LaRC LASER FLUORESCENCE ANEMOMETER
LFA TRAVERSE SYSTEM
MOTOR DRIVE SYSTEM TO TCS8 ENCODER SIGNALS**

DRAWING DATE JULY 3, 1991	FILE NAME LANGLEY ENCODER
DESIGN ENGINEER TODD A. AMBUR	DRAWING NUMBER LFA-TAA-003

**COMPLERE INC.
MOTOR DRIVE SYSTEM**

**LINEAR INDUSTRIES
LIMIT SWITCHES**



SIGNAL	DESCRIPTION
HOME	Home Switch Signal
CW	End of Travel Limit Signal Clockwise
CCW	End of Travel Limit Signal Counter Clockwise
COMMON	Logical Ground


COMPLERE INC.

**NASA LaRC LASER FLUORESCENCE ANEMOMETER
LFA TRAVERSE SYSTEM
MOTOR DRIVE SYSTEM TO LINEAR INDUSTRIES LIMIT SWITCHES**

DRAWING DATE JULY 3, 1991	FILE NAME LANGLEY LIMIT SWITCH
DESIGN ENGINEER TODD A. AMBUR	DRAWING NUMBER LFA-TAA-004

COMPLERE INC.
MOTOR DRIVE SYSTEM

TCS8



C1 - CHANNEL 1
C2 - CHANNEL 2
C3 - CHANNEL 3
C4 - CHANNEL 4

COMPLERE INC.

NASA LaRC LASER FLUORESCENCE ANEMOMETER
LFA TRAVERSE SYSTEM
MDS TO TCS8 ENCODER SIGNALS

DRAWING DATE

JULY 3, 1991

FILE NAME

ENCODER SIGNALS

DESIGN ENGINEER


TODD A. AMBUR

DRAWING NUMBER

LFA-TAA-005

HP 9000-330

TCS8
COM1/COM
2




SIGNAL	DESCRIPTION
TxD	Transmit Data
RxD	Receive Data
RTS	Ready to Send
CTS	Clear to Send
DCD	Data Carrier Detect
DTR	Data Terminal Ready
DSR	Data Set Ready
RI	Ring Indicator
GND	Logical Ground
NC	No Connection

COMPLERE INC.

NASA LaRC LASER FLUORESCENCE ANEMOMETER
LFA TRAVERSE SYSTEM
HP SERIES 9000 MODEL 330 TO TCS8

DRAWING DATE JULY 3, 1991	FILE NAME HOST SERIAL
DESIGN ENGINEER TODD A. AMBUR	DRAWING NUMBER LFA-TAA-006



SIGNAL	DESCRIPTION
TxD	Transmit Data
RxD	Receive Data
RTS	Ready to Send
CTS	Clear to Send
GND	Logical Ground
NC	No Connection

COMPLERE INC.

NASA LaRC LASER FLUORESCENCE ANEMOMETER
LFA TRAVERSE SYSTEM
TCS8 TO MDS SERIAL

DRAWING DATE	FILE NAME
JULY 3, 1991	TCS8 SERIAL
DESIGN ENGINEER	DRAWING NUMBER
TODD A. AMBUR	LFA-TAA-007

Appendix B Laser Velocimeter Data Acquisition System

The LVDAS acquires simultaneous digital data, analog data, and time information data. The data are sampled, multiplexed, buffered, and then transferred to the facility's host computer for further data reduction, analysis, and presentation.

Four 16 bit word parallel input ports are provided to accept the digital output of LV counter processors and/or other instrumentation.

New applications in laser velocimetry have brought about the need for a more advanced laser velocimeter data acquisition system. These new applications require high data rates that are not hindered by on-line time dependent data sorting and real time graphic data presentation. The new Laser Velocimeter Data Acquisition System (LVDAS) was designed specifically to meet these advanced requirements.

High data acquisition rates are achieved by providing a separate latched input for each laser velocimeter digital input and a separate converter for each laser fluorescence analog input. The system will allow for a data acquisition rates of approximately 100,000 samples per second simultaneously on each of the laser velocimeter and laser fluorescence inputs.

A 32 bit time of day (TOD) 10MHz counter is used to tag arrival times to acquired digital LDV data as they become available on each of four digital inputs. When a data valid sync pulse is sensed for a particular channel, the LVDAS latches the current TOD into a 32 bit time of arrival register (TOA). A separate TOA register is available for each digital input so that particle arrival times of measured velocity information U, V, W can be monitored for coincidence. The latched times of arrivals have a resolution of 100 ns and maximum time of over 7 minutes.

All of the acquired digital velocity data with corresponding time of arrival data can be processed and stored. However, if coincident data is required, then the arrival time of the various channels can be conditionally accepted if they all occur within a finite window of time. These coincident events can then be assigned interarrival times which represent elapsed time since the previous event.

The coincidence control logic allows for 3 channel coincidence. The coincidence time is adjustable from 0.1 μ s to 1 s. In addition to the laser velocimeter inputs, three additional data words are generated internally. They are the inter arrival time, the coincidence time, and status words. The inter arrival and coincidence time is provided by a clock whose resolution and maximum time is 100 ns and 500 seconds respectively. The status word contains information about coincidence and the order in which the laser velocimeter data arrived.

During data acquisition, it is important that the user obtain some visual feedback about the data being acquired. This is necessary so that the user can make informed decisions about both the quality and quantity of data received. The user is either reassured about the quality of the data or can make alterations and improvements in technique while on line. To help achieve this, the instantaneous velocities are used to generate real time histograms from which probability density distributions are determined for all velocity components.

Additionally, the laser velocimeter data acquisition system has the capability of reducing the raw laser velocimeter data. Each laser velocimeter output contains the information required to

calculate the instantaneous velocities. From the instantaneous velocity determinations, the average velocities, turbulence levels, and the turbulence cross correlations are all be calculated.

The coincidence control logic will allow for up to 4 channel coincidence. The coincidence time can be adjustable to any resolution or duration within the capability of the time of arrival registers. When coincident criteria are met, the analog inputs can be sampled and converted to provide concurrent data with the digital data. A single time of arrival is latched for all four of the analog to digital inputs since they are all sampled and converted simultaneously. A final time of arrival is latched for external events. These might be derived from such sources as oscillating models or model surfaces, rotating helicopter blades, rotating engine fans, or flow sensors.

All digital Macrodyne data, optional digital data, analog to digital data, and time of arrival data can be sent by the LVDAS to other computers via two serial and two parallel input/output ports. One parallel port will be used for the HP series 9000 model 330 computer while the other can be used by the facility host computer. The serial ports can be used by PC type computers such as IBMs or MACs. Software has been developed for on-line data acquisition and display. A program listing is enclosed.

```

100 Main: -----! LDVWT
110 !
120 !           NASA LANGLEY RESEARCH CENTER           Property of COMPLERE INC.
130 !           16 BY 24 INCH WATER TUNNEL             Proprietary software
140 !                                           Copyright February 25, 1992
150 !           LASER FLUORESCENCE ANEMOMETER           Developed by: T. Kevin McDevitt
160 !
170 ! -----!
180 !
190 ! PROGRAM DESCRIPTION:
200 !
210 !     This program provides the capability to acquire simultaneous Laser Doppler Velocimeter (LDV), Laser
220 !     Fluorescence Anemometer (LFA), and Analog Voltage Data at user selectable traverse controlled probe volume
230 !     positions within the water tunnel flow.
240 !     The LVDAS (Laser Velocimeter Data Acquisition System) is used to sample the LDV, LFA, and Analog Voltage
250 !     data simultaneously with a coincidence criterion being applied to LDV incoming data. The LVDAS also generates
260 !     interarrival times and coincidence time.
270 !     The measured LDV data provides the necessary frequency information from which three components of flow
280 !     velocities can be determined. These velocities are measured directly in "LASER" coordinates. Coordinate
290 !     system transformations are applied to these measured velocities to obtain velocities in "TUNNEL" and "MODEL"
300 !     coordinates.
310 !     The TCS8 (Traverse Control System) is used to precisely move the LDV probe volume within the tunnel and
320 !     about the model. The TCS8 provides three axes plus one auxiliary axis of traverse capability for both the
330 !     transmitting and receiving side optical packages. The Tx and Rx side traverses can be moved independently to
340 !     achieve laser alignment or they can be moved together to maintain laser alignment.
350 !     The TCS8 will give the traverse positions in TCS8 coordinates where one inch of commanded movement will
360 !     yield one inch of movement on the traverse slides. However, this will not yield one inch of movement of the
370 !     probe volume crossover point within the water filled tunnel because of the differences of refraction in air,
380 !     glass, and water. Therefore, coordinate system transformations are applied to TCS8 positions to obtain
390 !     positions in "TUNNEL" and "MODEL" coordinates.
400 !     During data acquisition, real time histograms will be displayed of the LDV and analog data. After the data
410 !     has been acquired, the averages, standard deviations, and shear stresses will be calculated and displayed in
420 !     profile plots where the data is plotted versus traverse position. The reduced data is also sent to the
430 !     printer in tabular form. The reduced data as well as the raw data are stored along with the tunnel conditions
440 !     on the hard disc for archival purposes and also to allow for further data reduction, data plotting, or data
450 !     transfer to other computers.
460 !
470 ! PROGRAM OPERATION:
480 !
490 !     The following power up sequences should be completed before this program is run:
500 !     1. Turn on the "Motor Drive System" boxes.
510 !     2. Turn on the "TCS8" traverse control system.
520 !     3. Turn on the "LVDAS" Laser Velocimeter Data Acquisition System.
530 !     4. Turn on the HP series 9000 model 330 computer.
540 !     This program will automatically be loaded and run when the computer is turned on. If it is not loaded then
550 !     you can type in the following commands to load and then run it.
560 !     LOAD "LDVWT:,1400,0,0"
570 !     RUN
580 !     When the program is ready for user operation, it will display three things on the CRT. These are the main
590 !     menu, TCS8 traverse positions, and new sets of histogram & profile graphs. If they do not appear on the CRT
600 !     then you should perform the following actions to reinitialize the systems.
610 !     1. Press shift reset on the HP series 9000 model 330 computers keyboard.
620 !     2. Press reset on the back of the TCS8.
630 !     3. Press reset on the front (or back) of the LVDAS.
640 !     4. LOAD "LDVWT:,1400,0,0"
650 !     5. RUN
660 !
670 ! PROGRAM VARIABLES:
680 !
690 ! Mass Storage Variables:
700 !
710 !     System$    Tells the program where to read/store system data related files.
720 !     Data$      Tells the program where to read/store raw & reduced data related files.
730 !     File$      File name for tunnel conditions data or raw & reduced data.
740 !
750 ! Menu Variables:
760 !
770 !     Menu$(*)   String array where each element describes its corresponding menu subroutine's function.
780 !     Menu       Used as an index to the string array Menu$(*). Indicates which of the menus has been
790 !               selected as the current menu.
800 !     Key        Used as an index to the string array Menu$(*). Indicates which one of eight menu
810 !               subroutines in the menu is to be executed.
820 !     Busy       Tells the Menu Status subprogram to display the current menu selection in inverse video.
830 !     Ready      Tells the Menu Status subprogram to display the current menu selection in normal text.
840 !
850 ! Traverse Position Variables:
860 !
870 !     Tcs1(*)    TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TCS coordinates.
880 !     Tcs2(*)    TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TCS coordinates.
890 !     Tun(*)     Traverse positions (X,Y,Z) in TUNNEL coordinates.

```

```

900      !      Mod(*)      Traverse positions (X,Y,Z) in MODEL coordinates.
910      !
920      ! Auto Move Traverse Position Variables:
930      !
940      !      Pos(*)      Array of preprogrammed auto move positions.
950      !      Pname$(*)   Names for the variables in Pos(*).
960      !      Pimage$(*)  Image formats for the variables in Pos(*).
970      !      Punits$(*)  Units for the variables in Pos(*).
980      !      Npos       Number of preprogrammed auto move positions in Pos(*).
990      !      Paxis      Specifies which axis is to be traversed for the profile. Also defines axis for plots.
1000     !
1010     ! Traverse Coordinate System Transformation Variables:
1020     !
1030     !      Index(*)     Array of indexes of refraction for air, glass, and water.
1040     !                  Index(1) : Index of refraction for Air.
1050     !                  Index(2) : Index of refraction for Glass.
1060     !                  Index(3) : Index of refraction for Water.
1070     !      Theta      Tx Side Off Axis Angle.
1080     !      Fs         Focal length for sending side onaxis and offaxis lenses.
1090     !      Fr         Focal length for receiving side offaxis lens.
1100     !      Bs         Beam spacings for sending side onaxis and offaxis beam pairs.
1110     !      Br         Beam spacing for receiving side offaxis.
1120     !      Ts         Angle of offaxis sending side beam pair.
1130     !      Tr         Angle of offaxis receiving side beam pair.
1140     !      Ta         Sending side offaxis auxiliary traverse angle.
1150     !      Tcs2tun1(*) Sending side coordinate system transformation matrix for converting Tcs1(*) to Tun(*).
1160     !      Tun2tcs1(*) Sending side coordinate system transformation matrix for converting Tun(*) to Tcs1(*).
1170     !      Tcs2tun2(*) Receiving side coordinate system transformation matrix for converting Tcs2(*) to Tun(*).
1180     !      Tun2tcs2(*) Receiving side coordinate system transformation matrix for converting Tun(*) to Tcs2(*).
1190     !      Tun2mod(*)  Coordinate system transformation matrix for converting Tun(*) to Mod(*).
1200     !      Mod2tun(*)  Coordinate system transformation matrix for converting Mod(*) to Tun(*).
1210     !
1220     ! Velocity Coordinate System Transformation Variables:
1230     !
1240     !      Index(*)     Array of indexes of refraction for air, glass, and water.
1250     !                  Index(1) : Index of refraction for Air.
1260     !                  Index(2) : Index of refraction for Glass.
1270     !                  Index(3) : Index of refraction for Water.
1280     !      Theta1(*)   Angles between LASER & TUNNEL UVW laser beams in Air (N=Index1).
1290     !      Theta3(*)   Angles between LASER & TUNNEL UVW laser beams in Water (N=Index3).
1300     !      Tun2ldv(*)  Coordinate system transformation matrix for converting from TUNNEL to LASER.
1310     !      Ldv2tun(*)  Coordinate system transformation matrix for converting from LASER to TUNNEL.
1320     !
1330     ! Traverse & Velocity Coordinate System Transformation Variables:
1340     !
1350     !      Alpha(*)     Angles of attack, yaw, and roll.
1360     !                  Alpha(1) : Angle of Attack.
1370     !                  Alpha(2) : Angle of Yaw.
1380     !                  Alpha(3) : Angle of Roll.
1390     !      Mod2tun(*)  Coordinate system transformation matrix for converting positions & velocities from MODEL to TUNNEL.
1400     !      Tun2mod(*)  Coordinate system transformation matrix for converting positions & velocities from TUNNEL to MODEL.
1410     !
1420     ! Tunnel Condition Variables:
1430     !
1440     !      Array(*)     Array of tunnel conditions, laser parameters, graph scales, etc.
1450     !      Name$(*)    Names for the variables in Array(*).
1460     !      Image$(*)   Image formats for the variables in Array(*).
1470     !      Units$(*)   Units for the variables in Array(*).
1480     !
1490     ! Misc. Tunnel Condition Variables:
1500     !
1510     !      Date        Date.
1520     !      Time        Time.
1530     !      Run         Run Number.
1540     !      File        File Number.
1550     !      Mach        Mach Number.
1560     !      Temp        Room Temperature (deg. F).
1570     !      Uedge       Freestream Velocity.
1580     !      Ujet_ue     Jet exit velocity normalized by Uedge.
1590     !
1600     ! LVDAS Variables:
1610     !
1620     !      Table(*)     Lookup table of frequencies.
1630     !      Atime        The maximum desired acquisition time (seconds).
1640     !      Ctime        The maximum desired coincidence time (seconds).
1650     !      At_exp       Exponent for interarrival times.
1660     !      Ct_exp       Exponent for coincidence times.
1670     !      Nreads       Number of desired samples.
1680     !      Nsam         Number of acquired samples.
1690     !      Coin(*)     Coincidence criteria.

```

```

1700      !      Cmask      Coincidence mask for U,V,W selection.
1710      !      Raw(*)     Array of raw data acquired from the LVDAS.
1720      !
1730      ! Instantaneous Velocity and Voltage Variables:
1740      !
1750      !      Ui(*)       Read from LVDAS as the instantaneous U frequency data, then converted into U velocities.
1760      !      Vi(*)       Read from LVDAS as the instantaneous V frequency data, then converted into V velocities.
1770      !      Wi(*)       Read from LVDAS as the instantaneous W frequency data, then converted into W velocities.
1780      !      Ai(*)       Read from LVDAS as the instantaneous A voltage data.
1790      !      Bi(*)       Read from LVDAS as the instantaneous B voltage data.
1800      !      Ii(*)       Read from LVDAS as the raw interarrival time data, then converted into interarrival times.
1810      !      Ci(*)       Read from LVDAS as the raw coincidence time data, then converted into coincidence times.
1820      !      Valid(*)    Validation words. Initially all ones, then some set to zero during histogram clipping.
1830      !
1840      ! Histogram Clipping Variables:
1850      !
1860      !      Umin         The minimum acceptable U frequency (MHz).
1870      !      Umax         The maximum acceptable U frequency (MHz).
1880      !      Vmin         The minimum acceptable V frequency (MHz).
1890      !      Vmax         The maximum acceptable V frequency (MHz).
1900      !      Wmin         The minimum acceptable W frequency (MHz).
1910      !      Wmax         The maximum acceptable W frequency (MHz).
1920      !      Clip         Clip: 1 turn histogram clipping on; 0 turns it off.
1930      !
1940      ! Frequency to Velocity Conversion Variables:
1950      !
1960      !      Beam_spc(*)    Beam spacing at lens.
1970      !      Focl_len(*)   Focal length.
1980      !      Beam_sep(*)   Beam separation angle in degrees (full angle).
1990      !      Wave_len(*)   Wave length.
2000      !      Frng_spc(*)  Fringe Spacings.
2010      !      Brq_frq(*)   Bragg Frequencies.
2020      !      Mix_frq(*)   Mixing Frequencies.
2030      !      Mea_sgn(*)    Measured Frequencies' Signs.
2040      !      Brq_sgn(*)   Bragg Frequencies' Signs.
2050      !      Mix_sgn(*)   Mixing Frequencies' Signs.
2060      !
2070      ! Summation Variables:
2080      !
2090      !      Sum(1,1)       Summation of all of the valid Ui.
2100      !      Sum(2,1)       Summation of all of the valid Vi.
2110      !      Sum(3,1)       Summation of all of the valid Wi.
2120      !      Sum(4,1)       Summation of all of the valid Ai.
2130      !      Sum(5,1)       Summation of all of the valid Bi.
2140      !      Sum(6,1)       Summation of all of the valid Ii.
2150      !      Sum(7,1)       Summation of all of the valid Ci.
2160      !      Sum(1,2)       Summation of all of the valid Ui*Ui.
2170      !      Sum(2,2)       Summation of all of the valid Vi*Vi.
2180      !      Sum(3,2)       Summation of all of the valid Wi*Wi.
2190      !      Sum(4,2)       Summation of all of the valid Ai*Ai.
2200      !      Sum(5,2)       Summation of all of the valid Bi*Bi.
2210      !      Sum(6,2)       Summation of all of the valid Ii*Ii.
2220      !      Sum(7,2)       Summation of all of the valid Ci*Ci.
2230      !      Sum(1,3)       Summation of all of the valid Ui*Vi.
2240      !      Sum(2,3)       Summation of all of the valid Vi*Wi.
2250      !      Sum(3,3)       Summation of all of the valid Wi*Ui.
2260      !      Sum(4,3)       Summation of all of the valid Ai*Bi.
2270      !      Sum(5,3)       Summation of all of the valid Ui*Ai.
2280      !      Sum(6,3)       Summation of all of the valid Vi*Ai.
2290      !      Sum(7,3)       Summation of all of the valid Wi*Ai.
2300      !      N(*)          Number of valid samples for the above summations.
2310      !
2320      ! Reduced Data Variables:
2330      !
2340      !      U            Average U frequency or velocity.
2350      !      V            Average V frequency or velocity.
2360      !      W            Average W frequency or velocity.
2370      !      A            Average A voltage.
2380      !      B            Average B voltage.
2390      !      I            Average interarrival time.
2400      !      C            Average coincidence time.
2410      !      U1           Standard deviation for U frequency or velocity.
2420      !      V1           Standard deviation for V frequency or velocity.
2430      !      W1           Standard deviation for W frequency or velocity.
2440      !      A1           Standard deviation for A voltage.
2450      !      B1           Standard deviation for B voltage.
2460      !      I1           Standard deviation for interarrival time.
2470      !      C1           Standard deviation for coincidence time.
2480      !      Ulvl         Velocity:Velocity Shear Stress.
2490      !      Vlwl         Velocity:Velocity Shear Stress.

```

```

2500      !      Wlul      Velocity:Velocity Shear Stress.
2510      !      Albl      Voltage :Voltage Shear Stress.
2520      !      Ulal      Velocity:Voltage Shear Stress.
2530      !      Vlal      Velocity:Voltage Shear Stress.
2540      !      Wlal      Velocity:Voltage Shear Stress.
2550      !
2560      ! Data Plotting Symbol Variables:
2570      !
2580      !      Symbols(*)   Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
2590      !      Symbol(*)   Array of coordinates which when connected produce a distinct geometric symbol.
2600      !      Dot(*)      Array of coordinates which produce a dot. The dot symbol is added to all symbols.
2610      !      Noc          The number of coordinates in a symbol.
2620      !      S            Used to index the Symbols array.
2630      !
2640      ! Histogram and Profile Graph Variables:
2650      !
2660      !      Wndw(*)       Array containing the plot's scales.
2670      !      Vwprt(*)     Array containing the plot's CRT position.
2680      !      Xdiv(*)      Array containing the number of X divisions for the plot's X axis.
2690      !      Ydiv(*)      Array containing the number of Y divisions for the plot's Y axis.
2700      !      Xlabel$(*)  String array containing labels for the X axis.
2710      !      Ylabel$(*)  String array containing labels for the Y axis.
2720      !      Title$(*)   String array containing labels for the Plots.
2730      !      Ximage$(*)   String array containing image formats for the X axis labeling.
2740      !      Yimage$(*)   String array containing image formats for the Y axis labeling.
2750      !      Legend$(*)   String array containing labels for each symbol in a profile plot.
2760      !      G            Used as an index to the above arrays. Specifies one of nine plots.
2770      !      Gsave(*)     Used to save the entire graphics contents of the CRT.
2780      !
2790      !
2800      ! OPTION BASE 1
2810      ! COM /Data/ INTEGER Raw(1000,10),Valid(1000),REAL Table(0:32766),U1(1000),V1(1000),W1(1000),A1(1000),
2820      !      Bi(1000),Ii(1000),Ci(1000)
2830      ! COM /Array/ Name$(100,4)[10],Image$(100,4)[10],Units$(100,4)[10],REAL Array(100,4)
2840      ! COM /Pos/ Pname$(25,1)[10],Pimage$(25,1)[10],Punits$(25,1)[10],REAL Pos(25,1),Npos
2850      ! COM /Graph/ Wndw(9,4),Vwprt(9,4),Xdiv(9),Ydiv(9),Xlabel$(9)[80],Ylabel$(9)[80],Title$(9)[80],
2860      !      Ximage$(9)[80],Yimage$(9)[80],Legend$(9,5)[80]
2870      ! COM Run,File,Paxis
2880      ! DIM Menu$(5,8)[80],System$(20),Data$(20),File$(50),L$(160)
2890      ! INTEGER Gsave(1280,1024),At_exp,Ct_exp,Cmask,Nsam,N(10,3)
2900      ! REAL Atime,Ctime,Sum(10,3),Symbols(5,0:20,3),Apos,Bpos
2910      ! DIM Tcs2tun1(4,4),Tun2tcs1(4,4),Tun2mod(3,3),Tun2ldv(3,3),Tun(4),Tcs1(4)
2920      ! DIM Tcs2tun2(4,4),Tun2tcs2(4,4),Mod2tun(3,3),Ldv2tun(3,3),Mod(4),Tcs2(4)
2930      ! DIM Beam_spc(3),Focl_len(3),Mea_sgn(3),Mix_frq(3),Mix_sgn(3),Frng_spc(3),Thetal(3,3)
2940      ! DIM Beam_sep(3),Wave_len(3),Brg_frq(3),Brg_sgn(3),Index(3),Coin(3),Theta3(3,3),Alpha(3)
2950      ! Perform trigonometric operations in degrees.
2960      ! DEG
2970      ! Clear the CRT and direct printed output to it.
2980      ! CLEAR SCREEN
2990      ! GCLEAR
3000      ! PRINTER IS CRT
3010      ! Perform any necessary setup and initialization routines.
3020      ! GOSUB Lvds_set_up      ! Initialize the HP to LVDAS interface.
3030      ! GOSUB File_set_up    ! Select mass storage device for system & data files.
3040      ! GOSUB Tcs8_set_up    ! Initialize the HP to TCS8 interface.
3050      ! GOSUB Menu_set_up    ! Initialize the user driven menus and display the main menu.
3060      ! GOSUB Grph_set_up    ! Initialize the CRT and plot the nine empty plots for profiles and histograms.
3070      ! Here:               ! The main program, while continually displaying the time of day, will wait hear for menu key selection.
3080      ! Date=TIMEDATE
3090      ! Time=Date
3100      ! DISP TIME$(Time),DATE$(Date)
3110      ! GOTO Here
3120      ! On_key:
3130      ! ON KEY 1 GOSUB Key1    ! If the user function key #1 is ever pressed then execute the "Key1" subroutine.
3140      ! ON KEY 2 GOSUB Key2    ! If the user function key #2 is ever pressed then execute the "Key2" subroutine.
3150      ! ON KEY 3 GOSUB Key3    ! If the user function key #3 is ever pressed then execute the "Key3" subroutine.
3160      ! ON KEY 4 GOSUB Key4    ! If the user function key #4 is ever pressed then execute the "Key4" subroutine.
3170      ! ON KEY 5 GOSUB Key5    ! If the user function key #5 is ever pressed then execute the "Key5" subroutine.
3180      ! ON KEY 6 GOSUB Key6    ! If the user function key #6 is ever pressed then execute the "Key6" subroutine.
3190      ! ON KEY 7 GOSUB Key7    ! If the user function key #7 is ever pressed then execute the "Key7" subroutine.
3200      ! ON KEY 8 GOSUB Key8    ! If the user function key #8 is ever pressed then execute the "Key8" subroutine.
3210      ! RETURN
3220      ! Subroutine Key1,Key2,Key3,Key4,Key5,Key6,Key7,Key8 descriptions:
3230      ! When one of the special user function keys is pressed, the main program will execute one the
3240      ! following eight subroutines. Each of these subroutines performs essentially the same basic
3250      ! function in that it subsequently executes one of the menu subroutines. The particular menu
3260      ! subroutine to be executed will depend on the current menu selected and the current key pressed.
3270      ! Before the selected menu subroutine is executed, the corresponding menu entry at the top of
3280      ! the CRT is redisplayed in inverse video. This indicates that the menu selection has been
3290      ! acknowledged and that any resultant actions are still in progress. When the highlighted menu
3300      ! subroutine has completed the current TCS8 traverse positions will be read and updated on the CRT
3310      ! display. The corresponding menu entry displayed at the top of the CRT is redisplayed in normal

```

```

3280      !      text to indicate the completion of the menu subroutine.  The user can then select another special
3290      !      function key.
3300      !      Variables:
3310      !      Menu      Indicates which of the menus has been selected as the current menu.
3320      !      Key       Indicates which one of eight menu subroutines in the menu is to be executed.
3330      !      Menu$(*)  String array where each element describes its corresponding menu subroutine's function.
3340      !      Busy      Tells the Menu Status subroutine to display the current menu selection in inverse video.
3350      !      Ready     Tells the Menu Status subroutine to display the current menu selection in normal text.
3360 Key1:  Key=1
3370      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3380      ON Menu GOSUB M1k1,M2k1,M3k1,M4k1,M5k1,M6k1,M7k1
3390      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3400      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3410      RETURN
3420 Key2:  Key=2
3430      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3440      ON Menu GOSUB M1k2,M2k2,M3k2,M4k2,M5k2,M6k2,M7k2
3450      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3460      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3470      RETURN
3480 Key3:  Key=3
3490      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3500      ON Menu GOSUB M1k3,M2k3,M3k3,M4k3,M5k3,M6k3,M7k3
3510      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3520      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3530      RETURN
3540 Key4:  Key=4
3550      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3560      ON Menu GOSUB M1k4,M2k4,M3k4,M4k4,M5k4,M6k4,M7k4
3570      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3580      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3590      RETURN
3600 Key5:  Key=5
3610      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3620      ON Menu GOSUB M1k5,M2k5,M3k5,M4k5,M5k5,M6k5,M7k5
3630      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3640      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3650      RETURN
3660 Key6:  Key=6
3670      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3680      ON Menu GOSUB M1k6,M2k6,M3k6,M4k6,M5k6,M6k6,M7k6
3690      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3700      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3710      RETURN
3720 Key7:  Key=7
3730      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3740      ON Menu GOSUB M1k7,M2k7,M3k7,M4k7,M5k7,M6k7,M7k7
3750      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3760      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3770      RETURN
3780 Key8:  Key=8
3790      CALL Menu_status(Menu,Key,Busy,Menu$(*))
3800      ON Menu GOSUB M1k8,M2k8,M3k8,M4k8,M5k8,M6k8,M7k8
3810      CALL Menu_status(Menu,Key,Ready,Menu$(*))
3820      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
3830      RETURN
3840 Menu1:  ! Descriptions of the "Main Menu" subroutines M1K1,...,M1K8:
3850      !      The eight subroutines M1K1,...,M1K8 together implement the "Main Menu".  The following will be
3860      !      displayed at the top left of the CRT display when the "Main Menu" is selected:
3870      !
3880      !      M1K1: Laser Alignment
3890      !      M1K2: Pre Run
3900      !      M1K3: Post Run (Dump Graphics)
3910      !      M1K4: Set Auto Move Positions
3920      !      M1K5: Move traverse
3930      !      M1K6: Take data
3940      !      M1K7: Auto move and take
3950      !      M1K8: Display Histograms
3960      !
3970      !      M1K1 will change the current active menu from the "Main Menu" to the "Laser Alignment Menu".
3980      !      M1K2 will change the current active menu from the "Main Menu" to the "Pre Run Menu".  M1K3 will
3990      !      transfer the graphics contents of the CRT to the printer.  This provides a hard copy of the profile
4000      !      plots.  M1K4 has the user enter predefined traverse positions for a profile plot.  M1K5 moves the
4010      !      traverse to a user selectable position.  M1K6 acquires LVDAS data at the current TCS8 traverse
4020      !      position.  M1K7 acquires LVDAS data at each of the pre programed TCS8 traverse positions set up by
4030      !      M1K4.  M1K8 repeatedly displays five channels of real time histograms until the user presses any
4040      !      key on the keyboard.
4050      !
4060 M1k1:  ! Change the current active menu from the "Main Menu" to the "Laser Alignment Menu".
4070      Menu=2

```



```

4080      CALL Menu_disp(Menu,Menu$(*))
4090      RETURN
4100 M1k2:  ! Change the current active menu from the "Main Menu" to the "Pre Run Menu".
4110      Menu=3
4120      CALL Menu_disp(Menu,Menu$(*))
4130      RETURN
4140 M1k3:  ! Transfer the graphics contents of the CRT to the printer. This provides a hard copy of the plots.
4150      KEY LABELS OFF          ! Turn off the key labels so that they won't be printed.
4160      PRINTER IS CRT;WIDTH 132
4170      DISP ""                ! Clear the CRT's display line so that they won't be printed.
4180      FOR L=1 TO 9            ! Clear the CRT's menu lines so that it won't be printed.
4190          PRINT TABXY(1,L);RPTS(" ",120)
4200      NEXT L
4210      PRINTER IS PRT
4220      PRINT USING "#,@"      ! Move to the top of the next page on the printer.
4230      DUMP GRAPHICS          ! Dump the entire CRT's contents to the printer.
4240      PRINT USING "#,@"      ! Move to the top of the next page on the printer.
4250      PRINTER IS CRT
4260      CALL Menu_disp(Menu,Menu$(*)) ! Redisplay the menus.
4270      RETURN
4280 M1k4:  ! Have the user enter predefined traverse positions for a profile plot.
4290      CALL Enter_value("number of traverse positions",Npos,"K")
4300      REDIM Pos(Npos,1),Pname$(Npos,1),Pimage$(Npos,1),Punits$(Npos,1)
4310      MAT Pimage$= ("M4D.4D")
4320      MAT Punits$= ("in")
4330      FOR K=1 TO Npos
4340          Pname$(K,1)="Pos"&VAL$(K)
4350      NEXT K
4360      GSTORE Gsave(*)
4370      CALL Change("VALUES",Pos(*),Pname$(*),Pimage$(*),Punits$(*))
4380      GLOAD Gsave(*)
4390      CALL Menu_disp(Menu,Menu$(*))
4400      RETURN
4410 M1k5:  ! Moves the traverse to a user selectable position.
4420      GOSUB Read_calc_fill
4430      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
4440      CALL Enter_value(CHRS(NUM("X")+Paxis-1),Mod(Paxis),"K")
4450 M1k5a: ON KBD CALL Do_nothing
4460      DISP "Moving"
4470      Movement=Mod(Paxis)
4480      CALL Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),"Tx & Rx","MODEL",
          "ABSOLUTE",Paxis,Movement)
4490      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
4500      GOSUB Calc
4510      GOSUB Fill
4520      DISP ""
4530      OFF KBD
4540      RETURN
4550 M1k6:  ! Acquire LVDAS data at the current TCS8 traverse position.
4560      ! DISP "Press any key to TAKE DATA"
4570      ! CALL Rt_histo(@Lvdas,Symbols(*),1)
4580      Cmask=Coin(1)*1+Coin(2)*2+Coin(3)*4
4590      Nsam=MIN(Nreads,1000)
4600      Date=TIMEDATE
4610      Time=Date
4620      CALL Lvdas_take(@Lvdas,Atime,Ctime,At_exp,Ct_exp,Cmask,Nsam)
4630      IF Nsam>1 THEN
4640          SELECT Paxis          ! Select the non-scanning axes for printing their position values at each point.
4650          CASE 3
4660              D$="X"
4670              F$="Y"
4680              Apos=Mod(1)
4690              Bpos=Mod(2)
4700          CASE 2
4710              D$="X"
4720              F$="Z"
4730              Apos=Mod(1)
4740              Bpos=Mod(3)
4750          CASE 1
4760              D$="Y"
4770              F$="Z"
4780              Apos=Mod(2)
4790              Bpos=Mod(3)
4800          END SELECT
4810          OUTPUT PRT USING "K,K";CHRS(27)&"&k2S"&CHRS(27)&"&19D",RPTS("=",140)
4820          PRINTER IS PRT
4830          Run$=VAL$(Run)
4840          File$=VAL$(File)
4850          PRINT " RUN "&Run$&" FILE "&File$
4860          A$=DATE$(TIMEDATE) ! Acquire the date and time for printing at each point.

```

```

4870 BS=TIMES(TIMEDATE)
4880 PRINT USING 4890;AS,BS,DS,Apos,FS,Bpos
4890 IMAGE 6X,11A,3X,8A,3X,A,"=",3D.4D,3X,A,"=",3D.4D
4900 PRINTER IS CRT
4910 CALL Data_reduce(At_exp,Ct_exp,Nsam)
4920 CALL Pt_histo(Symbols(*),Run,File,Mod(Paxis),Nsam)
4930 CALL Data_clip(Nsam,Umin,Umax,Vmin,Vmax,Wmin,Wmax)
4940 CALL Data_sum(Sum(*),N(*),Nsam)
4950 CALL Data_calc(N(*),Sum(*),U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4960 B=U/U ! Replacement of B, B1, and Alb1 is being made by velocity ratios since
4970 B1=V/U ! the second analog channel B is not being used
4980 Alb1=W/U
4990 CALL Data_print(Paxis,Mod(Paxis),Nsam,"MHZ",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,
    U1a1,V1a1,W1a1,Uedge)
5000 CALL Data_fconvert(Array(*))
5010 CALL Data_sum(Sum(*),N(*),Nsam)
5020 CALL Data_calc(N(*),Sum(*),U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
5030 B=U/Uedge
5040 B1=V/Uedge
5050 Alb1=W/Uedge
5060 CALL Data_print(Paxis,Mod(Paxis),Nsam,"LDV",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,
    U1a1,V1a1,W1a1,Uedge)
5070 CALL Data_trnsfrm(Ldv2tun(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1,U1a1,V1a1,W1a1)
5080 B=U/Uedge
5090 B1=V/Uedge
5100 Alb1=W/Uedge
5110 CALL Data_print(Paxis,Mod(Paxis),Nsam,"TUN",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,
    U1a1,V1a1,W1a1,Uedge)
5120 CALL Data_trnsfrm(Tun2mod(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1,U1a1,V1a1,W1a1)
5130 B=U/Uedge
5140 B1=V/Uedge
5150 Alb1=W/Uedge
5160 CALL Data_print(Paxis,Mod(Paxis),Nsam,"MOD",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,
    U1a1,V1a1,W1a1,Uedge)
5170 CALL Data_plot(Array(*),Symbols(*),6,Mod(Paxis),U,V,W,1/Uedge,N(1,1),N(2,1),N(3,1))
5180 CALL Data_plot(Array(*),Symbols(*),7,Mod(Paxis),U1,V1,W1,1/Uedge,N(1,2),N(2,2),N(3,2))
5190 CALL Data_plot(Array(*),Symbols(*),8,Mod(Paxis),U1v1,V1w1,W1u1,1/Uedge^2,N(1,3),N(2,3),N(3,3))
5200 CALL Data_plot(Array(*),Symbols(*),9,Mod(Paxis),A,A,A1,1,N(4,1),N(4,1),N(4,2))
5210 OUTPUT PRT USING "K,K";CHR$(27)&"&k2S"&CHR$(27)&"&l9D",RPT$("=",140)
5220 GOSUB Store_file
5230 File=File+1
5240 END IF
5250 RETURN
5260 M1k7: ! Acquire LVDAS data at each of the pre programed TCS8 traverse positions set up by M1k4.
5270 Quit=0
5280 ON KBD GOSUB Quit
5290 FOR J=1 TO Npos
5300 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
5310 Mod(Paxis)=Pos(J,1)
5320 GOSUB M1k5a
5330 GOSUB M1k6
5340 IF Quit THEN 5360
5350 NEXT J
5360 OFF KBD
5370 GOSUB On_key
5380 CALL Menu_disp(Menu,Menu$(*))
5390 RETURN
5400 M1k8: ! Repeatedly displays five channels of real time histograms until the user presses any key on the keyboard.
5410 DISP "Press any key to return to main menu"
5420 CALL Rt_histo(@Lvdas,Symbols(*),1)
5430 RETURN
5440 Menu2: ! Descriptions of the "Laser Alignment Menu" subroutines M2K1,...,M2K8:
5450 ! The eight subroutines M2K1,...,M2K8 together implement the "Laser Alignment Menu". The
5460 ! following will be displayed at the top left of the CRT display when the "Laser Alignment Menu" is
5470 ! selected:
5480 !
5490 ! M2K1: Return to main menu
5500 ! M2K2: Sides : Tx & Rx
5510 ! M2K3: Coordinates: MODEL
5520 ! M2K4: Mode : ABSOLUTE
5530 ! M2K5: Move X
5540 ! M2K6: Move Y
5550 ! M2K7: Move Z
5560 ! M2K8: Move A
5570 !
5580 ! M2K1 will change the current active menu from the "Laser Alignment Menu" to the "Main Menu".
5590 ! M2K2 selects whether the transmitting, receiving, or both sides of the traverse are to be moved.
5600 ! M2K3 selects the TCS, TUNNEL, or MODEL coordinate systems for traverse movements. M2K4
5610 ! specifies movements to be relative to the current position or to absolute positions. M2K5 has the
5620 ! user enter a movement for the X axis and then the movement is performed. M2K6 has the user enter

```

```

5630      !      a movement for the Y axis and then the movement is performed. M2K7 has the user enter a movement
5640      !      for the Z axis and then the movement is performed. M2K8 has the user enter a movement for the A
5650      !      axis and then the movement is performed.
5660      !
5670 M2k1:  ! Change the current active menu from the "Laser Alignment Menu" to the "Main Menu".
5680      Menu=1
5690      CALL Menu_disp(Menu,Menu$(*))
5700      RETURN
5710 M2k2:  ! Select whether the transmitting, receiving, or both sides of the traverse are to be moved.
5720      SELECT TRIMS(Menu$(Menu,Key)[20])
5730      CASE "Tx & Rx"
5740          Menu$(Menu,Key)[20]="Tx"
5750      CASE "Tx"
5760          Menu$(Menu,Key)[20]="Rx"
5770      CASE "Rx"
5780          Menu$(Menu,Key)[20]="Tx & Rx"
5790      END SELECT
5800      CALL Menu_disp(Menu,Menu$(*))
5810      RETURN
5820 M2k3:  ! Selects the TCS, TUNNEL, or MODEL coordinate systems for traverse movements.
5830      SELECT TRIMS(Menu$(Menu,Key)[20])
5840      CASE "MODEL"
5850          Menu$(Menu,Key)[20]="TUNNEL"
5860      CASE "TUNNEL"
5870          Menu$(Menu,Key)[20]="TCS"
5880      CASE "TCS"
5890          Menu$(Menu,Key)[20]="MODEL"
5900      END SELECT
5910      CALL Menu_disp(Menu,Menu$(*))
5920      RETURN
5930 M2k4:  ! Specifies movements to be relative to the current's position or to absolute positions.
5940      SELECT TRIMS(Menu$(Menu,Key)[20])
5950      CASE "ABSOLUTE"
5960          Menu$(Menu,Key)[20]="RELATIVE"
5970      CASE "RELATIVE"
5980          Menu$(Menu,Key)[20]="ABSOLUTE"
5990      END SELECT
6000      CALL Menu_disp(Menu,Menu$(*))
6010      RETURN
6020 M2k5:  !      The subroutines M2K5 thru M2K8 all execute the same code. The code will have the user enter a
6030 M2k6:  !      movement for the X,Y,Z, or A depending on what the value of "Key" is. The user specified movement
6040 M2k7:  !      for the selected axis will then be performed.
6050 M2k8:  !
6060      Side$=TRIMS(Menu$(Menu,2)[20])
6070      Coord$=TRIMS(Menu$(Menu,3)[20])
6080      Mode$=TRIMS(Menu$(Menu,4)[20])
6090      CALL Enter_value(Mode$&" Movement",Movement,"4D.5D")
6100      ON KBD CALL Do_nothing
6110      DISP "Moving"
6120      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
6130      CALL Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),Side$,Coord$,Mode$,
        Key-4,Movement)
6140      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
6150      DISP ""
6160      OFF KBD
6170      RETURN
6180 Menu3:  ! Descriptions of the "Pre Run Menu" subroutines M3K1,...,M3K8:
6190      !      The eight subroutines M3K1,...,M3K8 together implement the "Pre Run Menu". The following will
6200      !      be displayed at the top left of the CRT display when the "Pre Run Menu" is selected:
6210      !
6220      !      M3K1: Return to MAIN menu
6230      !      M3K2: Enter Run & File Numbers
6240      !      M3K3: Enter Number of Samples
6250      !      M3K4: Select Traverse Axis for Profile
6260      !      M3K5: Print Coordinate Transformation Matrices
6270      !      M3K6: Setup Graphics
6280      !      M3K7: Tunnel Conditions
6290      !      M3K8: Traverse
6300      !
6310      !      M3K1 will change the current active menu from the "Pre Run Menu" to the "Main Menu". M3K2 has
6320      !      the user enter a the Run and File numbers. A new run number should be assigned to each profile
6330      !      while a new file number is assigned to each set of data. M3K3 has the user enter the desired
6340      !      number of samples. M3K4 has the user select which axis to traverse in for the profiles. M3K5
6350      !      prints the coordinate system transformation matrices for both traverse positions and velocities.
6360      !      M3K6 creates a new set of empty plots for new profiles. M3K7 will change the current active menu
6370      !      from the "Pre Run Menu" to the "Tunnel Conditions Menu". M3K8 will change the current active menu
6380      !      from the "Pre Run Menu" to the "Traverse Menu".
6390      !
6400 M3k1:  ! Change the current active menu from the "Pre Run Menu" to the "Main Menu".
6410      Menu=1

```

```

6420      CALL Menu_disp(Menu,Menu$(*))
6430      RETURN
6440 M3k2:  ! Have the user enter a the Run and File numbers.
6450      CALL Enter_value("Run",Run,"3D.2D")
6460      CALL Enter_value("File",File,"3D")
6470      RETURN
6480 M3k3:  ! Have the user enter the desired number of samples.
6490      CALL Enter_value("Number of Samples ",Nreads,"K")
6500      RETURN
6510 M3k4:  ! Have the user select which axis to traverse in for the profiles.
6520      CALL Enter_string("Traverse Axis for Profiles ",Paxis$,"K")
6530      SELECT Paxis$
6540      CASE "X"
6550          Paxis=1
6560      CASE "Y"
6570          Paxis=2
6580      CASE "Z"
6590          Paxis=3
6600      CASE "A"
6610          Paxis=4
6620      CASE ELSE
6630          GOTO M3k4
6640      END SELECT
6650      GOSUB Fill
6660      RETURN
6670 M3k5:  ! Prints the coordinate system transformation matrices for both traverse positions and velocities.
6680      GOSUB Read_calc_fill
6690      OUTPUT PRT USING "#,2/"
6700      OUTPUT PRT USING "20X,K,/,/";"TRAVERSE COORDINATE TRANSFORMATION MATRICES"
6710      OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/);"Transmitting side TCS8 to TUNNEL",Tcs2tun1(*)
6720      OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/);"Receiving side TCS8 to TUNNEL",Tcs2tun2(*)
6730      OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/);"Transmitting side TUNNEL to TCS8",Tun2tcs1(*)
6740      OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/);"Receiving side TUNNEL to TCS8",Tun2tcs2(*)
6750      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"TUNNEL to MODEL",Tun2mod(*)
6760      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"MODEL to TUNNEL",Mod2tun(*)
6770      OUTPUT PRT USING "20X,K,/,/";"VELOCITY COORDINATE TRANSFORMATION MATRICES"
6780      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"LASER to TUNNEL",Ldv2tun(*)
6790      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"TUNNEL to LASER",Tun2ldv(*)
6800      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"TUNNEL to MODEL",Tun2mod(*)
6810      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"MODEL to TUNNEL",Mod2tun(*)
6820      OUTPUT PRT USING "#,0"
6830      RETURN
6840 M3k6:  ! Display a new set of plots for new profiles.
6850      CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
6860      RETURN
6870 M3k7:  ! Change the current active menu from the "Pre Run Menu" to the "Tunnel Conditions Menu".
6880      Menu=4
6890      CALL Menu_disp(Menu,Menu$(*))
6900      RETURN
6910 M3k8:  ! Change the current active menu from the "Pre Run Menu" to the "Traverse Menu".
6920      Menu=5
6930      CALL Menu_disp(Menu,Menu$(*))
6940      RETURN
6950 Menu4:  ! Descriptions of the "Tunnel Conditions Menu" subroutines M4K1,...,M4K8:
6960      !
6970      !     The eight subroutines M4K1,...,M4K8 together implement the "Tunnel Conditions Menu". The
6980      !     following will be displayed at the top left of the CRT display when the "Tunnel Conditions Menu" is
6990      !     selected:
7000      !
7010      !           M4K1: Return to PRE RUN menu
7020      !           M4K2: Load Tunnel Conditions
7030      !           M4K3: Save Tunnel Conditions
7040      !           M4K4: Print Tunnel Conditions
7050      !           M4K5: Enter Tunnel Condition Data
7060      !           M4K6: Enter Tunnel Condition Names
7070      !           M4K7: Enter Tunnel Condition Units
7080      !           M4K8: Enter Tunnel Condition Images
7090      !
7100      !           M4K1 will change the current active menu from the "Tunnel Conditions Menu" to the "Pre Run
7110      !           Menu". M4K2 loads the old tunnel conditions from a file on the disk. M4K3 saves the current
7120      !           tunnel conditions to a file on the disk. M4K2 & M4K3 load and save default tunnel conditions from
7130      !           the file "ARRAY" on the hard disk. The default values are not related to any particular run number.
7140      !           M4K4 sends the current tunnel conditions to the printer. M4K5 has the user enter values for the
7150      !           tunnel condition variables. M4K6 has the user enter names for the tunnel condition variables.
7160      !           M4K7 has the user enter units for the tunnel condition variables. M4K8 has the user enter image
7170      !           formats for the tunnel condition variables.
7180 M4k1:  ! Change the current active menu from the "Tunnel Conditions Menu" to the "Pre Run Menu".
7190      Menu=3
7200      CALL Menu_disp(Menu,Menu$(*))
7210      RETURN

```

```

7220 M4k2:      ! Load the old tunnel conditions from a file on the disk. This loads the default values.
7230      GOSUB Read_array
7240      GOSUB Read_calc_fill
7250      RETURN
7260 M4k3:      ! Save the current tunnel conditions to a file on the disk. This updates the default values on the disk.
7270      GOSUB Read_calc_fill
7280      GOSUB Save_array
7290      RETURN
7300 M4k4:      ! Print the current tunnel conditions.
7310      GOSUB Read_calc_fill
7320      GOSUB Print_header
7330      RETURN
7340 M4k5:      ! Have the user enter values for the tunnel condition variables.
7350      GSTORE Gsave(*)
7360      GOSUB Read_calc_fill
7370      CALL Change("VALUES",Array(*),Name$(*),Image$(*),Units$(*))
7380      GOSUB Read_calc_fill
7390      GLOAD Gsave(*)
7400      RETURN
7410 M4k6:      ! Have the user enter names for the tunnel condition variables.
7420      GSTORE Gsave(*)
7430      GOSUB Read_calc_fill
7440      CALL Change("NAMES",Array(*),Name$(*),Image$(*),Units$(*))
7450      GOSUB Read_calc_fill
7460      GLOAD Gsave(*)
7470      RETURN
7480 M4k7:      ! Have the user enter units for the tunnel condition variables.
7490      GSTORE Gsave(*)
7500      GOSUB Read_calc_fill
7510      CALL Change("UNITS",Array(*),Name$(*),Image$(*),Units$(*))
7520      GOSUB Read_calc_fill
7530      GLOAD Gsave(*)
7540      RETURN
7550 M4k8:      ! Have the user enter image formats for the tunnel condition variables.
7560      GSTORE Gsave(*)
7570      GOSUB Read_calc_fill
7580      CALL Change("IMAGES",Array(*),Name$(*),Image$(*),Units$(*))
7590      GOSUB Read_calc_fill
7600      GLOAD Gsave(*)
7610      RETURN
7620 Menu5:      ! Descriptions of the "Traverse Menu" subroutines M5K1,...,M5K8:
7630      !           The eight subroutines M5K1,...,M5K8 together implement the "Traverse Menu". The following will
7640      !           be displayed at the top left of the CRT display when the "Traverse Menu" is selected:
7650      !
7660      !           M5K1: Return to PRE RUN menu
7670      !           M5K2: View & Set TCS8 Positions
7680      !           M5K3: View & Set TCS8 Units
7690      !           M5K4: View & Set TCS8 Revolution
7700      !           M5K5: View & Set TCS8 Velocity
7710      !           M5K6: View & Set TCS8 Acceleration
7720      !           M5K7:
7730      !           M5K8:
7740      !
7750      !           M5K1 will change the current active menu from the "Traverse Menu" to the "Pre Run Menu". M5K2
7760      !           reads from the TCS8 the current positions and lets the user change them. The new positions are
7770      !           then sent to the TCS8. M5K3 reads from the TCS8 the current counts per unit length (inches) and
7780      !           lets the user change them. The new counts per unit length are then sent to the TCS8. M5K4 reads
7790      !           from the TCS8 the current counts per revolution and lets the user change them. The new counts per
7800      !           revolution are then sent to the TCS8. M5K5 reads from the TCS8 the current velocities and lets the
7810      !           user change them. The new velocities are then sent to the TCS8. M5K6 reads from the TCS8 the
7820      !           current accelerations and lets the user change them. The new accelerations are then sent to the
7830      !           TCS8. M5K7 does nothing. M5K8 does nothing.
7840      !
7850 M5k1:      ! Change the current active menu from the "Traverse Menu" to the "Pre Run Menu".
7860      Menu=3
7870      CALL Menu_disp(Menu,Menu$(*))
7880      RETURN
7890 M5k2:      ! Read current TCS8 positions, have the user update them, & then send the new values to the TCS8.
7900      CALL Tcs8set("P",@Tcs8)      ! View and set TCS8 Positions.
7910      GRAPHICS ON
7920      CALL Menu_disp(Menu,Menu$(*))
7930      RETURN
7940 M5k3:      ! Read current TCS8 counts per inch, have the user update them, & then send the new values to the TCS8.
7950      CALL Tcs8set("U",@Tcs8)      ! View and set TCS8 counts per Unit length.
7960      GRAPHICS ON
7970      CALL Menu_disp(Menu,Menu$(*))
7980      RETURN
7990 M5k4:      ! Read current TCS8 counts per revolution, have the user update them, & then send new values to the TCS8.
8000      CALL Tcs8set("R",@Tcs8)      ! View and set TCS8 counts per Revolution.
8010      GRAPHICS ON

```

```

8020      CALL Menu_disp(Menu,Menu$(*))
8030      RETURN
8040 M5k5:  ! Read current TCS8 velocities, have the user update them, & then send the new values to the TCS8.
8050      CALL Tcs8set("V",@Tcs8)          ! View and set TCS8 Velocities.
8060      GRAPHICS ON
8070      CALL Menu_disp(Menu,Menu$(*))
8080      RETURN
8090 M5k6:  ! Read current TCS8 accelerations, have the user update them, & then send the new values to the TCS8.
8100      CALL Tcs8set("A",@Tcs8)          ! View and set TCS8 Accelerations.
8110      GRAPHICS ON
8120      CALL Menu_disp(Menu,Menu$(*))
8130      RETURN
8140 M5k7:  ! This subroutine does nothing.
8150 M5k8:  ! This subroutine does nothing.
8160 Quit:  Quit=1 ! Quit will be set during a multiple position scan (see M1K7) if any key is pressed on the
8170      RETURN ! keyboard during the scan. This indicates that the scan should be terminated.
8180 Lvds_set_up: ! This subroutine initializes the HP to LVDAS high speed parallel interface. A communications path named
8190      ! "@Lvds" is opened. Also, this subroutine creates the raw data to frequency conversion look up table.
8200      CALL Lvdsas_init(@Lvdsas)
8210      CALL Table(Table(*))
8220      RETURN
8230 File_set_up: ! This subroutine reads the initialization files from the disk. System$ tells the program where to read
8240      ! system related files while Data$ tells the program where to read and store raw and reduced data.
8250      System$=":",1400,0,0"
8260      Data$=":",1400,0,1"
8270      LOAD KEY "KEYS"&System$
8280      GOSUB Read_array
8290      GOSUB Read_calc_fill
8300      GOSUB Save_array
8310      CLEAR SCREEN
8320      RETURN
8330 Tcs8_set_up: ! This subroutine initialized the HP to TCS8 serial interface. The communications path "@Tcs8" is opened.
8340      CALL Tcs8init(@Tcs8)
8350      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
8360      GOSUB Calc
8370      GOSUB Fill
8380      RETURN
8390 Grph_set_up: ! This subroutine defines the graphics symbols for plotting data points, clears and initializes the CRT,
8400      ! and displays a new empty set of graphs for histogram and profile plotting.
8410      CALL Read_symbols(Symbols(*))
8420      CALL Crt_init
8430      CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
8440      RETURN
8450 Menu_set_up: ! This subroutine defines the menu descriptors for all of the menus. The current menu is set to the "Main
8460      ! Menu" and its menu is displayed at the top of the screen.
8470      CALL Menu_read(Menu$(*))
8480      CALL Menu_disp(Menu,Menu$(*))
8490      GOSUB On_key
8500      Busy=0
8510      Ready=1
8520      RETURN
8530 Print_header: ! This subroutine prints a header on the printers paper. The "header" is a formatted list of all of the
8540      ! tunnel conditions, laser parameters, and graph scales.
8550      PRINTER IS PRT;WIDTH 144
8560      PRINT USING "#,5(K)";CHR$(27)&"&k2S"&CHR$(27)&"&l9D"
8570      CALL Array_print(Array(*),Name$(*),Image$(*),Units$(*))
8580      PRINT USING "#,0,5(K)";CHR$(27)&"E"
8590      PRINTER IS CRT
8600      RETURN
8610 Read_calc_fill: ! This subroutine extracts (reads) the tunnel conditions from the Array(*). These values can be used to
8620      ! calculate other tunnel conditions. The original tunnel conditions along with any calculated tunnel
8630      ! conditions are then put back (filled) into the Array(*).
8640      GOSUB Read
8650      GOSUB Calc
8660      GOSUB Fill
8670      RETURN
8680 Store_header: ! This subroutine stores the header Array(*) and other arrays onto the disk. There will be one header
8690      ! file for each run number. For example, if the run number equal 1, then the data will be stored in a
8700      ! file named "R1". This file will include an extensive list of tunnel conditions, laser parameters, graph
8710      ! scales, traverse positions, coordinate system transformation matrices, etc.
8720      DISP "Storing Header"
8730      ! Set File$ equal to the file name for the header file. Each run number will have a different file name.
8740      File$="R"&VAL$(Run)&Data$
8750      ! Check if the file already exists. If it does then, ask the user if he wants to overwrite the old file.
8760      ON ERROR GOTO 8960
8770      ASSIGN @Data TO File$
8780      OFF ERROR
8790      FOR K=1 TO 10
8800          WAIT .2
8810          BEEP

```

```

8820      NEXT K
8830      INPUT "Over Write old file? (Y or N) ",L$
8840      SELECT L$(1,1)
8850      CASE "Y","y"          ! If the user wants to overwrite the old file, then purge the old file.
8860          ASSIGN @Data TO *
8870          PURGE File$
8880          GOTO 8960
8890      CASE "N","n"          ! If the user doesn't want to overwrite the old file, then have a new run# entered.
8900          CALL Enter_value("Run",Run,"3D.2D")
8910          CALL Enter_value("File",File,"3D")
8920          GOTO Store_header
8930      CASE ELSE
8940          GOTO Store_header
8950      END SELECT
8960      OFF ERROR
8970      Fsize=INT((3200+4000*3+128*4+72*4)/256*1.05+1)      ! Calculate the headers file size.
8980      CREATE BDAT File$,Fsize      ! Create the header's file.
8990      ASSIGN @Data TO File$      ! Open the header's file.
9000      OUTPUT @Data;Array(*),Name$(*),Image$(*),Units$(*)
9010      OUTPUT @Data;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
9020      OUTPUT @Data;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
9030      ASSIGN @Data TO *      ! Close the header's file.
9040      RETURN
9050 Store_file:      ! This subroutine stores the header Array(*), the raw data, and the reduced data onto the disk. There
9060      ! will be one data file for each data set. For example, if the run and file numbers equal 7 and 5
9070      ! respectively, then the data will be stored in a file named "R7F5".
9080      GOSUB Calc      ! Use the tunnel conditions to calculate and/or update other tunnel conditions.
9090      GOSUB Fill      ! Fill Array(*) with the original tunnel conditions along with the updated tunnel conditions.
9100      IF File=1 THEN GOSUB Store_header
9110      DISP "Storing Data"
9120      File$="R"&VAL$(Run)&"F"&VAL$(File)&Data$
9130      ! Check if the file already exists. If it does, then ask the user if he wants to overwrite the old file.
9140      ON ERROR GOTO 9340
9150      ASSIGN @Data TO File$
9160      OFF ERROR
9170      FOR K=1 TO 10
9180          WAIT .2
9190          BEEP
9200      NEXT K
9210      INPUT "Over Write old file? (Y or N) ",L$
9220      SELECT L$(1,1)
9230      CASE "Y","y"          ! If the user wants to overwrite the old file, then purge the old file.
9240          ASSIGN @Data TO *
9250          PURGE File$
9260          GOTO 9340
9270      CASE "N","n"          ! If the user doesn't want to overwrite the old file, then have a new run# entered.
9280          CALL Enter_value("Run",Run,"3D.2D")
9290          CALL Enter_value("File",File,"3D")
9300          GOTO Store_file
9310      CASE ELSE
9320          GOTO Store_file
9330      END SELECT
9340      OFF ERROR
9350      Fsize=INT((3200+Nsam*10*2+60+240)/256*1.05+1)      ! Calculate the data's file size.
9360      CREATE BDAT File$,Fsize      ! Create the data's file.
9370      ASSIGN @Data TO File$      ! Open the data's file.
9380      OUTPUT @Data;Array(*),Raw(*),N(*),Sum(*)
9390      ASSIGN @Data TO *      ! Close the data's file.
9400      RETURN
9410 Read_array:      ! This subroutine reads the header Array(*) off of the disk from a file named "ARRAY". The file will the
9420      ! have default values for the tunnel conditions, laser parameters, graph scales, etc. This file is not
9430      ! meant to be attached to any run number or profile scan. It is used to provide default values for the
9440      ! program so that the user will not have to enter a rather lengthy list of tunnel conditions.
9450      ON ERROR GOTO 9550
9460      ! If the file already exists, then read the Array(*) from the disk.
9470      ASSIGN @File TO "ARRAY"&System$
9480      ENTER @File;Array(*),Name$(*),Image$(*),Units$(*)
9490      ENTER @File;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
9500      ENTER @File;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
9510      ASSIGN @File TO *
9520      OFF ERROR
9530      RETURN
9540      ! If the file doesn't exist then create the file, read in default data, and store the Array(*) on disk.
9550      OFF ERROR
9560      ASSIGN @File TO *
9570      ON ERROR GOTO 9590
9580      PURGE "ARRAY"&System$
9590      OFF ERROR
9600      CALL Array_init(Name$(*),Array(*),Image$(*),Units$(*))
9610      CREATE BDAT "ARRAY"&System$,50

```

```

9620      GOSUB Save_array
9630      RETURN
9640 Save_array: ! This subroutine saves the header Array(*) onto the disk in a file named "ARRAY". The file will then have
9650      ! default values for the tunnel conditions, laser parameters, graph scales, etc. This file is not meant to
9660      ! be attached to any run number or profile scan. It is used to provide default values for the program so
9670      ! that the user will not have to enter a rather lengthy list of tunnel conditions.
9680      ASSIGN @File TO "ARRAY"&System$
9690      OUTPUT @File;Array(*),Name$(*),Image$(*),Units$(*)
9700      OUTPUT @File;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
9710      OUTPUT @File;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
9720      ASSIGN @File TO *
9730      RETURN
9740 Fill: ! This subroutine fills the Array(*) with the current tunnel conditions, laser parameters, and histogram
9750      ! & profile scales.
9760      Array(1,1)=Date ! Date.
9770      Array(1,2)=Mach ! Mach Number.
9780      Array(1,4)=Alpha(1) ! Angle of Attack.
9790      Array(2,1)=Time ! Time.
9800      Array(2,2)=Temp ! Room Temperature (deg. F).
9810      Array(2,4)=Alpha(2) ! Angle of Yaw.
9820      Array(3,1)=Run ! Run Number.
9830      Array(3,2)=Uedge ! Freestream Velocity.
9840      Array(3,4)=Alpha(3) ! Angle of Roll.
9850      Array(4,1)=File ! File Number.
9860      Array(4,2)=Ujet_ue ! Jet exit velocity normalized by Uedge.
9870      Array(4,4)=Theta ! Tx Side Off Axis Angle.
9880      MAT Array(11:14,1)= Mod ! Probe volume positions in MODEL coordinates.
9890      MAT Array(11:14,2)= Tun ! Probe volume positions in TUNNEL coordinates.
9900      MAT Array(11:14,3)= Tcs1 ! Tx side traverse positions in Tcs8 coordinates.
9910      MAT Array(11:14,4)= Tcs2 ! Rx side traverse positions in Tcs8 coordinates.
9920      MAT Array(21,1:3)= Index ! Index of refraction of for laser light (eg: Nair,Nglass,Nwater).
9930      MAT Array(22:24,1:3)= Theta1 ! Angles between LASER & TUNNEL UVW laser beams in Air (N=Index1).
9940      MAT Array(25:27,1:3)= Theta3 ! Angles between LASER & TUNNEL UVW laser beams in Water (N=Index3).
9950      MAT Array(31,1:3)= Beam_spc ! Beam spacing at lens.
9960      MAT Array(32,1:3)= Focl_len ! Focal length.
9970      MAT Array(33,1:3)= Beam_sep ! Beam separation angle in degrees (full angle).
9980      MAT Array(34,1:3)= Wave_len ! Wave length.
9990      MAT Array(35,1:3)= Frng_spc ! Fringe spacing.
10000     MAT Array(36,1:3)= Brg_frq ! Bragg frequency.
10010     MAT Array(37,1:3)= Mix_frq ! Mixing frequency.
10020     MAT Array(38,1:3)= Mea_sgn ! Sign of measured frequency in velocity equation.
10030     MAT Array(39,1:3)= Brg_sgn ! Sign of bragg frequency in velocity equation.
10040     MAT Array(40,1:3)= Mix_sgn ! Sign of mixing frequency in velocity equation.
10050     MAT Array(41,1:3)= Coin ! Coincidence criteria
10060     Array(42,1)=Umin ! Frequency minimum for U calculation.
10070     Array(42,2)=Vmin ! Frequency minimum for V calculation.
10080     Array(42,3)=Wmin ! Frequency minimum for W calculation.
10090     Array(43,1)=Umax ! Frequency maximum for U calculation.
10100     Array(43,2)=Vmax ! Frequency maximum for V calculation.
10110     Array(43,3)=Wmax ! Frequency maximum for W calculation.
10120     Array(51,1)=Nreads ! Number of desired samples.
10130     Array(52,1)=Nsam ! Number of acquired samples.
10140     Array(51,2)=Atime ! Acquisition time.
10150     Array(52,2)=Ctime ! Coincidence time.
10160     Array(51,3)=At_exp ! Acquisition time exponent.
10170     Array(52,3)=Ct_exp ! Coincidence time exponent.
10180     Array(51,4)=Paxis ! Axis for plots.
10190     Array(52,4)=Clip ! Clip: 1 turn histogram clipping on; 0 turns it off.
10200     RETURN
10210 Read: ! This subroutine extracts (reads) the current tunnel conditions, laser parameters, and histogram
10220     ! & profile scales from the Array(*).
10230     Date=TIMEDATE ! Date.
10240     Mach=Array(1,2) ! Mach Number.
10250     Alpha(1)=Array(1,4) ! Angle of Attack.
10260     Time=Date ! Time.
10270     Temp=Array(2,2) ! Room Temperature (deg. F).
10280     Alpha(2)=Array(2,4) ! Angle of Yaw.
10290     Uedge=Array(3,2) ! Freestream Velocity.
10300     Alpha(3)=Array(3,4) ! Angle of Roll.
10310     Ujet_ue=Array(4,2) ! Jet exit velocity normalized by Uedge.
10320     Theta=Array(4,4) ! Tx Side Off Axis Angle.
10330     MAT Mod= Array(11:14,1) ! Probe volume positions in MODEL coordinates.
10340     MAT Tun= Array(11:14,2) ! Probe volume positions in TUNNEL coordinates.
10350     MAT Tcs1= Array(11:14,3) ! Tx side traverse positions in Tcs8 coordinates.
10360     MAT Tcs2= Array(11:14,4) ! Rx side traverse positions in Tcs8 coordinates.
10370     MAT Index= Array(21,1:3) ! Index of refraction of for laser light (eg: Nair,Nglass,Nwater).
10380     MAT Theta1= Array(22:24,1:3) ! Angles between LASER & TUNNEL UVW laser beams in Air (N=Index1).
10390     MAT Theta3= Array(25:27,1:3) ! Angles between LASER & TUNNEL UVW laser beams in Water (N=Index3).
10400     MAT Beam_spc= Array(31,1:3) ! Beam spacing at lens.
10410     MAT Focl_len= Array(32,1:3) ! Focal length.

```



```

10420 MAT Beam_sep= Array(33,1:3) ! Beam separation angle in degrees (full angle).
10430 MAT Wave_len= Array(34,1:3) ! Wave length.
10440 MAT Frng_spc= Array(35,1:3) ! Fringe spacing.
10450 MAT Brg_frq= Array(36,1:3) ! Bragg frequency.
10460 MAT Mix_frq= Array(37,1:3) ! Mixing frequency.
10470 MAT Mea_sgn= Array(38,1:3) ! Sign of measured frequency in velocity equation.
10480 MAT Brg_sgn= Array(39,1:3) ! Sign of bragg frequency in velocity equation.
10490 MAT Mix_sgn= Array(40,1:3) ! Sign of mixing frequency in velocity equation.
10500 MAT Coin= Array(41,1:3) ! Coincidence criteria.
10510 Umin=Array(42,1) ! Frequency minimum for U calculation.
10520 Vmin=Array(42,2) ! Frequency minimum for V calculation.
10530 Wmin=Array(42,3) ! Frequency minimum for W calculation.
10540 Umax=Array(43,1) ! Frequency maximum for U calculation.
10550 Vmax=Array(43,2) ! Frequency maximum for V calculation.
10560 Wmax=Array(43,3) ! Frequency maximum for W calculation.
10570 Nreads=Array(51,1) ! Number of desired samples.
10580 Nsam=Array(52,1) ! Number of acquired samples.
10590 Atime=Array(51,2) ! Acquisition time.
10600 Ctime=Array(52,2) ! Coincidence time.
10610 At_exp=Array(51,3) ! Acquisition time exponent.
10620 Ct_exp=Array(52,3) ! Coincidence time exponent.
10630 Paxis=Array(51,4) ! Axis for plots.
10640 Clip=Array(52,4) ! Clip: 1 turn histogram clipping on; 0 turns it off.
10650 RETURN
10660 Calc: ! This subroutine uses the current tunnel conditions and laser parameters to calculate and/or update other
10670 ! tunnel conditions and laser parameters.
10680 FOR K=1 TO 3
10690 IF K=2 THEN
10700 Beam1=Theta+ATN(Beam_spc(K)/2/Focl_len(K)) ! Angles of the off axis beam pair in air.
10710 Beam2=Theta-ATN(Beam_spc(K)/2/Focl_len(K))
10720 ELSE
10730 Beam1=0+ATN(Beam_spc(K)/2/Focl_len(K)) ! Angles of the on axis beam pairs in air.
10740 Beam2=0-ATN(Beam_spc(K)/2/Focl_len(K))
10750 END IF
10760 Beam1=ASN(Index(1)/Index(3)*SIN(Beam1)) ! Angle of the beam pairs in water.
10770 Beam2=ASN(Index(1)/Index(3)*SIN(Beam2))
10780 Beam_sep(K)=Beam1-Beam2 ! Beam pair separation angle.
10790 Frng_spc(K)=Wave_len(K)/Index(3)/(2*SIN(Beam_sep(K)/2))/1000 ! Fringe spacing in water (um).
10800 NEXT K
10810 MAT Array(33,1:3)= Beam_sep ! Beam separation angle in degrees (full angle).
10820 MAT Array(35,1:3)= Frng_spc ! Fringe spacing.
10830 ! Calculate the TCS to TUNNEL (and visa versa) traverse coordinate system transformation matrices.
10840 Fs=Focl_len(1) ! Focal length of sending side lenses (inches).
10850 Fr=Focl_len(1) ! Focal length of receiving side lenses (inches).
10860 Bs=Beam_spc(1) ! Beam spacing at sending lenses (inches).
10870 Br=3 ! Receiving side lens diameter (inches).
10880 Ts=Theta ! Sending side off axis angle (degrees).
10890 Tr=17.05 ! Receiving side off axis angle (degrees).
10900 ! Ta is the offaxis sending side auxiliary rotation angle (degrees).
10910 CALL Ctm_tcs(Tcs2tun1(*),Tcs2tun2(*),Tun2tcs1(*),Tun2tcs2(*),Fs,Fr,Bs,Br,Index(*),Ts,Tr,Ta)
10920 ! Calculate the LASER to TUNNEL (and visa versa) velocity coordinate system transformation matrices.
10930 CALL Refract(Index(*),Theta1(*),Theta3(*))
10940 CALL Ctm_ldv(Theta3(*),Tun2ldv(*),Ldv2tun(*))
10950 ! Calculate the TUNNEL to MODEL (and visa versa) coordinate system transformation matrices.
10960 CALL Ctm_mod(Alpha(*),Mod2tun(*),Tun2mod(*))
10970 ! Define the coincidence mask depending on the value of Coin(*).
10980 Cmask=Coin(1)*1+Coin(2)*2+Coin(3)*4
10990 ! Define Paxis$ depending on the value of Paxis.
11000 SELECT Paxis
11010 CASE 1
11020 Paxis$="X"
11030 CASE 2
11040 Paxis$="Y"
11050 CASE 3
11060 Paxis$="Z"
11070 CASE 4
11080 Paxis$="A"
11090 CASE ELSE
11100 Paxis=2
11110 Paxis$="Y"
11120 GOSUB M3k4
11130 END SELECT
11140 ! If the Run number or File number have not been defined then have the user enter their values.
11150 IF Run=0 OR File=0 THEN
11160 CALL Enter_value("Run Number ",Run,"3D.2D")
11170 CALL Enter_value("File Number ",File,"3D")
11180 GOTO 11150
11190 END IF
11200 RETURN
11210 END

```

```

11220 Do_nothing: SUB Do_nothing
11230 ! Description:
11240 ! This subprogram is called when the keys on the keyboard are pressed during TCS8 traverse
11250 ! movements. This is done so that any STOP, PAUSE, or RESET keys will be ignored. This prevents
11260 ! stopping the program while the HP and TCS8 are communicating with each other. Otherwise, they
11270 ! might get out of sync while communicating resulting in system hang ups.
11280 ! Variables:
11290 ! KS String used to flush the keyboard buffer.
11300 KS=KBD$
11310 SUBEND
11320 Menu: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
11330 Menu_read: SUB Menu_read(Menu$(*))
11340 ! Description:
11350 ! This subprogram reads in the menu descriptors for each entry of the five menus.
11360 ! Variables:
11370 ! Menu Used as an index to the string array Menu$(*).
11380 ! Key Used as an index to the string array Menu$(*).
11390 ! Menu$(*) String array where each element describes its corresponding menu subroutine's function.
11400 ! L$ String use to read in the menu descriptor from the data statements.
11410 OPTION BASE 1
11420 DIM L$(80)
11430 ! Fill all of the menu entry's descriptions with "MxKx".
11440 FOR Menu=1 TO SIZE(Menu$,1)
11450 FOR Key=1 TO 8
11460 Menu$(Menu,Key)="M"&VAL$(Menu)&"K"&VAL$(Key)&":"
11470 NEXT Key
11480 NEXT Menu
11490 ON ERROR GOTO 11570 ! The following while loop will get error#36 when the data statements run out.
11500 ! For each menu and key, enter the menu entry's description.
11510 WHILE 1=1
11520 READ L$
11530 Menu=VAL(L$[2,2])
11540 Key=VAL(L$[4,4])
11550 Menu$(Menu,Key)=L$
11560 END WHILE
11570 SUBEXIT
11580 DATA "M1K1: Laser Alignment"
11590 DATA "M2K1: Return to main menu"
11600 DATA "M2K2: Sides : Tx & Rx"
11610 DATA "M2K3: Coordinates: MODEL"
11620 DATA "M2K4: Mode : ABSOLUTE"
11630 DATA "M2K5: Move X"
11640 DATA "M2K6: Move Y"
11650 DATA "M2K7: Move Z"
11660 DATA "M2K8: Move A"
11670 DATA "M1K2: Pre Run"
11680 DATA "M3K1: Return to MAIN menu"
11690 DATA "M3K2: Enter Run & File Numbers"
11700 DATA "M3K3: Enter Number of Samples"
11710 DATA "M3K4: Select Traverse Axis for Profile"
11720 DATA "M3K5: Print Coordinate Transformation Matrices"
11730 DATA "M3K6: Setup Graphics"
11740 DATA "M3K7: Tunnel Conditions"
11750 DATA "M4K1: Return to PRE RUN menu"
11760 DATA "M4K2: Load Tunnel Conditions"
11770 DATA "M4K3: Save Tunnel Conditions"
11780 DATA "M4K4: Print Tunnel Conditions"
11790 DATA "M4K5: Enter Tunnel Condition Data"
11800 DATA "M4K6: Enter Tunnel Condition Names"
11810 DATA "M4K7: Enter Tunnel Condition Units"
11820 DATA "M4K8: Enter Tunnel Condition Images"
11830 DATA "M3K8: Traverse"
11840 DATA "M5K1: Return to PRE RUN menu"
11850 DATA "M5K2: View & Set TCS8 Positions"
11860 DATA "M5K3: View & Set TCS8 Units"
11870 DATA "M5K4: View & Set TCS8 Revolution"
11880 DATA "M5K5: View & Set TCS8 Velocity"
11890 DATA "M5K6: View & Set TCS8 Acceleration"
11900 DATA "M1K3: Post Run (Dump Graphics)"
11910 DATA "M1K4: Set Auto Move Positions"
11920 DATA "M1K5: Move traverse"
11930 DATA "M1K6: Take data"
11940 DATA "M1K7: Auto move and take"
11950 DATA "M1K8: Display Histograms"
11960 SUBEND
11970 Menu_disp: SUB Menu_disp(Menu,Menu$(*))
11980 ! Description:
11990 ! This subprogram displays the current menu at the top of the CRT.
12000 ! Variables:
12010 ! Menu Used as an index to the string array Menu$(*).

```

```

12020      !      Key      Used as an index to the string array Menu$(*).
12030      !      Menu$(*) String array where each element describes its corresponding menu subroutine's function.
12040      PRINTER IS CRT
12050      PRINT CHR$(128);          ! Turn off inverse video if it is on.
12060      IF Menu=0 THEN Menu=1
12070      FOR Key=1 TO 8
12080          Menu$(Menu,Key)=Menu$(Menu,Key)&RPT$( " ",50-LEN(Menu$(Menu,Key)))
12090          PRINT TABXY(1,Key);Menu$(Menu,Key) [3]
12100      NEXT Key
12110      SUBEND
12120 Menu_status: SUB Menu_status(Menu,Key, Pen,Menu$(*))
12130      !      Description:
12140      !      This subprogram displays the current menu selection in normal or inverse video. The inverse
12150      !      video text style indicates that the subroutine for the current menu selection is busy. The
12160      !      normal text style indicates that the subroutine for the current menu selection is has completed.
12170      !      Variables:
12180      !      Menu      Indicates which of the menus has been selected as the current menu.
12190      !      Key       Indicates which one of eight menu subroutines in the menu is to be executed.
12200      !      Pen       Indicates Busy/Ready Status. Pen=0 for busy. Pen=1 for ready.
12210      !      Menu$(*) String array where each element describes its corresponding menu subroutine's function.
12220      PRINTER IS CRT
12230      PRINT TABXY(1,Key);CHR$(129-Pen);Menu$(Menu,Key) [3];CHR$(128)
12240      WAIT .1
12250      SUBEND
12260 Enter:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
12270 Enter_value: SUB Enter_value(Name$,Value,Image$)
12280      !      Description:
12290      !      This subprogram displays the current value of a variable and then has the user enter its new
12300      !      value. The old value will be kept if the RETURN key is pressed and no data is entered.
12310      !      Variables:
12320      !      Name$      Name of the variable.
12330      !      Image$     Image format of the variable. Used for printing the variable with a format.
12340      !      Value     Contains the initial value and then the updated value for the variable.
12350      IF Name$="Date" OR Name$="Time" THEN SUBEXIT
12360      DISP CHR$(129);
12370      DISP USING 12380;Name$
12380      IMAGE #,"Old ",K,"="
12390      IF Image$<>"" THEN DISP USING "#,"&Image$;Value
12400      IF Image$="" THEN DISP USING "#,K";Value
12410      DISP USING 12420;Name$
12420      IMAGE #,"      Enter new ",K
12430      INPUT " ? ",Value
12440      DISP CHR$(128);
12450      SUBEND
12460 Enter_string: SUB Enter_string(Name$,Value$,Image$)
12470      !      Description:
12480      !      This subprogram displays the current value of a string variable and then has the user enter its
12490      !      new value. The old value will be kept if the RETURN key is pressed and no data is entered.
12500      !      Variables:
12510      !      Name$      Name of the variable.
12520      !      Value$     Contains the initial value and then the updated value for the string variable.
12530      DISP CHR$(129);
12540      DISP USING 12550;Name$
12550      IMAGE #,"Old ",K,"="
12560      DISP USING "#,"&Image$;Value$
12570      DISP USING 12580;Name$
12580      IMAGE #,"      Enter new ",K
12590      INPUT " ? ",Value$
12600      DISP CHR$(128);
12610      SUBEND
12620 Array:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
12630 Array_init: SUB Array_init(Name$(*),Array(*),Image$(*),Units$(*))
12640      !      Description:
12650      !      This subprogram reads in default data for each of the variable's names, values, image formats,
12660      !      and units. These variables include, but are not limited to, the tunnel conditions, laser
12670      !      parameters, graph scales, traverse positions, and coordinate system transformation matrices.
12680      !      Variables:
12690      !      Array(*)    Array of tunnel conditions, laser parameters, graph scales, etc.
12700      !      Name$(*)   Names for the variables in Array(*).
12710      !      Image$(*)  Image formats for the variables in Array(*).
12720      !      Units$(*)  Units for the variables in Array(*).
12730      !      X         Used as an index to the above arrays and string arrays.
12740      !      Y         Used as an index to the above arrays and string arrays.
12750      !      Before     Number of digits before the decimal point in the image format.
12760      !      After      Number of digits after the decimal point in the image format.
12770      ON ERROR GOTO 12950
12780      READ Y
12790      FOR X=1 TO SIZE(Name$,2)
12800          READ Name$(Y,X),Array(Y,X),Image$(Y,X),Units$(Y,X)
12810          SELECT Image$(Y,X)

```

```

12820 CASE "0"
12830     Image$(Y,X)="9D"
12840 CASE "1" TO "7"
12850     After=VAL(Image$(Y,X))
12860     Before=8-After
12870     Image$(Y,X)=VAL$(Before)&"D."&VAL$(After)&"D"
12880 CASE "K"
12890 CASE "N"
12900 CASE ELSE
12910     Image$(Y,X)="9D"
12920 END SELECT
12930 NEXT X
12940 GOTO 12780
12950 SUBEXIT
12960 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
12970 DATA 1, Date , 0,0,"", Mach , 0,4,"", "" , 0,0,"", Alpha1 , 0,4,3
12980 DATA 2, Time , 0,0,"", Temp , 68.5,4,3F, "" , 0,0,"", Alpha2 , 0,4,3
12990 DATA 3, Run , 0,2,"", Uedge ,.0762,4,m/s, "" , 0,0,"", Alpha3 , 0,4,3
13000 DATA 4, File , 0,0,"", Ujet/Ue , 0,4,m/s, "" , 0,0,"", Theta , 45,4,3
13010 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13020 DATA 11, Xmod , 0,4,in , Xtun , 0,4,in , Xltcs , 0,4,in , X2tcs , 0,4,in
13030 DATA 12, Ymod , 0,4,in , Ytun , 0,4,in , Yltcs , 0,4,in , Y2tcs , 0,4,in
13040 DATA 13, Zmod , 0,4,in , Ztun , 0,4,in , Zltcs , 0,4,in , Z2tcs , 0,4,in
13050 DATA 14, Amod , 0,4,in , Atun , 0,4,in , Altcs , 0,4,in , A2tcs , 0,4,in
13060 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13070 DATA 21, Index1 ,1.000,3,"", Index2 ,1.430,3,"", Index3 ,1.333,3,"", "" , 0,0,""
13080 DATA 22, Theta1AU, 0,4,3 , Theta1AV, 90,4,3 , Theta1AW, 90,4,3 , "" , 0,0,""
13090 DATA 23, Theta1BU, 45,4,3 , Theta1BV, 135,4,3 , Theta1BW, 90,4,3 , "" , 0,0,""
13100 DATA 24, Theta1CU, 90,4,3 , Theta1CV, 90,4,3 , Theta1CW, 0,4,3 , "" , 0,0,""
13110 DATA 25, Theta3AU, 0,4,3 , Theta3AV, 90,4,3 , Theta3AW, 90,4,3 , "" , 0,0,""
13120 DATA 26, Theta3BU, 45,4,3 , Theta3BV, 135,4,3 , Theta3BW, 90,4,3 , "" , 0,0,""
13130 DATA 27, Theta3CU, 90,4,3 , Theta3CV, 90,4,3 , Theta3CW, 0,4,3 , "" , 0,0,""
13140 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13150 DATA 31, UBeamSpc,2.362,3,in , VBeamSpc,2.362,3,in , WBeamSpc,2.362,3,in , "" , 0,0,""
13160 DATA 32, UFocLen,19.413,3,in , VFocLen,19.413,3,in , WFocLen,19.413,3,in , "" , 0,0,""
13170 DATA 33, UBeamSep,0.000,3,3 , VBeamSep,0.000,3,3 , WBeamSep,0.000,3,3 , "" , 0,0,""
13180 DATA 34, UWaveLen,476.5,3,nm , VWaveLen,514.5,3,nm , WWaveLen,488.0,3,nm , "" , 0,0,""
13190 DATA 35, UFrngSpc,00.00,3,um , VFrngSpc,00.00,3,um , WFrngSpc,00.00,3,um , "" , 0,0,""
13200 DATA 36, Ubrag ,40.00,4,MHz, Vbrag ,40.00,4,MHz, Wbrag ,40.00,4,MHz, "" , 0,0,""
13210 DATA 37, Umix ,39.90,4,MHz, Vmix ,39.90,4,MHz, Wmix ,39.90,4,MHz, "" , 0,0,""
13220 DATA 38, UmeaSgn ,+1,0,"", VmeaSgn ,+1,0,"", WmeaSgn ,+1,0,"", "" , 0,0,""
13230 DATA 39, UbrgSgn , -1,0,"", VbrgSgn , -1,0,"", WbrgSgn , -1,0,"", "" , 0,0,""
13240 DATA 40, UmixSgn ,+1,0,"", VmixSgn ,+1,0,"", WmixSgn ,+1,0,"", "" , 0,0,""
13250 DATA 41, U coin , 1,0,"", V coin , 1,0,"", W coin , 1,0,"", "" , 0,0,""
13260 DATA 42, UFreqMin, -99,4,MHz, VFreqMin, -99,4,MHz, WFreqMin, -99,4,MHz, "" , 0,0,""
13270 DATA 43, UFreqMax, 99,4,MHz, VFreqMax, 99,4,MHz, WFreqMax, 99,4,MHz, "" , 0,0,""
13280 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13290 DATA 51, Nreads , 1000,0,"", Atime , 30,6,s , AExp , 10,0,"", Paxis , 2,0,""
13300 DATA 52, Nsam , 1000,0,"", Ctime , 1E-2,6,s , CExp , 5,0,"", Clip , 0,0,""
13310 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13320 DATA 61, Xmin1 , 0.00,1,"", Xmax1 , 1.00,1,"", Ymin1 , 0,0,"", Ymax1 , 100,0,""
13330 DATA 62, Xmin2 , 0.00,1,"", Xmax2 , 1.00,1,"", Ymin2 , 0,0,"", Ymax2 , 100,0,""
13340 DATA 63, Xmin3 , 0.00,1,"", Xmax3 , 1.00,1,"", Ymin3 , 0,0,"", Ymax3 , 100,0,""
13350 DATA 64, Xmin4 , 0.00,1,"", Xmax4 , 2.00,1,"", Ymin4 , 0,0,"", Ymax4 , 100,0,""
13360 DATA 65, Xmin5 , 0.00,1,"", Xmax5 , 2.00,1,"", Ymin5 , 0,0,"", Ymax5 , 100,0,""
13370 DATA 66, Xmin6 , 0,1,"", Xmax6 , 3,1,"", Ymin6 , -1.5,2,"", Ymax6 , 1.5,2,""
13380 DATA 67, Xmin7 , 0,1,"", Xmax7 , 5,1,"", Ymin7 , -1.5,2,"", Ymax7 , 1.5,2,""
13390 DATA 68, Xmin8 , -.025,3,"", Xmax8 , .025,3,"", Ymin8 , -1.5,2,"", Ymax8 , 1.5,2,""
13400 DATA 69, Xmin9 , -.1,2,"", Xmax9 , .1,2,"", Ymin9 , -1.5,2,"", Ymax9 , 1.5,2,""
13410 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13420 DATA 71, Xmin1 , 935,0,pxl, Xmax1 , 1235,0,pxl, Ymin1 , 725,0,pxl, Ymax1 , 825,0,pxl
13430 DATA 72, Xmin2 , 935,0,pxl, Xmax2 , 1235,0,pxl, Ymin2 , 585,0,pxl, Ymax2 , 685,0,pxl
13440 DATA 73, Xmin3 , 935,0,pxl, Xmax3 , 1235,0,pxl, Ymin3 , 445,0,pxl, Ymax3 , 545,0,pxl
13450 DATA 74, Xmin4 , 935,0,pxl, Xmax4 , 1235,0,pxl, Ymin4 , 305,0,pxl, Ymax4 , 405,0,pxl
13460 DATA 75, Xmin5 , 935,0,pxl, Xmax5 , 1235,0,pxl, Ymin5 , 165,0,pxl, Ymax5 , 265,0,pxl
13470 DATA 76, Xmin6 , 75,0,pxl, Xmax6 , 325,0,pxl, Ymin6 , 525,0,pxl, Ymax6 , 825,0,pxl
13480 DATA 77, Xmin7 , 425,0,pxl, Xmax7 , 675,0,pxl, Ymin7 , 525,0,pxl, Ymax7 , 825,0,pxl
13490 DATA 78, Xmin8 , 75,0,pxl, Xmax8 , 325,0,pxl, Ymin8 , 165,0,pxl, Ymax8 , 465,0,pxl
13500 DATA 79, Xmin9 , 425,0,pxl, Xmax9 , 675,0,pxl, Ymin9 , 165,0,pxl, Ymax9 , 465,0,pxl
13510 ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
13520 DATA 81, Xdiv1 , 5,0,"", Ydiv1 , 4,0,"", Xdiv6 , 6,0,"", Ydiv6 , 6,0,""
13530 DATA 82, Xdiv2 , 5,0,"", Ydiv2 , 4,0,"", Xdiv7 , 5,0,"", Ydiv7 , 6,0,""
13540 DATA 83, Xdiv3 , 5,0,"", Ydiv3 , 4,0,"", Xdiv8 , 2,0,"", Ydiv8 , 6,0,""
13550 DATA 84, Xdiv4 , 5,0,"", Ydiv4 , 4,0,"", Xdiv9 , 4,0,"", Ydiv9 , 6,0,""
13560 DATA 85, Xdiv5 , 5,0,"", Ydiv5 , 4,0,"", "" , 0,0,"", "" , 0,0,""
13570 SUBEND
13580 Array_print: SUB Array_print(Array(*),Name$(*),Image$(*),Units$(*))
13590 ! Description:
13600 ! This subprogram prints the values of each of the variables with their names, image formats, and
13610 ! units. These variables include, but are not limited to, the tunnel conditions, laser

```

```

13620 ! parameters, and graph scales.
13630 ! Variables:
13640 ! Array(*) Array of tunnel conditions, laser parameters, graph scales, etc.
13650 ! Name$(*) Names for the variables in Array(*).
13660 ! Image$(*) Image formats for the variables in Array(*).
13670 ! Units$(*) Units for the variables in Array(*).
13680 ! X Used as in index to the above arrays and string arrays.
13690 ! Y Used as in index to the above arrays and string arrays.
13700 PRINT USING "#,5/"
13710 FOR Y=1 TO SIZE(Array,1)
13720 MAT SEARCH Array(Y,*),#LOC(<>0);L1
13730 MAT SEARCH Name$(Y,*),#LOC(<>"");L2
13740 IF L1+L2=0 AND L3=0 THEN 13980
13750 L3=L1+L2
13760 PRINT USING "#,28X"
13770 FOR X=1 TO SIZE(Array,2)
13780 SELECT Name$(Y,X)
13790 CASE ""
13800 PRINT USING "#,28X"
13810 CASE "Date"
13820 L$=DATES(Array(Y,X))
13830 L$=L$[1,2]&L$[4,6]&L$[8,11]
13840 PRINT USING "#,8A,A,9A,X,3A,6X";TRIMS(Name$(Y,X)), "=",L$,Units$(Y,X)
13850 CASE "Time"
13860 L$=" "&TIMES(Array(Y,X))
13870 PRINT USING "#,8A,A,9A,X,3A,6X";TRIMS(Name$(Y,X)), "=",L$,Units$(Y,X)
13880 CASE ELSE
13890 IF Image$(Y,X)="" THEN Image$(Y,X)="9D"
13900 ON ERROR GOTO 13930
13910 PRINT USING "#,8A,A,"&Image$(Y,X)&","X,3A,6X";TRIMS(Name$(Y,X)), "=",Array(Y,X),Units$(Y,X)
13920 GOTO 13950
13930 OFF ERROR
13940 PRINT USING "#,8A,A,K,X,3A,6X";TRIMS(Name$(Y,X)), "=",Array(Y,X),Units$(Y,X)
13950 END SELECT
13960 NEXT X
13970 PRINT
13980 NEXT Y
13990 SUBEND
14000 Change: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
14010 Change: SUB Change(Type$,Array(*),Name$(*),Image$(*),Units$(*))
14020 ! Description:
14030 ! This subprogram displays on the CRT the values of each of the variables with their names,
14040 ! image formats, and units. The user can select one of the variables and enter a new value,
14050 ! name, image format, or units. The user selects the particular variable by using the
14060 ! left, right, up, and down cursor keys. The selected variable will appear in inverse video.
14070 ! When it is not selected, it will appear in normal text. When the user has selected the
14080 ! appropriate variable he should then press the "Select" key on the keyboard. Then, depending on
14090 ! the value of Type$ he will be asked to enter a new value, name, image format, or units. To
14100 ! exit the change variables mode press the "Escape" key.
14110 ! There are three types of data that are passed to the subprogram. The first type of data
14120 ! includes, but is not limited to, the tunnel conditions, laser parameters, and graph scales.
14130 ! With this first type the user is allowed to enter new variable values, names, image formats, and
14140 ! units. The second type of data is the "Auto Move and Take" data. These data are for the pre
14150 ! programed traverse positions used in a profile scan. The third type of data is the "View and
14160 ! Set TCS8 parameters" data acquired from and then sent back to the TCS8.
14170 ! Variables:
14180 ! Array(*) Array whose values, names, image formats, or units are to be modified.
14190 ! Name$(*) Names for the variables in Array(*).
14200 ! Image$(*) Image formats for the variables in Array(*).
14210 ! Units$(*) Units for the variables in Array(*).
14220 ! Type$ Indicates which type of data is to be entered.
14230 ! Type$="VALUES" has the user enter a new value for the selected variable.
14240 ! Type$="NAMES" has the user enter a new name for the selected variable.
14250 ! Type$="IMAGES" has the user enter a new image format for the selected variable.
14260 ! Type$="UNITS" has the user enter a new units for the selected variable.
14270 ! X Used as in index to the above arrays and string arrays.
14280 ! Y Used as in index to the above arrays and string arrays.
14290 PRINTER IS CRT
14300 FOR Y=1 TO SIZE(Array,1)
14310 FOR Y1=Y TO SIZE(Array,1)
14320 FOR X=1 TO SIZE(Array,2)
14330 IF Name$(Y1,X)<>"" THEN 14380
14340 NEXT X
14350 NEXT Y1
14360 CLEAR SCREEN
14370 SUBEXIT
14380 FOR Y2=Y1 TO SIZE(Array,1)
14390 FOR X=1 TO SIZE(Array,2)
14400 IF Name$(Y2,X)<>"" THEN 14430
14410 NEXT X

```

```

14420          GOTO 14440
14430      NEXT Y2
14440      FOR Y2=Y2 TO SIZE(Array,1)
14450          FOR X=1 TO SIZE(Array,2)
14460              IF Name$(Y2,X) <> "" THEN 14490
14470              NEXT X
14480          NEXT Y2
14490          Y2=Y2-1
14500          CLEAR SCREEN
14510          CALL Display(Type$,Y1,Y2,Array(*),Name$(*),Image$(*),Units$(*))
14520          Done=0
14530          X=1
14540          Y=Y1
14550          ON KBD ALL,15 GOSUB Kbd
14560 Wait:      IF NOT Done THEN Wait
14570          OFF KBD
14580          CLEAR SCREEN
14590          Y=Y2
14600      NEXT Y
14610      SUBEXIT
14620 Kbd:      CALL Update(Type$,X,Y,Y1,Y2,Done,Array(*),Name$(*),Image$(*),Units$(*))
14630      RETURN
14640  SUBEND
14650 Display:  SUB Display(Type$,Y1,Y2,Array(*),Name$(*),Image$(*),Units$(*))
14660      ! Description:
14670      !     This subprogram displays on the CRT the values of each of variables with their names, image
14680      !     formats, and units.
14690      ! Variables:
14700      !     Array(*)   Array whose values, names, image formats, or units are to be modified.
14710      !     Name$(*)   Names for the variables in Array(*).
14720      !     Image$(*)  Image formats for the variables in Array(*).
14730      !     Units$(*)  Units for the variables in Array(*).
14740      !     Type$      Indicates which type of data is to be entered.
14750      !                 Type$="VALUES" has the user enter a new value for the selected variable.
14760      !                 Type$="NAMES"  has the user enter a new name for the selected variable.
14770      !                 Type$="IMAGES" has the user enter a new image format for the selected variable.
14780      !                 Type$="UNITS"  has the user enter a new units for the selected variable.
14790      !     X          Used as in index to the above arrays and string arrays.
14800      !     Y          Used as in index to the above arrays and string arrays.
14810      FOR Y=Y1 TO Y2
14820          FOR X=1 TO SIZE(Array,2)
14830              CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
14840          NEXT X
14850      NEXT Y
14860      CALL Select(Type$,1,Y1,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
14870  SUBEND
14880 Select:  SUB Select(Type$,X,Y,Y1,Y2,C,Array(*),Name$(*),Image$(*),Units$(*))
14890      ! Description:
14900      !     This subprogram displays on the CRT the value of one variable along with its names, image
14910      !     format, and units.
14920      ! Variables:
14930      !     Array(*)   Array whose values, names, image formats, or units are to be modified.
14940      !     Name$(*)   Names for the variables in Array(*)
14950      !     Image$(*)  Image formats for the variables in Array(*)
14960      !     Units$(*)  Units for the variables in Array(*)
14970      !     Type$      Indicates which type of data is to be entered.
14980      !                 Type$="VALUES" has the user enter a new value for the selected variable.
14990      !                 Type$="NAMES"  has the user enter a new name for the selected variable.
15000      !                 Type$="IMAGES" has the user enter a new image format for the selected variable.
15010      !                 Type$="UNITS"  has the user enter a new units for the selected variable.
15020      !     X          Used as in index to the above arrays and string arrays.
15030      !     Y          Used as in index to the above arrays and string arrays.
15040      PRINT CHR$(128+C);TABXY(26*X-24,15+Y-Y1+1);
15050      PRINT RPT$( " ",23);TABXY(26*X-24,15+Y-Y1+1);
15060      IF Name$(Y,X)="" AND Array(Y,X)=0 THEN 15260
15070      Img$=Image$(Y,X)
15080      Unt$=Units$(Y,X)
15090      IF Image$(Y,X)="" THEN Img$="K"
15100      IF Units$(Y,X)="" THEN Unt$=" "
15110      SELECT Type$
15120      CASE "VALUES"
15130          SELECT Name$(Y,X)
15140          CASE "Date"
15150          CASE "Time"
15160          CASE ELSE
15170              PRINT USING "#,10A,A,"&Img$&","X,3A";Name$(Y,X),":",Array(Y,X),Unt$
15180          END SELECT
15190      CASE "NAMES"
15200          PRINT USING "#,10A,A,8A";Name$(Y,X),":",Name$(Y,X)
15210      CASE "UNITS"

```

```

15220      PRINT USING "#,10A,A,8A";Name$(Y,X),":",Units$(Y,X)
15230      CASE "IMAGES"
15240      PRINT USING "#,10A,A,8A";Name$(Y,X),":",Image$(Y,X)
15250      END SELECT
15260      PRINT CHR$(128);
15270      SUBEND
15280 Update: SUB Update(Type$,X,Y,Y1,Y2,Done,Array(*),Name$(*),Image$(*),Units$(*))
15290      ! Description:
15300      ! This subprogram scrolls through the variables displayed on the CRT and has the user enter
15310      ! updated values. The user can select one of the variables and enter a new value, name, image
15320      ! format, or units. The user selects the particular variable by using the left, right, up, down
15330      ! cursor keys. This subprogram will only have been called after a keyboard key has been pressed.
15340      ! If a cursor key has been pressed then the previously selected variable will be redisplayed in
15350      ! normal text and the new selected variable will appear in inverse video text. When the user has
15360      ! selected the appropriate variable he will have pressed the "Select" key on the keyboard. Then,
15370      ! depending on the value of the Type$ he will be asked to enter a new value, name, image format,
15380      ! or units. To exit the change variables mode the user will have pressed the "Escape" key.
15390      ! Variables:
15400      ! Array(*) Array of tunnel conditions, laser parameters, graph scales, etc.
15410      ! Name$(*) Names for the variables in Array(*).
15420      ! Image$(*) Image formats for the variables in Array(*).
15430      ! Units$(*) Units for the variables in Array(*).
15440      ! Type$ Indicates which type of data is to be entered.
15450      ! Type$="VALUES" has the user enter a new value for the selected variable.
15460      ! Type$="NAMES" has the user enter a new name for the selected variable.
15470      ! Type$="IMAGES" has the user enter a new image format for the selected variable.
15480      ! Type$="UNITS" has the user enter a new units for the selected variable.
15490      ! X Used as in index to the above arrays and string arrays.
15500      ! Y Used as in index to the above arrays and string arrays.
15510      DISABLE
15520      K$=KBD$
15530      IF K$="" THEN 15990
15540      SELECT NUM(K$(1,1))
15550      CASE 27 ! ESC
15560      Done=1
15570      CASE 255
15580      CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
15590      SELECT NUM(K$(2,2))
15600      CASE 73,80 ! Break,Stop
15610      PAUSE
15620      CASE 124 ! Menu
15630      Done=1
15640      CASE 38 ! Select
15650      CALL Select(Type$,X,Y,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
15660      SELECT Type$
15670      CASE "VALUES"
15680      IF Name$(Y,X)="" THEN CALL Enter_string("Name for "&Name$(Y,X),Name$(Y,X),"K")
15690      IF Image$(Y,X)="" THEN CALL Enter_string("Image for "&Name$(Y,X),Image$(Y,X),"K")
15700      CALL Enter_value(Name$(Y,X),Array(Y,X),Image$(Y,X))
15710      CASE "NAMES"
15720      CALL Enter_string("Name for "&Name$(Y,X),Name$(Y,X),"K")
15730      CASE "UNITS"
15740      CALL Enter_string("Units for "&Name$(Y,X),Units$(Y,X),"K")
15750      CASE "IMAGES"
15760      CALL Enter_string("Image for "&Name$(Y,X),Image$(Y,X),"K")
15770      END SELECT
15780      CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
15790      IF X=SIZE(Array,2) THEN Y=Y+1
15800      X=X+1
15810      CASE 60 ! Left
15820      X=X-1
15830      CASE 62 ! Right
15840      X=X+1
15850      CASE 94 ! Up
15860      Y=Y-1
15870      CASE 86 ! Down
15880      Y=Y+1
15890      CASE 92 ! First
15900      X=1
15910      Y=1
15920      END SELECT
15930      X=(X-1) MOD SIZE(Array,2)+1
15940      Y=(Y-1+1-1) MOD (Y2-Y1+1)+Y1
15950      IF X<1 THEN X=SIZE(Array,2)
15960      IF Y<Y1 THEN Y=Y2
15970      CALL Select(Type$,X,Y,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
15980      END SELECT
15990      ENABLE
16000      SUBEXIT
16010      SUBEND

```

```

16020 Misc:
16030 Convert2words: SUB Convert2words(Real,INTEGER High,Low)
16040      ! Description:
16050      !   This subprogram converts a single real precision variable into two 16 bit words. The initial
16060      !   real precision variables is converted in to a 32 bit integer and then separated into high and
16070      !   low 16 bit integers. The most significant 16 bits will be in the "High" variable while the
16080      !   least significant 16 bits will be placed the the "Low" variable. The main purpose of this
16090      !   subprogram is to provide a means to send a 32 bit integer to the LVDAS over the 16 bit high
16100      !   speed interface.
16110      ! Variables:
16120      !   Real   Initial real precision value for the variable.
16130      !   Hex$   Hex value of "Real". String length will be 8 bytes for 32 bits.
16140      !   High   Most significant 16 bits of integerized "Real".
16150      !   Low    Least significant 16 bits of integerized "Real".
16160      Hex$=DVAL$(Real,16)
16170      High=IVAL(Hex$[1,4],16)
16180      Low=IVAL(Hex$[5,8],16)
16190      SUBEND
16200 Error:
16210      SUB Error
16220      ! Description:
16230      !   This subprogram will print an error message when ever a program error occurs. The error message
16240      !   will be displayed at the top of the CRT and also printed on the printers paper. Such errors
16250      !   might occur when data to be printed will not fit in the image formats. Other errors will also
16260      !   generate a displayed and printed error message.
16270      BEEP
16280      DISP ERRMS
16290      OUTPUT PRT;ERRMS
16300      PRT=VAL(SYSTEM$( "PRINTER IS" ))
16310      PRINTER IS CRT
16320      PRINT TABXY(95,1);ERRMS
16330      PRINTER IS PRT
16340      ERROR SUBEXIT
16350      SUBEND
16360 Scale:
16370      SUB Scale(G)
16380      ! Description:
16390      !   This subprogram selects one of nine histogram or profile plots. The plot's area of the CRT is
16400      !   selected and scaled to the appropriate scales.
16410      !   OPTION BASE 1
16420      !   COM /Graph/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),
16430      !   Legend$(*)
16440      !   VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
16450      !   WINDOW Wndw(G,1),Wndw(G,2),Wndw(G,3),Wndw(G,4)
16460      SUBEND
16470 Table:
16480      SUB Table(Table(*))
16490      ! Description:
16500      !   This subprogram is used to create a lookup table array. The lookup table array facilitates
16510      !   the rapid conversion of raw encoded Macrodyne data into a usable frequency. Once the table
16520      !   has been filled, then the raw Macrodyne data can be used as an index to the table array.
16530      ! Variables:
16540      !   Table(*)   Lookup table of frequencies.
16550      !   Mantissa(*) The 10 bit mantissa part of the raw Macrodyne data (0..1023).
16560      !   Fringes    The 1 bit Fringe Count part of the raw Macrodyne data (0:16, 1:8 fringes).
16570      !   Exponent    The 4 bit Exponent part of the raw Macrodyne data.
16580      !   Time(*)     An array of measurement times for a given number of Fringes and Exponent.
16590      !   Freq(*)     An array of measured frequencies for a given number of Fringes and Exponent.
16600      !   Bin        Used to index Mantissa(*).
16610      !   Min         Used as a subrange index for Table(*).
16620      !   Max         Used as a subrange index for Table(*).
16630      !   OPTION BASE 1
16640      !   REAL Mantissa(0:1023),Time(0:1023),Freq(0:1023)
16650      !   If the last entry in the table is not zero then the table has already been created.
16660      !   IF Table(32766) THEN SUBEXIT
16670      !   FOR Bin=0 TO 1023      ! Fill Mantissa array.
16680      !       Mantissa(Bin)=Bin
16690      !   NEXT Bin
16700      !   Mantissa(0)=1
16710      !   Min=0
16720      !   FOR Fringes=0 TO 1      ! 0 indicates 16 fringes while 1 indicates 8 fringes.
16730      !       FOR Exponent=0 TO 15
16740      !           Max=Min+1023
16750      !           IF Max=32767 THEN      ! Maximum size of an array is 32766.
16760      !               Max=32766
16770      !               REDIM Mantissa(0:1022),Time(0:1022),Freq(0:1022)
16780      !           END IF
16790      !           DISP Fringes,Exponent
16800      !           !MAT Time= Mantissa*(2^(Exponent-1)/5000000000)      ! Use this line with new macrodynes.
16810      !           MAT Time= Mantissa*(2^(Exponent-3)/5000000000)      ! Use this line with old macrodynes.
16820      !           MAT Freq= (2^(4-Fringes))/Time
16830      !           MAT Freq= Freq/(1000000)

```



```

16810             MAT Table(Min:Max)= Freq
16820             Min=Min+1024
16830             NEXT Exponent
16840             NEXT Fringes
16850         SUBEND
16860 Lvdas:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
16870 Lvdas_init: SUB Lvdas_init(@Lvdas)
16880             ! Description:
16890             !       This subprogram is used to initialize the HP98622-66501 Rev B 16-bit General Purpose
16900             !       Input Output (GPIO) interface. The subprogram also opens the LVDAS path on the HP computer
16910             !       for command and data transfer. The I/O path is given the name "@Lvdas". Data transferred
16920             !       from the HP to the LVDAS will use the "OUTPUT @Lvdas" statement. Data transferred to the HP
16930             !       from LVDAS will use the "ENTER @Lvdas" statement.
16940             !       The I/O path has a select code of 12 and is initialized to perform unformatted word
16950             !       transfers without any end of line designations. The DIP switches on the HP98622-66501 Rev B
16960             !       printed circuit board need to be set as shown below:
16970             !       DIP switches for INT LVL      : Bit1=0   Bit0=0
16980             !       DIP switches for Select Code : Bit4=0   Bit3=1   Bit2=1   Bit1=0   Bit0=0
16990             !       DIP switches for DI15to08 clk: RDY =1   BSY =0   RD  =1
17000             !       DIP switches for DI07to00 clk: RDY =1   BSY =0   RD  =1
17010             !       DIP switches for Hndsk Levels: DOUT=0   DIN  =0   HSHK=1   PSTS=0   PFLG=0   PCTL=1
17020             ASSIGN @Lvdas TO 12;WORD,FORMAT OFF,EOL ""
17030             OUTPUT @Lvdas USING "#,AA","HP"
17040         SUBEND
17050 Lvdas_take: SUB Lvdas_take(@Lvdas,Atime,Ctime,INTEGER At_exp,Ct_exp,Cmask,Nsam)
17060             ! Description:
17070             !       This subprogram samples the two analog, three digital, and two external trigger channels
17080             !       from the LVDAS. The HP sends a "CS" to sample the LVDAS data with coincidence. Following the
17090             !       "CS" the HP sends the LVDAS an additional eight words to specify the acquisition and
17100             !       coincidence times, the interarrival and coincidence time exponents, the coincidence mask, and
17110             !       the number of desired samples. After the desired number of samples is acquired or the desired
17120             !       acquisition time expires then the LVDAS sends to the HP an updated number of samples (Nsam).
17130             !       The updated Nsam may be less than the original Nsam if the desired acquisition time expires
17140             !       before the desired Nsam samples are realized.
17150             ! Variables:
17160             !       Atime  The maximum desired acquisition time (seconds).
17170             !       Ctime  The maximum desired coincidence time (seconds).
17180             !       At1    The upper word of integer of 10000000*Atime.
17190             !       At2    The lower word of integer of 10000000*Atime.
17200             !       Ct1    The upper word of integer of 10000000*Ctime.
17210             !       Ct2    The lower word of integer of 10000000*Ctime.
17220             !       At_exp Exponent for interarrival times.
17230             !       Ct_exp Exponent for coincidence times.
17240             !       Nsam   Number of desired samples.
17250             !       Cmask  Coincidence Mask for U,V,W selection.
17260             !       Raw(*) Array of raw data acquired LVDAS data.
17270             OPTION BASE 1
17280             COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
17290             INTEGER At1,At2,Ct1,Ct2
17300             DISP "Taking Data"
17310             CALL Convert2words(Atime*10000000,At1,At2)
17320             CALL Convert2words(Ctime*10000000,Ct1,Ct2)
17330             OUTPUT @Lvdas USING "AA,8(W)";"CS",At1,At2,Ct1,Ct2,At_exp,Ct_exp,Cmask,Nsam
17340             ENTER @Lvdas USING "#,W";Nsam
17350             IF Nsam=0 THEN SUBEXIT
17360             REDIM Raw(1:Nsam,1:10)
17370             ENTER @Lvdas USING "#,W";Raw(*)
17380         SUBEND
17390 Data:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
17400 Data_reduce: SUB Data_reduce(INTEGER At_exp,Ct_exp,Nsam)
17410             ! Description:
17420             !       This subprogram separates the ten by Nsam Raw(*) data array into multiple one by Nsam
17430             !       arrays. The frequency arrays U1,V1,W1 are extracted from columns 6,7,8 of the Raw data array.
17440             !       The voltage arrays A1,B1 are extracted from columns 9,10 of the Raw data array. The
17450             !       interarrival time array I1 is extracted from columns 1 of the Raw data array. The coincidence
17460             !       time array C1 is extracted from columns 2 of the Raw data array. The validation word array
17470             !       Valid is extracted from columns 5 of the Raw data array. If i'th sample acquired contains
17480             !       valid data, then Valid(i) will be equal to one, and zero otherwise. All values for the Valid
17490             !       array are initially set to one by the LVDAS.
17500             !       The raw data from arrays U1,V1,W1 are converted into frequencies by using their initial
17510             !       values as indexes to the frequency look up table array Table(*). The raw data from arrays
17520             !       A1,B1 are converted into voltages by multiplying their initial values by 5 volts over 2^15.
17530             !       The raw data from array I1 are converted into interarrival times by multiplying their initial
17540             !       values by 2^At_exp over 10 to get us. The raw data from array C1 are converted into
17550             !       coincidence times by multiplying their initial values by 2^Ct_exp over 10 to get us.
17560             ! Variables:
17570             !       Table(*) Lookup table of frequencies.
17580             !       Raw(*)   Array of raw data acquired LVDAS data.
17590             !       U1(*)   Array of extracted raw U frequency data.
17600             !       V1(*)   Array of extracted raw V frequency data.

```

```

17610      !      Wi(*)      Array of extracted raw W frequency data.
17620      !      Ai(*)      Array of extracted raw A voltage data.
17630      !      Bi(*)      Array of extracted raw B voltage data.
17640      !      Ii(*)      Array of extracted raw interarrival time data.
17650      !      Ci(*)      Array of extracted raw coincidence time data.
17660      !      Valid(*)   Array of extracted raw validation words.
17670      !      At_exp     Exponent of interarrival times.
17680      !      Ct_exp     Exponent of coincidence times.
17690      !      Nsam       Number of samples acquired.
17700      OPTION BASE 1
17710      COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),Ui(*),Vi(*),Wi(*),Ai(*),Bi(*),Ii(*),Ci(*)
17720      REDIM Ui(Nsam),Vi(Nsam),Wi(Nsam),Ai(Nsam),Bi(Nsam),Ii(Nsam),Ci(Nsam),Valid(Nsam)
17730      DISP "Reducing Data"
17740      MAT Ii= Raw(*,1)
17750      MAT Ci= Raw(*,2)
17760      MAT Valid= Raw(*,5)
17770      MAT Ui= Raw(*,6)
17780      MAT Vi= Raw(*,7)
17790      MAT Wi= Raw(*,8)
17800      MAT Ai= Raw(*,9)
17810      MAT Bi= Raw(*,10)
17820      FOR K=1 TO Nsam
17830          Ui(K)=Table(Ui(K))
17840          Vi(K)=Table(Vi(K))
17850          Wi(K)=Table(Wi(K))
17860      NEXT K
17870      MAT Ai= Ai*(5/32768)
17880      MAT Bi= Bi*(5/32768)
17890      MAT Ii= Ii*(2^At_exp/10)
17900      MAT Ci= Ci*(2^Ct_exp/10)
17910      MAT Ui= Ui . Valid
17920      MAT Vi= Vi . Valid
17930      MAT Wi= Wi . Valid
17940      MAT Ai= Ai . Valid
17950      MAT Bi= Bi . Valid
17960      MAT Ii= Ii . Valid
17970      MAT Ci= Ci . Valid
17980      SUBEND
17990 Data_clip: SUB Data_clip(INTEGER Nsam,REAL Umin,Umax,Vmin,Vmax,Wmin,Wmax)
18000      !      Description:
18010      !      This subprogram compares each of the instantaneous U,V, and W frequencies with user
18020      !      selectable minimum and maximum frequencies. If the instantaneous value is less than the
18030      !      desired minimum, then the validation word is set to zero. Also, if the instantaneous value is
18040      !      greater than the desired maximum, then the validation word is set to zero. The setting of the
18050      !      validation words to zero will have the net effect of discarding the data samples from the data
18060      !      set. In other words, the data is weighted as zero for the average, sdv, normal and shear
18070      !      stress calculations.
18080      !      Variables:
18090      !      Nsam       Number of samples acquired.
18100      !      Ui(*)      Array of instantaneous U frequencies (MHz).
18110      !      Vi(*)      Array of instantaneous V frequencies (MHz).
18120      !      Wi(*)      Array of instantaneous W frequencies (MHz).
18130      !      Valid(*)   Array of sample validation words.
18140      !      Umin       The minimum acceptable U frequency (MHz).
18150      !      Umax       The maximum acceptable U frequency (MHz).
18160      !      Vmin       The minimum acceptable V frequency (MHz).
18170      !      Vmax       The maximum acceptable V frequency (MHz).
18180      !      Wmin       The minimum acceptable W frequency (MHz).
18190      !      Wmax       The maximum acceptable W frequency (MHz).
18200      OPTION BASE 1
18210      COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),Ui(*),Vi(*),Wi(*),Ai(*),Bi(*),Ii(*),Ci(*)
18220      DISP "Clipping Histograms"
18230      FOR K=1 TO Nsam
18240          MAT SEARCH Ui(*),LOC(<Umin);L,K
18250          IF L<Nsam THEN Valid(L)=0
18260          K=L
18270      NEXT K
18280      FOR K=1 TO Nsam
18290          MAT SEARCH Vi(*),LOC(>Umax);L,K
18300          IF L<Nsam THEN Valid(L)=0
18310          K=L
18320      NEXT K
18330      FOR K=1 TO Nsam
18340          MAT SEARCH Vi(*),LOC(<Vmin);L,K
18350          IF L<Nsam THEN Valid(L)=0
18360          K=L
18370      NEXT K
18380      FOR K=1 TO Nsam
18390          MAT SEARCH Vi(*),LOC(>Vmax);L,K
18400          IF L<Nsam THEN Valid(L)=0

```

```

18410         K=L
18420     NEXT K
18430     FOR K=1 TO Nsam
18440         MAT SEARCH Wl(*),LOC(<Wmin);L,K
18450         IF L<Nsam THEN Valid(L)=0
18460         K=L
18470     NEXT K
18480     FOR K=1 TO Nsam
18490         MAT SEARCH Wl(*),LOC(>Wmax);L,K
18500         IF L<Nsam THEN Valid(L)=0
18510         K=L
18520     NEXT K
18530     MAT Ui= Ui . Valid
18540     MAT Vi= Vi . Valid
18550     MAT Wi= Wi . Valid
18560     MAT Ai= Ai . Valid
18570     MAT Bi= Bi . Valid
18580     MAT Ii= Ii . Valid
18590     MAT Ci= Ci . Valid
18600 SUBEND
18610 Data_fconvert: SUB Data_fconvert(Array(*))
18620     ! Description:
18630     ! This subprogram takes the frequency values from the arrays Ui,Vi,Wi and replaces them with
18640     ! velocities after doing the frequency to velocity conversion.
18650     ! Variables:
18660     ! Array(*) An array containing relevant LDV laser and tunnel condition parameters
18670     ! Frng_spc(*) Fringe Spacings extracted from Array(*).
18680     ! Brg_frq(*) Bragg Frequencies extracted from Array(*).
18690     ! Mix_frq(*) Mixing Freqs. extracted from Array(*).
18700     ! Mea_sgn(*) Measured Freq's. Signs extracted from Array(*)
18710     ! Brg_sgn(*) Bragg Freq's. Signs extracted from Array(*).
18720     ! Mix_sgn(*) Mixing Freq's. Signs extracted from Array(*).
18730     ! Ui(*) Array of instantaneous U data.
18740     ! Vi(*) Array of instantaneous V data.
18750     ! Wi(*) Array of instantaneous W data.
18760     ! Equations:
18770     ! The following equations are used to convert the frequencies to velocities
18780     ! Velocity = Fs * Ftotal
18790     ! Ftotal = MeaSgn*Fmeas+BrgSgn*Fbrag+MixSgn*Fmix
18800     OPTION BASE 1
18810     COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),Ui(*),Vi(*),Wi(*),Ai(*),Bi(*),Ii(*),Ci(*)
18820     DIM Frng_spc(3),Brg_frq(3),Mix_frq(3),Mea_sgn(3),Brg_sgn(3),Mix_sgn(3)
18830     DISP "Converting Data"
18840     MAT Frng_spc= Array(35,1:3)
18850     MAT Brg_frq= Array(36,1:3)
18860     MAT Mix_frq= Array(37,1:3)
18870     MAT Mea_sgn= Array(38,1:3)
18880     MAT Brg_sgn= Array(39,1:3)
18890     MAT Mix_sgn= Array(40,1:3)
18900     MAT Ui= Ui*(Mea_sgn(1))
18910     MAT Vi= Vi*(Mea_sgn(2))
18920     MAT Wi= Wi*(Mea_sgn(3))
18930     MAT Ui= Ui+(Brg_sgn(1)*Brg_frq(1)+Mix_sgn(1)*Mix_frq(1))
18940     MAT Vi= Vi+(Brg_sgn(2)*Brg_frq(2)+Mix_sgn(2)*Mix_frq(2))
18950     MAT Wi= Wi+(Brg_sgn(3)*Brg_frq(3)+Mix_sgn(3)*Mix_frq(3))
18960     MAT Ui= Ui*(Frng_spc(1))
18970     MAT Vi= Vi*(Frng_spc(2))
18980     MAT Wi= Wi*(Frng_spc(3))
18990 SUBEND
19000 Data_sum: SUB Data_sum(Sum(*),INTEGER N(*),Nsam)
19010     ! Description:
19020     ! This subprogram performs the summations on the instantaneous LDV and analog data. Data
19030     ! will be weighted as zero in the summations if the value of the validation word is set to zero.
19040     ! Intermediate arrays will be made so that summations of the products of the LDV and analog data
19050     ! can be determined.
19060     ! Variables:
19070     ! Nsam Number of samples acquired.
19080     ! Valid(*) Array of sample validation words.
19090     ! Ui(*) Array of instantaneous U frequency or velocity samples.
19100     ! Vi(*) Array of instantaneous V frequency or velocity samples.
19110     ! Wi(*) Array of instantaneous W frequency or velocity samples.
19120     ! Ai(*) Array of instantaneous A voltage samples.
19130     ! Bi(*) Array of instantaneous B voltage samples.
19140     ! Ii(*) Array of interarrival times.
19150     ! Ci(*) Array of coincidence times.
19160     ! Uu(*) Instantaneous product of the instantaneous Ui & Ui.
19170     ! Vv(*) Instantaneous product of the instantaneous Vi & Vi.
19180     ! Ww(*) Instantaneous product of the instantaneous Wi & Wi.
19190     ! Aa(*) Instantaneous product of the instantaneous Ai & Ai.
19200     ! Bb(*) Instantaneous product of the instantaneous Bi & Bi.

```

```

19210      !      I2(*)      Instantaneous product of the instantaneous Ii & Ii.
19220      !      C2(*)      Instantaneous product of the instantaneous Ci & Ci.
19230      !      Uv(*)      Instantaneous product of the instantaneous Ui & Vi.
19240      !      Vw(*)      Instantaneous product of the instantaneous Vi & Wi.
19250      !      Wu(*)      Instantaneous product of the instantaneous Wi & Ui.
19260      !      Ab(*)      Instantaneous product of the instantaneous Ai & Bi.
19270      !      Ua(*)      Instantaneous product of the instantaneous Ui & Ai.
19280      !      Va(*)      Instantaneous product of the instantaneous Vi & Ai.
19290      !      Wa(*)      Instantaneous product of the instantaneous Wi & Ai.
19300      !      Sum(1,1)    Summation of the array Ui.
19310      !      Sum(2,1)    Summation of the array Vi.
19320      !      Sum(3,1)    Summation of the array Wi.
19330      !      Sum(4,1)    Summation of the array Ai.
19340      !      Sum(5,1)    Summation of the array Bi.
19350      !      Sum(6,1)    Summation of the array Ii.
19360      !      Sum(7,1)    Summation of the array Ci.
19370      !      Sum(1,2)    Summation of the array Uu.
19380      !      Sum(2,2)    Summation of the array Vv.
19390      !      Sum(3,2)    Summation of the array Ww.
19400      !      Sum(4,2)    Summation of the array Aa.
19410      !      Sum(5,2)    Summation of the array Bb.
19420      !      Sum(6,2)    Summation of the array I2.
19430      !      Sum(7,2)    Summation of the array C2.
19440      !      Sum(1,3)    Summation of the array Uv.
19450      !      Sum(2,3)    Summation of the array Vw.
19460      !      Sum(3,3)    Summation of the array Wu.
19470      !      Sum(4,3)    Summation of the array Ab.
19480      !      Sum(5,3)    Summation of the array Ua.
19490      !      Sum(6,3)    Summation of the array Va.
19500      !      Sum(7,3)    Summation of the array Wa.
19510      !      OPTION BASE 1
19520      !      COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
19530      !      REAL Uu(1000),Vv(1000),Ww(1000),Aa(1000),Bb(1000),I2(1000),C2(1000)
19540      !      REAL Uv(1000),Vw(1000),Wu(1000),Ab(1000),Ua(1000),Va(1000),Wa(1000)
19550      !      REDIM Uu(Nsam),Vv(Nsam),Ww(Nsam),Aa(Nsam),Bb(Nsam),I2(Nsam),C2(Nsam)
19560      !      REDIM Uv(Nsam),Vw(Nsam),Wu(Nsam),Ab(Nsam),Ua(Nsam),Va(Nsam),Wa(Nsam)
19570      !      DISP "Summing Data"
19580      !
19590      !      MAT Uu= U1 . U1
19600      !      MAT Vv= V1 . V1
19610      !      MAT Ww= W1 . W1
19620      !      MAT Aa= A1 . A1
19630      !      MAT Bb= B1 . B1
19640      !      MAT Uv= U1 . V1
19650      !      MAT Vw= V1 . W1
19660      !      MAT Wu= W1 . U1
19670      !      MAT Ab= A1 . B1
19680      !      MAT Ua= U1 . A1
19690      !      MAT Va= V1 . A1
19700      !      MAT Wa= W1 . A1
19710      !      MAT I2= I1 . I1
19720      !      MAT C2= C1 . C1
19730      !
19740      !      Sum(1,1)=SUM(U1)
19750      !      Sum(2,1)=SUM(V1)
19760      !      Sum(3,1)=SUM(W1)
19770      !      Sum(4,1)=SUM(A1)
19780      !      Sum(5,1)=SUM(B1)
19790      !      Sum(6,1)=SUM(I1)
19800      !      Sum(7,1)=SUM(C1)
19810      !      Sum(1,2)=SUM(Uu)
19820      !      Sum(2,2)=SUM(Vv)
19830      !      Sum(3,2)=SUM(Ww)
19840      !      Sum(4,2)=SUM(Aa)
19850      !      Sum(5,2)=SUM(Bb)
19860      !      Sum(6,2)=SUM(I2)
19870      !      Sum(7,2)=SUM(C2)
19880      !      Sum(1,3)=SUM(Uv)
19890      !      Sum(2,3)=SUM(Vw)
19900      !      Sum(3,3)=SUM(Wu)
19910      !      Sum(4,3)=SUM(Ab)
19920      !      Sum(5,3)=SUM(Ua)
19930      !      Sum(6,3)=SUM(Va)
19940      !      Sum(7,3)=SUM(Wa)
19950      !      MAT N= (SUM(Valid))
19960      !      SUBEND
19970      !      Data_calc: SUB Data_calc(INTEGER N(*),REAL Sum(*),U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,
19980      !      !      Description:
19990      !      !      This subprogram uses the summations on the instantaneous LDV and analog data as well as the

```

```

20000      !      summations of the products of the LDV and analog data. The subprogram takes these summations
20010      !      and calculates the averages, standard deviations, and shear stresses.
20020      !      Variables:
20030      !      Nsam      The number of valid samples.
20040      !      Sum(1,1)  Summation of the array U1.
20050      !      Sum(2,1)  Summation of the array V1.
20060      !      Sum(3,1)  Summation of the array W1.
20070      !      Sum(4,1)  Summation of the array A1.
20080      !      Sum(5,1)  Summation of the array B1.
20090      !      Sum(6,1)  Summation of the array I1.
20100      !      Sum(7,1)  Summation of the array C1.
20110      !      Sum(1,2)  Summation of the array Uu.
20120      !      Sum(2,2)  Summation of the array Vv.
20130      !      Sum(3,2)  Summation of the array Ww.
20140      !      Sum(4,2)  Summation of the array Aa.
20150      !      Sum(5,2)  Summation of the array Bb.
20160      !      Sum(6,2)  Summation of the array I2.
20170      !      Sum(7,2)  Summation of the array C2.
20180      !      Sum(1,3)  Summation of the array Uv.
20190      !      Sum(2,3)  Summation of the array Vw.
20200      !      Sum(3,3)  Summation of the array Wu.
20210      !      Sum(4,3)  Summation of the array Ab.
20220      !      Sum(5,3)  Summation of the array Ua.
20230      !      Sum(6,3)  Summation of the array Va.
20240      !      Sum(7,3)  Summation of the array Wa.
20250      !      U        Average U frequency or velocity.
20260      !      V        Average V frequency or velocity.
20270      !      W        Average W frequency or velocity.
20280      !      A        Average A voltage.
20290      !      B        Average B voltage.
20300      !      I        Average interarrival time.
20310      !      C        Average coincidence time.
20320      !      U1       Standard deviation for U frequency or velocity.
20330      !      V1       Standard deviation for V frequency or velocity.
20340      !      W1       Standard deviation for W frequency or velocity.
20350      !      A1       Standard deviation for A voltage.
20360      !      B1       Standard deviation for B voltage.
20370      !      I1       Standard deviation for interarrival time.
20380      !      C1       Standard deviation for coincidence time.
20390      !      U1v1     Velocity:Velocity Shear Stress.
20400      !      V1w1     Velocity:Velocity Shear Stress.
20410      !      W1u1     Velocity:Velocity Shear Stress.
20420      !      A1b1     Voltage :Voltage Shear Stress.
20430      !      U1a1     Velocity:Voltage Shear Stress.
20440      !      V1a1     Velocity:Voltage Shear Stress.
20450      !      W1a1     Velocity:Voltage Shear Stress.
20460      DISP "Calculating Results"
20470      Nsam=N(1,1)
20480      IF Nsam>0 THEN
20490          U=Sum(1,1)/Nsam
20500          V=Sum(2,1)/Nsam
20510          W=Sum(3,1)/Nsam
20520          A=Sum(4,1)/Nsam
20530          B=Sum(5,1)/Nsam
20540          I=Sum(6,1)/Nsam
20550          C=Sum(7,1)/Nsam
20560          U1=SQR(ABS(Sum(1,2)/Nsam-U*U))
20570          V1=SQR(ABS(Sum(2,2)/Nsam-V*V))
20580          W1=SQR(ABS(Sum(3,2)/Nsam-W*W))
20590          A1=SQR(ABS(Sum(4,2)/Nsam-A*A))
20600          B1=SQR(ABS(Sum(5,2)/Nsam-B*B))
20610          I1=SQR(ABS(Sum(6,2)/Nsam-I*I))
20620          C1=SQR(ABS(Sum(7,2)/Nsam-C*C))
20630          U1v1=Sum(1,3)/Nsam-U*V
20640          V1w1=Sum(2,3)/Nsam-V*W
20650          W1u1=Sum(3,3)/Nsam-W*U
20660          A1b1=Sum(4,3)/Nsam-A*B
20670          U1a1=Sum(5,3)/Nsam-U*A
20680          V1a1=Sum(6,3)/Nsam-V*A
20690          W1a1=Sum(7,3)/Nsam-W*A
20700      ELSE
20710          U=0
20720          V=0
20730          W=0
20740          A=0
20750          B=0
20760          I=0
20770          C=0
20780          U1=0
20790          V1=0

```

```

20800      W1=0
20810      A1=0
20820      B1=0
20830      I1=0
20840      C1=0
20850      Ulv1=0
20860      V1w1=0
20870      W1u1=0
20880      Albl=0
20890      U1a1=0
20900      V1a1=0
20910      W1a1=0
20920      END IF
20930      SUBEND
20940 Data_trnsfrm: SUB Data_trnsfrm(REAL K(*),U,V,W,U1,V1,W1,Ulv1,V1w1,W1u1,U1a1,V1a1,W1a1)
20950      ! Description:
20960      ! This subprogram performs a coordinate system transformation on the averages, standard
20970      ! deviations, and shear stresses. The coordinate system transformation to be applied is passed
20980      ! through the "K3X3" array. If a LASER to TUNNEL coordinate system transformation is to be
20990      ! performed, then the array "Ldv2tun" array will be passed to the "K3X3" array. If a TUNNEL to
21000      ! MODEL coordinate system transformation is to be performed, then the array "Tun2mod" array will
21010      ! be passed to the "K3X3" array.
21020      ! Variables:
21030      ! U      Average U velocity.
21040      ! V      Average V velocity.
21050      ! W      Average W velocity.
21060      ! U1     Standard deviation for U velocity.
21070      ! V1     Standard deviation for V velocity.
21080      ! W1     Standard deviation for W velocity.
21090      ! Ulul   Velocity:Velocity Normal Stress.
21100      ! Ulv1   Velocity:Velocity Shear Stress.
21110      ! Ulw1   Velocity:Velocity Shear Stress.
21120      ! V1u1   Velocity:Velocity Shear Stress.
21130      ! V1v1   Velocity:Velocity Normal Stress.
21140      ! V1w1   Velocity:Velocity Shear Stress.
21150      ! W1u1   Velocity:Velocity Shear Stress.
21160      ! W1v1   Velocity:Velocity Shear Stress.
21170      ! W1w1   Velocity:Velocity Normal Stress.
21180      ! U1a1   Velocity:Voltage Shear Stress.
21190      ! V1a1   Velocity:Voltage Shear Stress.
21200      ! W1a1   Velocity:Voltage Shear Stress.
21210      ! R(*)   Original U,V,W.
21220      ! F(*)   Original Ulal,Vlal,Wlal.
21230      ! P(*)   Original stress terms Ulul,Ulv1,...,W1w1.
21240      ! K3X3   Coordinate system transformation matrix for average and Velocity:Voltage shear stress
21250      !          conversions.
21260      ! K9X9   Coordinate system transformation matrix for Velocity:Velocity normal and shear stress
21270      !          conversions.
21280      ! S(*)   Transformed U,V,W.
21290      ! H(*)   Transformed Ulal,Vlal,Wlal.
21300      ! Q(*)   Transformed stress terms Ulul,Ulv1,...,W1w1.
21310      OPTION BASE 1
21320      REAL R(3),S(3),F(3),H(3),P(9),Q(9),K3x3(3,3),K9x9(9,9)
21330      DISP "Transforming Results"
21340      ! Calculate Ulul,V1v1,W1w1 using U1,V1,W1.
21350      Ulul=U1*U1
21360      V1v1=V1*V1
21370      W1w1=W1*W1
21380      ! Set Ulw1,V1u1,W1v1 equal to W1u1,Ulv1,V1w1.
21390      Ulw1=W1u1
21400      V1u1=Ulv1
21410      W1v1=V1w1
21420      ! Fill the matrix R with U,V,W.
21430      R(1)=U
21440      R(2)=V
21450      R(3)=W
21460      ! Fill the matrix F with Ulal,Vlal,Wlal.
21470      F(1)=Ulal
21480      F(2)=Vlal
21490      F(3)=Wlal
21500      ! Fill the matrix P with Ulul,Ulv1,Ulw1,V1u1,V1v1,V1w1,W1u1,W1v1,W1w1.
21510      P(1)=Ulul
21520      P(2)=Ulv1
21530      P(3)=Ulw1
21540      P(4)=V1u1
21550      P(5)=V1v1
21560      P(6)=V1w1
21570      P(7)=W1u1
21580      P(8)=W1v1
21590      P(9)=W1w1

```

```

21600      ! Define the matrix K9x9 using products of the elements from then matrix K3x3.
21610      FOR X=1 TO 9
21620          FOR Y=1 TO 9
21630              Y1=((Y-1) DIV 3)+1
21640              X1=((X-1) DIV 3)+1
21650              Y2=((Y-1) MOD 3)+1
21660              X2=((X-1) MOD 3)+1
21670              K9x9(Y,X)=K3x3(Y1,X1)*K3x3(Y2,X2)
21680          NEXT Y
21690      NEXT X
21700      ! Transform matrix R to S using K3x3.
21710      MAT S= K3x3*R
21720      ! Transform matrix F to H using K3x3.
21730      MAT H= K3x3*F
21740      ! Transform matrix P to Q using K9x9.
21750      MAT Q= K9x9*P
21760      ! Extract the transformed U,V,W from the matrix S.
21770      U=S(1)
21780      V=S(2)
21790      W=S(3)
21800      ! Extract the transformed U1a1,V1a1,W1a1 from the matrix H.
21810      U1a1=H(1)
21820      V1a1=H(2)
21830      W1a1=H(3)
21840      ! Extract the transformed U1v1,U1w1,V1v1,V1w1,W1v1,W1w1 from the matrix Q.
21850      U1v1=Q(1)
21860      U1w1=Q(2)
21870      U1v1=Q(3)
21880      V1v1=Q(4)
21890      V1w1=Q(5)
21900      V1w1=Q(6)
21910      W1v1=Q(7)
21920      W1v1=Q(8)
21930      W1w1=Q(9)
21940      ! Calculate U1,V1,W1 using U1v1,V1v1,W1v1.
21950      U1=SQR(ABS(U1v1))
21960      V1=SQR(ABS(V1v1))
21970      W1=SQR(ABS(W1v1))
21980      ! Return transformed U,V,W,U1,V1,W1,U1v1,V1v1,W1v1,U1a1,V1a1,W1a1 to main program.
21990      SUBEND
22000      Print:
22010      Data_print: SUB Data_print(Axis,Pos,INTEGER Nsam,CS,REAL U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1v1,W1v1,A1b1,
      U1a1,V1a1,W1a1,Uedge)
22020      ! Description:
22030      ! This subprogram prints the averages, standard deviations, and shear & normal stress in
22040      ! tabular form. This subprogram may be called several times. The first call might print the
22050      ! reduced velocity data when their units are in frequency (MHz). Subsequent calls will print the
22060      ! reduced data when their units are in velocity (m/s). These subsequent calls will print the
22070      ! the data in one of three coordinate systems: LASER, TUNNEL, and MODEL.
22080      ! Variables:
22090      ! U Average U velocity.
22100      ! V Average V velocity.
22110      ! W Average W velocity.
22120      ! A Average W velocity.
22130      ! B Average B voltage.
22140      ! I Average interarrival time.
22150      ! C Average coincidence time.
22160      ! U1 Standard deviation for U velocity.
22170      ! V1 Standard deviation for V velocity.
22180      ! W1 Standard deviation for W velocity.
22190      ! A1 Standard deviation for A voltage.
22200      ! B1 Standard deviation for B voltage.
22210      ! I1 Standard deviation for interarrival time.
22220      ! C1 Standard deviation for coincidence time.
22230      ! U1v1 Velocity:Velocity Shear Stress.
22240      ! V1w1 Velocity:Velocity Shear Stress.
22250      ! W1u1 Velocity:Velocity Shear Stress.
22260      ! A1b1 Voltage :Voltage Shear Stress.
22270      ! U1a1 Velocity:Voltage Shear Stress.
22280      ! V1a1 Velocity:Voltage Shear Stress.
22290      ! W1a1 Velocity:Voltage Shear Stress.
22300      ! Axis$ Indicates one of the three axes X,Y,Z being traversed.
22310      ! Pos Current Traverse Position.
22320      ! Nsam Number of samples acquired.
22330      ! CS Indicates units and/or coordinate system of data printed.
22340      DISP "Printing Results"
22350      ON ERROR CALL Error
22360      PRINTER IS PRT;WIDTH 144
22370      Axis$=CHR$(NUM("X")+Axis-1)
22380      PRINT USING 22490;L$,Pos,U,U1,U1v1

```

```

22390 PRINT USING 22530;A,A1,A1b1,U1a1
22400 PRINT USING 22500;"N",Nsam,V,V1,V1w1
22410 PRINT USING 22540;B,B1,I1,V1a1
22420 PRINT USING 22510;C$[1,3],W,W1,W1u1
22430 PRINT USING 22550;C,I,C1,W1a1
22440 PRINT USING 22520;Uedge,U1/Uedge,V1/Uedge
22450 PRINT USING 22560;W1/Uedge,(U1v1)/Uedge^2,(V1w1)/Uedge^2,(W1u1)/Uedge^2
22460 IF C$<>"MOD" THEN PRINT
22470 PRINTER IS CRT
22480 OFF ERROR
22490 IMAGE #,8X,A,"=",3D.4D," U=",3D.5D," U1=",3D.5D," U1v1=",4D.6D
22500 IMAGE #,8X,A,"=",8D," V=",3D.5D," V1=",3D.5D," V1w1=",4D.6D
22510 IMAGE #,8X,3A,7X," W=",3D.5D," W1=",3D.5D," W1u1=",4D.6D
22520 IMAGE #,18X," UE=",3D.4D," U1/UE=",3D.5D," V1/UE=",3D.4D
22530 IMAGE " A =",3D.5D," A1 =",3D.5D," W/UE =",3D.5D," U1a1=",2D.6D
22540 IMAGE " U/UE=",3D.5D," V/UE =",3D.5D," IAT1=",9D," V1a1=",2D.6D
22550 IMAGE " CT=",9D," IAT=",9D," CT1 =",9D," W1a1=",2D.6D
22560 IMAGE " W1/UE=",1D.4D," U1v1/UE2=",2D.4D," V1w1/UE2=",2D.4D," W1u1/UE2=",2D.4D
22570 SUBEND
22580 Data_plot: SUB Data_plot(Array(*),Symbols(*),G,Y,P1,P2,P3,Scale,INTEGER N1,N2,N3)
22590 ! Description:
22600 ! This subprogram plots the averages, standard deviations, or shear & normal stress in
22610 ! the 4 profile plots on the CRT. This subprogram will typically be called 4 times. The first
22620 ! call will plot the average velocities normalized by Uedge. The second call will plot the
22630 ! normalized standard deviations of the velocities. The third call will plot the normalized
22640 ! velocity shear stresses. The forth and last call will plot the average and standard deviations
22650 ! of the analog data in the forth and last plot. Data points outside the plot boundaries will
22660 ! be plotted at the plot boundary.
22670 ! Variables:
22680 ! Array(*) Array containing the plot positions and scales.
22690 ! Symbols(*) Array of Symbol arrays. Each symbol array contains a distinct geometric symbol.
22700 ! G Indicates which plot that the normalized P1,P2,P3 will be plotted against Y in.
22710 ! Y Vertical position of the normalized data points in the plot.
22720 ! P1 Horizontal position of the 1st data point (P1 will be normalized by Scale).
22730 ! P2 Horizontal position of the 2nd data point (P2 will be normalized by Scale).
22740 ! P3 Horizontal position of the 3rd data point (P3 will be normalized by Scale).
22750 ! Scale The value that the horizontal positions will be normalized by.
22760 ! N1 The number of samples contributing to P1's value. P1 will be plotted if N1>0.
22770 ! N2 The number of samples contributing to P2's value. P2 will be plotted if N2>0.
22780 ! N3 The number of samples contributing to P3's value. P3 will be plotted if N3>0.
22790 ! Wndw(*) Array containing the plot's scales.
22800 ! Vwprt(*) Array containing the plot's CRT position.
22810 ! Symbol(*) Array containing a distinct geometric symbol.
22820 OPTION BASE 1
22830 DIM Wndw(4),Vwprt(4),Symbol(20,3)
22840 DISP "Plotting Results"
22850 AREA PEN -1
22860 PEN 1
22870 MAT Wndw= Array(60+G,*)
22880 MAT Vwprt= Array(70+G,*)
22890 VIEWPORT Vwprt(1)/10.23,Vwprt(2)/10.23,Vwprt(3)/10.23,Vwprt(4)/10.23
22900 WINDOW Wndw(1),Wndw(2),Wndw(3),Wndw(4)
22910 CLIP ON
22920 FOR I=0 TO 2
22930 IF I=0 AND N1=0 THEN 23120
22940 IF I=1 AND N2=0 THEN 23120
22950 IF I=2 AND N3=0 THEN 23120
22960 Sy=I+1
22970 Noc=Symbols(Sy,0,1)
22980 REDIM Symbol(Noc,3)
22990 MAT Symbol= Symbols(Sy,1:Noc,*)
23000 SELECT I
23010 CASE 0
23020 X=P1*Scale
23030 CASE 1
23040 X=P2*Scale
23050 CASE 2
23060 X=P3*Scale
23070 END SELECT
23080 Xm=MIN(MAX(X,Wndw(1)),Wndw(2))
23090 Ym=MIN(MAX(Y,Wndw(3)),Wndw(4))
23100 MOVE Xm,Ym
23110 SYMBOL Symbol(*),FILL,EDGE
23120 NEXT I
23130 SUBEND
23140 Tcs8:
23150 Tcs8init: SUB Tcs8init(@Tcs8)
23160 ! Description:
23170 ! This subprogram is used to initialize this computer's internal RS232 serial interface.
23180 ! The subprogram also opens the TCS8 path on the Hewlett Packard series 9000 model 3XX computer

```



```

23190      !      for command and data transfer. The I/O path is given the name "@Tcs8". Data transferred
23200      !      from the HP to the TCS8 will use the "OUTPUT @Tcs8" statement. Data transferred to the HP
23210      !      from TCS8 will use the "ENTER @Tcs8" statement.
23220      !      The I/O path has a select code of 9 and is initialized to perform unformatted byte
23230      !      transfers without any end of line designations.
23240      ASSIGN @Tcs8 TO 9;BYTE,FORMAT OFF,EOL ""
23250      CONTROL 9,0;1      ! Reset interface.
23260      CONTROL 9,3;9600    ! Select a baud rate of 9600.
23270      CONTROL 9,4;31     ! Select even parity, enable parity, 2 stop bits, 8 bits per character.
23280      CONTROL 9,12;IVAL("EF",16) ! Enable Carrier Detect. Disable Data Set Ready. Disable Clear To Send.
23290      CONTROL 9,13;9600    ! Default baud rate of 9600.
23300      CONTROL 9,14;31     ! Default character format: Even parity enabled, 2 stop, 8 bits/ char.
23310
23320 Tcs8set:      SUBEND
23330      SUB Tcs8set(Command$,@Tcs8)
23340      !      Description:
23350      !      This subprogram allows the user to view and then set the various initialization parameters
23360      !      of each channel of the TCS8. These parameters are the current position, counts per inch,
23370      !      counts per revolution, motor velocity, motor acceleration, plus and minus limit switches,
23380      !      home switch, and motor stall indication. All of these parameters can be viewed and set except
23390      !      the limit and home switches and the stall indication. They can only be viewed.
23400      !      Variables:
23410      !      Command$      A TCS8 command string which indicates which parameter we want to view & set.
23420      !      View(*)      Array of old TCS8 parameters viewed (received from TCS8). One for each channel.
23430      !      Set(*)      Array of new TCS8 parameters to be set (sent to TCS8). One for each channel.
23440      !      Name$(*)      String array of TCS8 parameter names.
23450      !      Image$(*)      String array of image formats.
23460      !      Units$(*)      String array of units.
23470      !      Channel      Indicates the TCS8 channel number. Used to index the above arrays.
23480      OPTION BASE 1
23490      DIM View(8,1),Set(8,2),Name$(8,1)[10],Image$(8,1)[10],Units$(8,1)[10]
23500      OUTPUT @Tcs8 USING "K,/","V"&Command$&"0"      ! Tell the TCS8 we want to View a parameter.
23510      ENTER @Tcs8 USING "8(K)";View(*)      ! Enter the parameter specified by Command$.
23520      ! Initialize the Name$,Image$,Units$ and Set arrays.
23530      READ Name$(*)
23540      MAT Image$= ("6D,3D")
23550      DATA X1,X2,Y1,Y2,Z1,Z2,A1,A2
23560      FOR Channel=1 TO 8
23570          Set(Channel,1)=Channel
23580          SELECT Command$
23590              CASE "P"      ! Command$="P" indicates we want to view the encoder Positions in inches.
23600                  Name$(Channel,1)=Name$(Channel,1)&" (pos)"
23610                  Units$(Channel,1)="in"
23620              CASE "U"      ! Command$="U" indicates we want to view the Units in counts per inch.
23630                  Name$(Channel,1)=Name$(Channel,1)&" (cpi)"
23640                  Units$(Channel,1)="cnt"
23650              CASE "R"      ! Command$="R" indicates we want to view the number counts per Revolution.
23660                  Name$(Channel,1)=Name$(Channel,1)&" (cpr)"
23670                  Units$(Channel,1)="cnt"
23680              CASE "V"      ! Command$="V" indicates we want to view the Velocity in revolution per second.
23690                  Name$(Channel,1)=Name$(Channel,1)&" (vel)"
23700                  Units$(Channel,1)="rev"
23710              CASE "A"      ! Command$="A" indicates we want to view the Acceleration in revolution per second^2.
23720                  Name$(Channel,1)=Name$(Channel,1)&" (acc)"
23730                  Units$(Channel,1)="rev"
23740              CASE "+"      ! Command$="+" indicates we want to view the current + direction limit switches.
23750                  Name$(Channel,1)=Name$(Channel,1)&" (+LS)"
23760                  Units$(Channel,1)=" "
23770              CASE "-"      ! Command$="-" indicates we want to view the current - direction limit switches.
23780                  Name$(Channel,1)=Name$(Channel,1)&" (-LS)"
23790                  Units$(Channel,1)=" "
23800              CASE "S"      ! Command$="S" indicates we want to view the current motor Stall indication status.
23810                  Name$(Channel,1)=Name$(Channel,1)&" (STALL)"
23820                  Units$(Channel,1)=" "
23830              CASE "H"      ! Command$="H" indicates we want to view the current Home limit switches.
23840                  Name$(Channel,1)=Name$(Channel,1)&" (HS)"
23850                  Units$(Channel,1)=" "
23860              END SELECT
23870          NEXT Channel
23880      ! The "Change" subprogram allows the user to see and then change the values of the viewed parameters
23890      CALL Change("VALUES",View(*),Name$(*),Image$(*),Units$(*))
23900      ! The "Set" parameters command is now sent to the TCS8
23910      SELECT Command$
23920          CASE "P","U","R","V","A"
23930              MAT Set(*,2)= View(*,1)
23940              OUTPUT @Tcs8 USING 23940;"S"&Command$,Set(*)
23950              IMAGE K,8(D,";",M6D,4D,""),/
23960          END SELECT
23970      SUBEND
23980 Tcs8read:      SUB Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
23990      !      Description:

```

```

23990      !      This subprogram reads the current TCS8 positions. The 8 positions are read in TCS
24000      !      coordinates with the units being in inches. Four of the eight positions (X1,Y1,Z1,A1) which
24010      !      are the transmitting side traverse positions are entered into the Tcs1 array. The other four
24020      !      positions (X2,Y2,Z2,A2) which are the receiving side traverse positions are entered into the
24030      !      Tcs2 array. The Tcs1 & Tcs2 arrays are converted from TCS to TUNNEL to MODEL coordinates.
24040      !      The current updated positions in the three coordinate systems are printed on the top of the
24050      !      CRT. They are also returned to the main program.
24060      !      Variables:
24070      !      Tcs1(*)      TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TCS coordinates.
24080      !      Tcs2(*)      TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TCS coordinates.
24090      !      Tun(*)       Traverse positions (X,Y,Z) in TUNNEL coordinates.
24100      !      Mod(*)       Traverse positions (X,Y,Z) in MODEL coordinates.
24110      !      Tcs2tun1(*) Coordinate system transformation matrix for converting Tcs1(*) to Tun(*).
24120      !      Tcs2tun2(*) Coordinate system transformation matrix for converting Tcs2(*) to Tun(*).
24130      !      Tun2mod(*)   Coordinate system transformation matrix for converting Tun(*) to Mod(*).
24140      OUTPUT @Tcs8 USING "K, /", "VP0"
24150      ENTER @Tcs8 USING "8(K)"; Tcs1(1), Tcs2(1), Tcs1(2), Tcs2(2), Tcs1(3), Tcs2(3), Tcs1(4), Tcs2(4)
24160      MAT Tun= Tcs2tun1*Tcs1
24170      REDIM Tun(1:3), Mod(1:3)
24180      MAT Mod= Tun2mod*Tun
24190      REDIM Tun(1:4), Mod(1:4)
24200      Mod(4)=0
24210      Tun(4)=0
24220      CALL Tcs8print (Mod(*), Tun(*), Tcs1(*), Tcs2(*))
24230      SUBEND
24240 Tcs8print: SUB Tcs8print (Mod(*), Tun(*), Tcs1(*), Tcs2(*))
24250      !      Description:
24260      !      This subprogram prints the current updated TCS8 positions at the top of the CRT. The
24270      !      positions are printed in TCS , TUNNEL , and MODEL coordinates.
24280      !      Variables:
24290      !      Tcs1(*)      TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TCS coordinates.
24300      !      Tcs2(*)      TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TCS coordinates.
24310      !      Tun(*)       Traverse positions (X,Y,Z) in TUNNEL coordinates.
24320      !      Mod(*)       Traverse positions (X,Y,Z) in MODEL coordinates.
24330      PRINT CHR$(128);
24340      PRINT TABXY(50,1); "
24350      PRINT TABXY(50,2); "      MOD      TUN      TCS1      TCS2      "
24360      PRINT TABXY(50,3); "
24370      PRINT TABXY(50,4);
24380      PRINT USING "#,K,4(M3D.4D),X"; " X:", Mod(1), Tun(1), Tcs1(1), Tcs2(1)
24390      PRINT TABXY(50,5);
24400      PRINT USING "#,K,4(M3D.4D),X"; " Y:", Mod(2), Tun(2), Tcs1(2), Tcs2(2)
24410      PRINT TABXY(50,6);
24420      PRINT USING "#,K,4(M3D.4D),X"; " Z:", Mod(3), Tun(3), Tcs1(3), Tcs2(3)
24430      PRINT TABXY(50,7);
24440      PRINT USING "#,K,4(M3D.4D),X"; " A:", Mod(4), Tun(4), Tcs1(4), Tcs2(4)
24450      PRINT TABXY(50,8); "
24460      SUBEND
24470 Tcs8move: SUB Tcs8move (@Tcs8, Mod(*), Tun(*), Tcs1(*), Tcs2(*), Mod2tun(*), Tun2tcs1(*), Tun2tcs2(*), Side$, Coor$, Mode$,
      K, Movement)
24480      !      Description:
24490      !      This subprogram allows for the movement of the probe volume and collecting optics in one of
24500      !      three coordinate systems. The three coordinate systems implemented are the TSC, TUNNEL and
24510      !      MODEL coordinate systems. Two movements modes are available. The first movement mode makes
24520      !      moves relative to the current position. The second movement mode makes moves to an absolute
24530      !      fixed position. Both the transmitting side and receiving side traverses can be moved in tandem
24540      !      or separately.
24550      !      Variables:
24560      !      Tcs1(*)      TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TCS coordinates.
24570      !      Tcs2(*)      TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TCS coordinates.
24580      !      Tun(*)       Traverse positions (X,Y,Z) in TUNNEL coordinates.
24590      !      Mod(*)       Traverse positions (X,Y,Z) in MODEL coordinates.
24600      !      Mod2tun(*)   Coordinate system transformation matrix for converting Tcs1(*) to Tun(*).
24610      !      Tun2tcs1(*)  Coordinate system transformation matrix for converting Tcs2(*) to Tun(*).
24620      !      Tun2tcs2(*)  Coordinate system transformation matrix for converting Tcs2(*) to Tun(*).
24630      !      Side$       Indicates which sides are to be moved:
24640      !                  Tx      : Transmitting side only.
24650      !                  Rx      : Receiving side only.
24660      !                  Tx & Rx : Both sides together.
24670      !      Coor$       Indicates which coordinate system the movement is to be made:
24680      !                  MODEL   : MODEL coordinates.
24690      !                  TUNNEL  : TUNNEL coordinates.
24700      !                  TCS    : TCS coordinates.
24710      !      Mode$       Indicates which movement mode is to be completed:
24720      !                  RELATIVE: Movements are relative to current positions.
24730      !                  ABSOLUTE: Movements are to absolute positions.
24740      !      K           Indicates which axis of the four axes is to be moved.
24750      !      Movement    Indicates the desired movement for the selected axis.
24760      !      I(*)         Array of viewed TCS8 "Initialized" parameters.
24770      !      C(*)         Array of viewed TCS8 "Currents On" parameters.

```

```

24780      OPTION BASE 1
24790      DIM L$(100)
24800      REAL Move(8,2),I(8),C(8)
24810      ! If all of the channels have not yet been initialized, then do so now.
24820      OUTPUT @Tcs8 USING "K,/","VIO"
24830      ENTER @Tcs8 USING "8(K)";I(*)
24840      IF SUM(I)<>8 THEN OUTPUT @Tcs8 USING "K,/","SIO"
24850      ! If all of the channels do not have their currents turned on, then do so now.
24860      OUTPUT @Tcs8 USING "K,/","VCO"
24870      ENTER @Tcs8 USING "8(K)";C(*)
24880      IF SUM(C)<>8 THEN OUTPUT @Tcs8 USING "K,/","SC0:1,"
24890      ! If the movement mode is to be RELATIVE, then clear all of the previously read positions.
24900      IF Mode$="RELATIVE" THEN
24910          MAT Mod= (0)
24920          MAT Tun= (0)
24930          MAT Tcs1= (0)
24940          MAT Tcs2= (0)
24950      END IF
24960      ! Set the new Tcs1(*) and Tcs2(*) position arrays.
24970      SELECT Coor$
24980      CASE "MODEL"
24990          Mod(K)=Movement
25000          REDIM Tun(1:3),Mod(1:3)
25010          MAT Tun= Mod2tun*Mod
25020          REDIM Tun(1:4),Mod(1:4)
25030          IF POS(Side$,"Tx") THEN MAT Tcs1= Tun2tcs1*Tun
25040          IF POS(Side$,"Rx") THEN MAT Tcs2= Tun2tcs2*Tun
25050      CASE "TUNNEL"
25060          Tun(K)=Movement
25070          IF POS(Side$,"Tx") THEN MAT Tcs1= Tun2tcs1*Tun
25080          IF POS(Side$,"Rx") THEN MAT Tcs2= Tun2tcs2*Tun
25090      CASE "TCS"
25100          IF POS(Side$,"Tx") THEN Tcs1(K)=Movement
25110          IF POS(Side$,"Rx") THEN Tcs2(K)=Movement
25120      END SELECT
25130      ! File the move array.
25140      FOR Channel=1 TO 8
25150          Move(Channel,1)=Channel
25160      NEXT Channel
25170      Move(1,2)=Tcs1(1)
25180      Move(2,2)=Tcs2(1)
25190      Move(3,2)=Tcs1(2)
25200      Move(4,2)=Tcs2(2)
25210      Move(5,2)=Tcs1(3)
25220      Move(6,2)=Tcs2(3)
25230      Move(7,2)=Tcs1(4)
25240      Move(8,2)=Tcs2(4)
25250      ! Initiate the start of the move.
25260      IF Mode$="ABSOLUTE" THEN OUTPUT @Tcs8 USING 25280;"MA",Move(*)
25270      IF Mode$="RELATIVE" THEN OUTPUT @Tcs8 USING 25280;"MR",Move(*)
25280      IMAGE K,8(D,";",S2D,5D,""),/
25290      ! The TCS8 will return the new updated positions only after the move is complete.
25300      ENTER @Tcs8 USING "8(K)";Tcs1(1),Tcs2(1),Tcs1(2),Tcs2(2),Tcs1(3),Tcs2(3),Tcs1(4),Tcs2(4)
25310      ! Turn off the motor drive currents.
25320      OUTPUT @Tcs8 USING "K,/","SC0:0,"
25330      SUBEND
25340      Ctm:
25350      Refract:
25360      SUB Refract(REAL Index(*),Thetal(*),Theta3(*))
25370          ! Description:
25380          ! This subprogram uses the laser beam angles outside the tunnel to compute the angles inside
25390          ! the water tunnel. This requires the knowledge of the indexes of refraction for the various
25400          ! media that the beams go through. The Mediums are air, glass, and water.
25410          ! Variables:
25420          ! Index(*) Array of indexes of refraction.
25430          ! Index(1) : Index of refraction for Air.
25440          ! Index(2) : Index of refraction for Glass.
25450          ! Index(3) : Index of refraction for Water.
25460          ! Thetal(*) Laser beam angles outside the water tunnel.
25470          ! Theta3(*) Laser beam angles inside the water tunnel.
25480      OPTION BASE 1
25490      ! Correct Theta for angles in water.
25500      MAT Theta3= Thetal
25510      IF Index(1)<>Index(3) THEN
25520          Theta3(2,1)=ASN(Index(1)/Index(3)*SIN(Thetal(2,1)))
25530          Theta3(2,2)=ASN(Index(1)/Index(3)*SIN(Thetal(2,2)))+90
25540      END IF
25550      SUBEND
25560      Ctm_ldv:
25570      SUB Ctm_ldv(Theta(*),Tun2ldv(*),Ldv2tun(*))
25580          ! Description:
25590          ! This subprogram computes directly the TUNNEL to LASER coordinate system transformation

```

```

25580      !      matrix "Tun2ldv(*)". However, the desired coordinate system transformation matrix "Ldv2tun" is
25590      !      required. It is the matrix inverse of "Tun2ldv".
25600      !      Variables:
25610      !      Theta(*)      Laser beam angles inside the water tunnel.
25620      !      Tun2ldv(*)    Coordinate system transformation matrix for converting from TUNNEL to LASER.
25630      !      Ldv2tun(*)    Coordinate system transformation matrix for converting from LASER to TUNNEL.
25640      OPTION BASE 1
25650      ! Tun2ldv converts TUNNEL coordinates to LASER coordinates.
25660      Tun2ldv(1,1)=COS(Theta(1,1))
25670      Tun2ldv(1,2)=COS(Theta(1,2))
25680      Tun2ldv(1,3)=COS(Theta(1,3))
25690      Tun2ldv(2,1)=COS(Theta(2,1))
25700      Tun2ldv(2,2)=COS(Theta(2,2))
25710      Tun2ldv(2,3)=COS(Theta(2,3))
25720      Tun2ldv(3,1)=COS(Theta(3,1))
25730      Tun2ldv(3,2)=COS(Theta(3,2))
25740      Tun2ldv(3,3)=COS(Theta(3,3))
25750      ! Ldv2tun converts LASER coordinates to TUNNEL coordinates.
25760      MAT Ldv2tun= INV(Tun2ldv)
25770      SUBEND
25780 Ctm_mod: SUB Ctm_mod(Alpha(*),Mod2tun(*),Tun2mod(*))
25790      ! Description:
25800      !      This subprogram computes directly the MODEL to TUNNEL coordinate system transformation
25810      !      matrix "Mod2tun(*)". However, the desired coordinate system transformation matrix "Tun2mod" is
25820      !      required. It is the matrix inverse of "Mod2tun".
25830      !      Variables:
25840      !      Alpha(*)      Angles of attack, yaw, and roll.
25850      !      T1(*)          Partial coordinate system transformation matrix for converting from MODEL to
25860      !                      TUNNEL coordinates. Takes into account a model at angle of attack.
25870      !      T2(*)          Partial coordinate system transformation matrix for converting from MODEL to
25880      !                      TUNNEL coordinates. Takes into account a model at angle of yaw.
25890      !      T3(*)          Partial coordinate system transformation matrix for converting from MODEL to
25900      !                      TUNNEL coordinates. Takes into account a model at angle of roll.
25910      !      Mod2tun(*)    Coordinate system transformation matrix for converting from MODEL to TUNNEL.
25920      !      Tun2mod(*)    Coordinate system transformation matrix for converting from TUNNEL to MODEL.
25930      OPTION BASE 1
25940      REAL T1(3,3),T2(3,3),T3(3,3),Temp(3,3)
25950      ! Define 1st coordinate transformation matrix for Mod2tun.
25960      ! Rotation in the x-y plane about the z-axis.
25970      ! Used when model is at an angle of attack.
25980      T1(1,1)=COS(Alpha(1))
25990      T1(1,2)=SIN(Alpha(1))
26000      T1(1,3)=0
26010      T1(2,1)=-SIN(Alpha(1))
26020      T1(2,2)=COS(Alpha(1))
26030      T1(2,3)=0
26040      T1(3,1)=0
26050      T1(3,2)=0
26060      T1(3,3)=1
26070      ! Define 2nd coordinate transformation matrix for Mod2tun.
26080      ! Rotation in the x-z plane about the y-axis.
26090      ! Used when model is at an angle of yaw.
26100      T2(1,1)=COS(-Alpha(2))
26110      T2(1,2)=0
26120      T2(1,3)=-SIN(-Alpha(2))
26130      T2(2,1)=0
26140      T2(2,2)=1
26150      T2(2,3)=0
26160      T2(3,1)=SIN(-Alpha(2))
26170      T2(3,2)=0
26180      T2(3,3)=COS(-Alpha(2))
26190      ! Define 3rd coordinate transformation matrix for Mod2tun.
26200      ! Rotation in the y-z plane about the x-axis.
26210      ! Used when model is at an angle of roll.
26220      T3(1,1)=1
26230      T3(1,2)=0
26240      T3(1,3)=0
26250      T3(2,1)=0
26260      T3(2,2)=COS(-Alpha(3))
26270      T3(2,3)=SIN(-Alpha(3))
26280      T3(3,1)=0
26290      T3(3,2)=-SIN(-Alpha(3))
26300      T3(3,3)=COS(-Alpha(3))
26310      ! Mod2tun converts MODEL coordinates to TUNNEL coordinates.
26320      MAT Temp= T2*T1
26330      MAT Mod2tun= T3*Temp
26340      ! Tun2mod converts TUNNEL coordinates to MODEL coordinates.
26350      MAT Tun2mod= INV(Mod2tun)
26360      SUBEND
26370 Ctm_tcs: SUB Ctm_tcs(Tcs2tun1(*),Tcs2tun2(*),Tun2tcs1(*),Tun2tcs2(*),Fs,Fr,Bs,Br,Index(*),Ts,Tr,Ta)

```

```

26380      ! Description:
26390      ! This subprogram computes the TUNNEL to TCS coordinate system transformation matrices
26400      ! "Tun2tcs1(*)" and "Tun2tcs2(*)". The coordinate system transformation matrices "Tcs2tun1" and
26410      ! "Tcs2tun2" are the matrix inverses of "Tun2tcs1(*)" and "Tun2tcs2(*)" respectively.
26420      ! Variables:
26430      ! Tcs2tun1(*) Sending side coordinate transformation matrix converting Tcs(*) to Tun(*).
26440      ! Tun2tcs1(*) Sending side coordinate transformation matrix converting Tun(*) to Tcs(*).
26450      ! Tcs2tun2(*) Receiving side coordinate transformation matrix converting Tcs(*) to Tun(*).
26460      ! Tun2tcs2(*) Receiving side coordinate transformation matrix converting Tun(*) to Tcs(*).
26470      ! Fs Focal length for sending side onaxis and offaxis lenses.
26480      ! Fr Focal length for receiving side offaxis lens.
26490      ! Bs Beam spacings for sending side onaxis and offaxis beam pairs.
26500      ! Br Beam spacing for receiving side offaxis.
26510      ! Index(*) Array of indexes of refraction for air, glass, and water.
26520      ! Ts Angle of offaxis sending side beam pair.
26530      ! Tr Angle of offaxis receiving side beam pair.
26540      ! Ta Sending side offaxis auxiliary traverse angle. Returned to main program.
26550      ! Xs_on,Ys_on Starting coordinates of onaxis sending side lens.
26560      ! Xs_offs,Ys_offs Starting coordinates of offaxis sending side lens.
26570      ! Xs_offr,Ys_offr Starting coordinates of offaxis receiving side lens.
26580      ! Xc,Yc The common point in air of two beam path equations.
26590      ! Ba,Bb The Y intercepts of two beam path equations.
26600      ! Theta(*) Array of angles in which each beam contacts the window.
26610      ! X(*) Array of X coordinates for the points in which each beam contacts the window.
26620      ! Yposition The Y coordinate of the point where all beams cross in the water.
26630      ! Y1,X2,Y2,X3,Y3 Temporary variables to hold the results of the first call to Findstart.
26640      ! Thickness The thickness of the window.
26650      ! Beam Subscript used while determining the X(*) array.
26660      OPTION BASE 1
26670      REAL Xs_on,Ys_on,Xs_offs,Ys_offs,Xs_offr,Ys_offr,Xc,Yc,Ba,Bb,Theta(6),X(6)
26680      REAL Yposition,Y1,X2,Y2,X3,Y3,Thickness
26690      INTEGER Beam
26700      Thickness=1.25
26710      Yposition=0
26720      GOSUB Findstart
26730      Y1=Ys_on
26740      X2=Xs_offs
26750      Y2=Ys_offs
26760      X3=Xs_offr
26770      Y3=Ys_offr
26780      Yposition=1
26790      GOSUB Findstart
26800      MAT Tun2tcs1= IDN
26810      MAT Tun2tcs2= IDN
26820      Tun2tcs1(2,2)=-Ys_on+Y1
26830      Tun2tcs1(4,2)=-SQRT((Xs_offs-X2)^2+(Ys_offs-(Ys_on-Y1+Y2))^2)
26840      Tun2tcs1(4,4)=0
26850      Ta=ATN((Xs_offs-X2)/(Ys_offs-(Ys_on-Y1+Y2)))
26860      Tun2tcs2(1,2)=Xs_offr-X3
26870      Tun2tcs2(2,2)=-Ys_offr+Y3
26880      Tun2tcs2(4,4)=0
26890      MAT Tcs2tun1= INV(Tun2tcs1)
26900      MAT Tcs2tun2= INV(Tun2tcs2)
26910      Tcs2tun1(4,2)=0
26920      Tcs2tun2(4,2)=0
26930      SUBEXIT
26940 Findstart:
26950      ! This subroutine finds the starting coordinates for the onaxis (Xs_on,Ys_on) and offaxis (Xs_offs,
26960      ! Ys_offs) sending side lenses and the offaxis (Xs_offr,Ys_offr) receiving side lens given the point
26970      ! at which all beams cross in the tunnel. The crossing point is given to be (0,Yposition). The
26980      ! method in which the starting coordinates are found involves solving simultaneously the equations for
26990      ! the path of each beam pair in air yielding the common point (Xc,Yc). Given the focal length of the
27000      ! lens, the starting coordinate can be calculated. The equation for each beam path in air is obtained
27010      ! by determining the angle and the point a beam contacts the window.
27020      !
27030      ! These six equations find the six angles.
27040      Theta(1)=ATN(Bs/(2*Fs))
27050      Theta(2)=-ATN(Bs/(2*Fs))
27060      Theta(3)=Ts+ATN(Bs/(2*Fs))
27070      Theta(4)=Ts-ATN(Bs/(2*Fs))
27080      Theta(5)=Tr+ATN(Br/(2*Fr))
27090      Theta(6)=Tr-ATN(Br/(2*Fr))
27100      ! This equation finds the X coordinate of the six points. The Y coordinate is equal to -Thickness.
27110      FOR Beam=1 TO 6
27120          X(Beam)=-Yposition*TAN(ASN(Index(1)/Index(3)*SIN(Theta(Beam))))-
27130              Thickness*TAN(ASN(Index(1)/Index(2)*SIN(Theta(Beam))))
27140      NEXT Beam
27150      ! Determine the Y intercepts for the onaxis beam paths.
27160      Ba=-Thickness-X(1)/TAN(Theta(1))
27170      Bb=-Thickness-X(2)/TAN(Theta(2))
27180      ! Solve for the common point.

```

```

27170      Xc=(Bb-Ba)/(1/TAN(Theta(1))-1/TAN(Theta(2)))
27180      Yc=Xc/TAN(Theta(2))+Bb
27190      ! Calculate the onaxis lens starting coordinate using the focal length and the onaxis angle (=0 deg).
27200      Xs_on=Xc-Fs*SIN(0)
27210      Ys_on=Yc-Fs*COS(0)
27220      ! Determine the Y intercepts for the offaxis sending side beam paths.
27230      Ba=-Thickness-X(3)/TAN(Theta(3))
27240      Bb=-Thickness-X(4)/TAN(Theta(4))
27250      ! Solve for the common point.
27260      Xc=(Bb-Ba)/(1/TAN(Theta(3))-1/TAN(Theta(4)))
27270      Yc=Xc/TAN(Theta(4))+Bb
27280      ! With the focal length and the offaxis angle calculate the starting coordinate
27290      ! of the offaxis sending side lens.
27300      Xs_offs=Xc-Fs*SIN(Ts)
27310      Ys_offs=Yc-Fs*COS(Ts)
27320      ! Determine the Y intercepts for the offaxis receiving side beam paths.
27330      Ba=-Thickness-X(5)/TAN(Theta(5))
27340      Bb=-Thickness-X(6)/TAN(Theta(6))
27350      ! Solve for the common point.
27360      Xc=(Bb-Ba)/(1/TAN(Theta(5))-1/TAN(Theta(6)))
27370      Yc=Xc/TAN(Theta(6))+Bb
27380      ! With the focal length and the offaxis angle calculate the starting coordinate
27390      ! of the offaxis receiving side lens.
27400      Xs_offr=Xc-Fr*SIN(Tr)
27410      Ys_offr=Yc-Fr*COS(Tr)
27420      RETURN
27430  SUBEND
27440 Graph: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
27450 Crt_init: SUB Crt_init
27460      ! Description:
27470      !       This subprogram initializes the CRT as the plotting device and clears both the alpha
27480      !       numerics and graphics part of the CRT.
27490      PLOTTER IS CRT,"INTERNAL"      ! Select the CRT as the plotting device.
27500      AREA PEN 0                    ! Select black for area fills.
27510      PEN 1                        ! Select white for line drawing and labeling.
27520      PRINTER IS CRT               ! Select the CRT as the printing device.
27530      PRINTALL IS CRT              ! Send ERROR and DISP message to CRT.
27540      KEY LABELS OFF              ! Hide the special function key labels for f1..f8.
27550  SUBEND
27560 Read_symbols: SUB Read_symbols(Symbols(*))
27570      ! Description:
27580      !       This subprogram defines 5 geometric symbols to be used with the SYMBOL statement. The
27590      !       symbols provided are as follows: Square,Octagon,Diamond, and Triangles (upwards & downwards
27600      !       pointing triangles). All of the symbols have a dot added to their center.
27610      ! Variables:
27620      !       Symbols(*) Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
27630      !       Symbol(*) Array of coordinates which when connected produce a distinct geometric symbol.
27640      !       Dot(*) Array of coordinates which produce a dot. The dot symbol is added to all symbols.
27650      !       Noc The number of coordinates in a symbol.
27660      !       S Used to index the Symbols array.
27670      OPTION BASE 1
27680      REAL Symbol(20,3),Dot(2,3)
27690      READ Dot(*)
27700      FOR S=1 TO 5
27710          READ Noc
27720          REDIM Symbol(Noc,3)
27730          READ Symbol(*)
27740          MAT Symbols(S,1:Noc,*)= Symbol
27750          MAT Symbols(S,Noc+1:Noc+2,*)= Dot
27760          Symbols(S,0,1)=Noc+2
27770      NEXT S
27780 Dot: DATA 4.5, 7.5,-2, 4.5, 7.5,-1
27790 Square: DATA 5, 0.5, 3.5,-2, 8.5, 3.5,-1, 8.5,11.5,-1, 0.5,11.5,-1, 0.5,3.5,-1
27800 Octagon: DATA 9, 0.5, 5.5,-2, 2.5, 3.5,-1, 6.5, 3.5,-1, 8.5, 5.5,-1, 8.5,9.5,-1, 6.5,11.5,-1,
27810          2.5,11.5,-1, 0.5,9.5,-1, 0.5,5.5,-1
27820 Diamond: DATA 5, -0.5, 7.5,-2, 4.5, 2.5,-1, 9.5, 7.5,-1, 4.5,12.5,-1, -0.5,7.5,-1
27830 Utriangle: DATA 4, 0.5, 4.5,-2, 8.5, 4.5,-1, 4.5,13.5,-1, 0.5, 4.5,-1
27840 Dtriangle: DATA 4, 0.5,10.5,-2, 8.5,10.5,-1, 4.5, 1.5,-1, 0.5,10.5,-1
27850 Setup_graph: SUB Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
27860      ! Description:
27870      !       This subprogram sets up nine empty plots on the CRT screen. Four plots are profile plots
27880      !       while the other five plots are histogram plots. The profile and histogram plots provided are
27890      !       as follows:
27900      !       Graph#      Type      Description
27910      !       1      Histogram #1      U frequency data in MHz.
27920      !       2      Histogram #2      V frequency data in MHz.
27930      !       3      Histogram #3      W frequency data in MHz.
27940      !       4      Histogram #4      Analog Channel #1 data in volts.
27950      !       5      Histogram #5      Analog Channel #2 data in volts.
27960      !       6      Profile Plot #1      Velocity Averages versus Traverse Position.

```

```

27960      !           7      Profile Plot #2      Velocity SDVs versus Traverse Position.
27970      !           8      Profile Plot #3      Velocity Shear Stresses versus Traverse Position.
27980      !           9      Profile Plot #4      Average voltages & SDVs versus Traverse Position.
27990      ! Variables:
28000      !      Array(*)      Array containing the plot positions and scales.
28010      !      Image$(*)      String array containing image formats for the axes labeling.
28020      !      Wndw(*)      Array containing the plot's scales.
28030      !      Vwprt(*)      Array containing the plot's CRT position.
28040      !      Xdiv(*)      Array containing the number of X divisions for the plot's X axis.
28050      !      Ydiv(*)      Array containing the number of Y divisions for the plot's Y axis.
28060      !      Xlabel$(*)      String array containing labels for the X axis.
28070      !      Ylabel$(*)      String array containing labels for the Y axis.
28080      !      Title$(*)      String array containing labels for the Plots.
28090      !      Ximage$(*)      String array containing image formats for the X axis labeling.
28100      !      Yimage$(*)      String array containing image formats for the Y axis labeling.
28110      !      Legend$(*)      String array containing labels for each symbol in a profile plot.
28120      !      Symbols(*)      Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
28130      !      G              Used as an index to the above arrays. Specifies one of nine plots.
28140      !      I              Used as an index to the Legend$ array.
28150      ! OPTION BASE 1
28160      ! COM /Graph/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),
      !      Legend$(*)
28170      ! MAT Wndw= Array(61:69,*)
28180      ! MAT Vwprt= Array(71:79,*)
28190      ! MAT Xdiv(1:5)= Array(81:85,1)
28200      ! MAT Xdiv(6:9)= Array(81:84,3)
28210      ! MAT Ydiv(1:5)= Array(81:85,2)
28220      ! MAT Ydiv(6:9)= Array(81:84,4)
28230      ! MAT Ximage$= Image$(61:69,1)
28240      ! MAT Yimage$= Image$(61:69,3)
28250      ! FOR G=1 TO 9
28260      !     READ G,Xlabel$(G)
28270      !     FOR I=1 TO SIZE(Legend$,2)
28280      !         READ Legend$(G,I)
28290      !     NEXT I
28300      !     SELECT G
28310      !     CASE 1 TO 5
28320      !         Ylabel$(G)=""
28330      !     CASE 6 TO 9
28340      !         Ylabel$(G)=CHR$(NUM("X")+Paxis-1)
28350      !     END SELECT
28360      !     CALL Set_up(G,Symbols(*))
28370      ! NEXT G
28380      ! SUBEXIT
28390      !      G, X axis Label      ,      Symbol #1...5 labels
28400      !      DATA 1, ""          ,      "", "", "", "", ""
28410      !      DATA 2, ""          ,      "", "", "", "", ""
28420      !      DATA 3, ""          ,      "", "", "", "", ""
28430      !      DATA 4, ""          ,      "", "", "", "", ""
28440      !      DATA 5, ""          ,      "", "", "", "", ""
28450      !      DATA 6, "Velocities / Uedge" ,      "U", "V", "W", "", ""
28460      !      DATA 7, "RMS / Uedge"      ,      "U1", "V1", "W1", "", ""
28470      !      DATA 8, "Shear Stress / Uedge^2" ,      "U1v1", "V1w1", "W1u1", "", ""
28480      !      DATA 9, "Fluorescence: C,C1 (volts)" ,      "C", "", "C1", "", ""
28490      ! SUBEND
28500      ! Set_up: SUB Set_up(G,Symbols(*))
28510      !      Description:
28520      !      This subprogram clears and then redraws one of nine empty plots on the CRT screen.
28530      !      Variables:
28540      !      Wndw(*)      Array containing the plot's scales.
28550      !      Vwprt(*)      Array containing the plot's CRT position.
28560      !      Xdiv(*)      Array containing the number of X divisions for the plot's X axis.
28570      !      Ydiv(*)      Array containing the number of Y divisions for the plot's Y axis.
28580      !      Xlabel$(*)      String array containing labels for the X axis.
28590      !      Ylabel$(*)      String array containing labels for the Y axis.
28600      !      Title$(*)      String array containing labels for the Plots.
28610      !      Ximage$(*)      String array containing image formats for the X axis labeling.
28620      !      Yimage$(*)      String array containing image formats for the Y axis labeling.
28630      !      Legend$(*)      String array containing labels for each symbol in a profile plot.
28640      !      Symbols(*)      Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
28650      !      G              Used as an index to the above arrays. Specifies one of nine plots.
28660      ! OPTION BASE 1
28670      ! COM /Graph/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),
      !      Legend$(*)
28680      ! DIM L$(80)
28690      ! ON ERROR CALL Error
28700      ! PLOTTER IS CRT,"INTERNAL"
28710      ! Define the pen numbers for the colors black and white.
28720      ! Black=-1
28730      ! White=1

```

```

28740 CSIZE 100*15/1023 ! Select a character labeling size of 15 pixels high.
28750 ! Define the values for the left,right,bottom,top ends of the horizontal and vertical scales.
28760 Xmin=Wndw(G,1)
28770 Xmax=Wndw(G,2)
28780 Ymin=Wndw(G,3)
28790 Ymax=Wndw(G,4)
28800 ! Define the values for the left,right,bottom,top pixel locations for the plot.
28810 Xpix1=Vwprt(G,1)
28820 Xpix2=Vwprt(G,2)
28830 Ypix1=Vwprt(G,3)
28840 Ypix2=Vwprt(G,4)
28850 ! Define the step size between grid lines, axis tick marks, and axis labels.
28860 Xstep=(Xmax-Xmin)/Xdiv(G)
28870 Ystep=(Ymax-Ymin)/Ydiv(G)
28880 ! Define the amount of scale X and Y which equals the size of one pixel (picture element).
28890 Xpixel=(Xmax-Xmin)/(Xpix2-Xpix1)
28900 Ypixel=(Ymax-Ymin)/(Ypix2-Ypix1)
28910 ! Clear the plots back ground & plot area and also draw the plots borders, grids, and axes.
28920 AREA PEN Black
28930 PEN White
28940 GOSUB Back_ground
28950 GOSUB Axes
28960 !GOSUB Grid
28970 GOSUB Plot_area
28980 ! Draw the X and Y axis labels.
28990 CLIP OFF
29000 GOSUB Ylabel
29010 GOSUB Xlabel
29020 ! Create a legend to define which symbol is used with which data.
29030 CALL Legend(G,Symbols(*))
29040 OFF ERROR
29050 SUBEXIT
29060 Back_ground: ! This subroutine clears the plot's background.
29070 VIEWPORT (Xpix1-75)/10.23,(Xpix2+25)/10.23,(Ypix1-33)/10.23,(Ypix2+6)/10.23
29080 WINDOW -1.E+9,1.E+9,-1.E+9,1.E+9
29090 MOVE 0,0
29100 WINDOW 0,1,0,1
29110 MOVE 0,0
29120 RECTANGLE 1,1,FILL
29130 RETURN
29140 Axes: ! This subroutine draws the plot's X and Y axes.
29150 VIEWPORT (Xpix1-1)/10.23,(Xpix2+1)/10.23,(Ypix1-6)/10.23,(Ypix1-1)/10.23
29160 WINDOW Xmin,Xmax,1,0
29170 AXES Xstep,2,Xmin,0,1,1,1
29180 VIEWPORT (Xpix1-1)/10.23,(Xpix2+1)/10.23,(Ypix2+1)/10.23,(Ypix2+6)/10.23
29190 WINDOW Xmin,Xmax,0,1
29200 AXES Xstep,2,Xmin,0,1,1,1
29210 VIEWPORT (Xpix1-6)/10.23,(Xpix1-1)/10.23,(Ypix1-1)/10.23,(Ypix2+1)/10.23
29220 WINDOW 1,0,Ymin,Ymax
29230 AXES 2,Ystep,0,Ymin,1,1,1
29240 VIEWPORT (Xpix2+1)/10.23,(Xpix2+6)/10.23,(Ypix1-1)/10.23,(Ypix2+1)/10.23
29250 WINDOW 0,1,Ymin,Ymax
29260 AXES 2,Ystep,0,Ymin,1,1,1
29270 RETURN
29280 Grid: ! This subroutine draws the plot's X and Y grid lines.
29290 VIEWPORT (Xpix1-1)/10.23,(Xpix2+1)/10.23,(Ypix1-1)/10.23,(Ypix2+1)/10.23
29300 WINDOW Xmin,Xmax,Ymin,Ymax
29310 LINE TYPE 4
29320 GRID Xstep,Ystep,Xmin,Ymin
29330 LINE TYPE 1
29340 RETURN
29350 Plot_area: ! This subroutine selects part of the CRT plot area and give it scales for the X and Y axes.
29360 VIEWPORT Xpix1/10.23,Xpix2/10.23,Ypix1/10.23,Ypix2/10.23
29370 WINDOW Xmin,Xmax,Ymin,Ymax
29380 RETURN
29390 Xlabel: ! This subroutine labels the X axis and also names the X axis.
29400 LOG 5
29410 FOR X=Xmin TO Xmax+Xstep/100 STEP Xstep
29420 MOVE X,Ymin-12*Ypixel
29430 OUTPUT L$ USING Ximage$(G);X
29440 LABEL TRIM$(L$)
29450 NEXT X
29460 MOVE (Xmin+Xmax)/2,Ymin-25*Ypixel
29470 LABEL Xlabel$(G)
29480 RETURN
29490 Ylabel: ! This subroutine labels the Y axis and also names the Y axis.
29500 LOG 8
29510 Len=0
29520 FOR Y=Ymin TO Ymax+Ystep/100 STEP Ystep
29530 MOVE Xmin-5*Xpixel,Y

```



```

29540         OUTPUT L$ USING Yimage$(G);Y
29550         LABEL TRIMS(L$)
29560         Len=MAX(Len,LEN(TRIMS(L$)))
29570     NEXT Y
29580     MOVE Xmin-(5+7*Len)*Xpixel,(Ymin+Ymax)/2
29590     LABEL Ylabel$(G)
29600     LOG 5
29610     RETURN
29620 SUBEND
29630 Legend: SUB Legend(G,Symbols(*) )
29640     ! Description:
29650     !     This subprogram produces a legend within one of the nine plots on the CRT screen.
29660     ! Variables:
29670     !     Wndw(*)      Array containing the plot's scales.
29680     !     Vwprt(*)     Array containing the plot's CRT position.
29690     !     Xdiv(*)     Array containing the number of X divisions for the plot's X axis.
29700     !     Ydiv(*)     Array containing the number of Y divisions for the plot's Y axis.
29710     !     Xlabel$(*)  String array containing labels for the X axis.
29720     !     Ylabel$(*)  String array containing labels for the Y axis.
29730     !     Title$(*)   String array containing labels for the Plots.
29740     !     Ximage$(*)  String array containing image formats for the X axis labeling.
29750     !     Yimage$(*)  String array containing image formats for the Y axis labeling.
29760     !     Legend$(*)  String array containing labels for each symbol in a profile plot.
29770     !     Symbols(*)  Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
29780     !     Symbol(*)   Array of coordinates which when connected produce a distinct geometric symbol.
29790     !     G           Used as an index to the above arrays. Specifies one of nine plots.
29800     !     S           Used to index the Legend$ array.
29810     !     Noc         The number of coordinates in a symbol.
29820     !     Len         Length of a Legend$ array element.
29830     OPTION BASE 1
29840     COM /Graph/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),
        Legend$(*)
29850     DIM Symbol(20,3)
29860     VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
29870     WINDOW Vwprt(G,1),Vwprt(G,2),Vwprt(G,3),Vwprt(G,4)
29880     ! Define the pen numbers for the colors black and white.
29890     Black=-1
29900     White=1
29910     ! Define the colors for symbol filling and edge drawing.
29920     AREA PEN Black
29930     PEN White
29940     CSIZE 100*15/1023      ! Select a character labeling size of 15 pixels high.
29950     LOG 2
29960     ! Calculate the maximum length of all of the symbol labels.
29970     Len=0
29980     FOR S=1 TO SIZE(Legend$,2)
29990         Len=MAX(LEN(Legend$(G,S)),Len)
30000     NEXT S
30010     ! For each symbol put up a sample symbol and its label.
30020     FOR S=1 TO SIZE(Legend$,2)
30030         IF LEN(Legend$(G,S))=0 THEN 30110
30040         Noc=Symbols(S,0,1)
30050         REDIM Symbol(Noc,3)
30060         MAT Symbol= Symbols(S,1:Noc,*)
30070         MOVE Vwprt(G,2)-7*Len-23,Vwprt(G,4)-15*S+5
30080         SYMBOL Symbol(*),FILL,EDGE
30090         MOVE Vwprt(G,2)-7*Len-10,Vwprt(G,4)-15*S+4
30100         LABEL Legend$(G,S)
30110     NEXT S
30120 SUBEND
30130 Histo:
30140 Rt_histo: SUB Rt_histo(@Lvdas,Symbols(*),Repeat)
30150     ! Description:
30160     !     This subprogram plots real time histograms within five of the nine plots on the CRT screen.
30170     !     The histogram data is acquired from the LVDAS over a specified acquisition time.
30180     ! Variables Defined in Main Program:
30190     !     Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
30200     ! Local Variables:
30210     !     Histo(*)     Array of bin numbers, old histogram bin heights, and new histogram bin heights.
30220     !     Nbins        Number of bins in the Histo(*).
30230     !     Bin          2^Bin is the bin width of individual histogram vertical bars.
30240     !     Min          Minimum value for histogram. Left side of histogram scale.
30250     !     Max          Maximum value for histogram. right side of histogram scale.
30260     !     F1           Upper 16bits of integerized Min.
30270     !     F2           Lower 16bits of integerized Min.
30280     !     A1           Upper 16bits of integerized histogram acquisition time.
30290     !     A2           Lower 16bits of integerized histogram acquisition time.
30300     !     Nnew         Number of samples in the most up to date histogram.
30310     !     Nold         Number of samples in the previous histogram.
30320     !     N(*)         Number of samples for each histogram of the five separate channels.

```

```

30330      ! Channel Used to select the LVDAS channel that will be sampled for a histogram.
30340      ! Kw      Converts Hz to MHz or raw data to volts.
30350      ! Ww      Window width of each vertical histogram bar.
30360      ! Old     Histogram height of previous histogram at a particular bin.
30370      ! New     Histogram height of current histogram at a particular bin.
30380      ! X1      Horizontal position of histogram rectangle.
30390      ! Y1      Vertical position of histogram rectangle.
30400      ! X2      Horizontal width of histogram rectangle.
30410      ! Y2      Vertical width of histogram rectangle.
30420      ! I       Used as an index to the Histo(*). Specifies one of Nbins bins.
30430      ! G       Used as an index to the graphics arrays. Specifies one of nine plots.
30440      OPTION BASE 1
30450      COM /Graph/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),
          Legend$(*)
30460      INTEGER Histo(1000,3),Nbins,F1,F2,A1,A2
30470      REAL Nnew,Nold,N(5)
30480      ! Clear all of the histogram data within the LVDAS.
30490      OUTPUT @Lvdas USING "AA","CA"
30500      ! Draw new plots for the five histograms.
30510      FOR Channel=1 TO 5
30520          CALL Set_up(Channel,Symbols(*))
30530      NEXT Channel
30540      ! Calculate the acquisition time. 0.1*10000000 will give an acquisition of 0.1 seconds.
30550      CALL Convert2words(.1*10000000,A1,A2)
30560      ! Enable the keyboard to terminate histogram plotting.
30570      ON KBD GOSUB Hdone
30580      REPEAT
30590          FOR Channel=1 TO 5
30600              G=Channel
30610              SELECT Channel
30620                  CASE 1,2,3
30630                      Kw=1000000
30640                      Min=Kw*Wndw(G,1)
30650                      Max=Kw*Wndw(G,2)
30660                      Bin=INT(LGT((Max-Min)/100)/LGT(2))+1
30670                      Ww=2^Bin
30680                      CALL Convert2words(Min,F1,F2)
30690                  CASE 4,5
30700                      Kw=32768/5
30710                      Min=Kw*Wndw(G,1)
30720                      Max=Kw*Wndw(G,2)
30730                      Bin=INT(LGT((Max-Min)/100)/LGT(2))+1
30740                      Ww=2^Bin
30750                      CALL Convert2words(Min,F1,F2)
30760                  CASE ELSE
30770                      GOTO 31100
30780              END SELECT
30790      Hsend:      ! Tell the LVDAS to Take a Histogram.
30800      OUTPUT @Lvdas USING "AA,6(W)","TH",F1,F2,Bin,A1,A2,Channel
30810      Henter:    ! Enter number of bins in the histogram.
30820      ENTER @Lvdas USING "#,W",Nbins
30830      ! Redimension the Histo(*) and then enter the histogram data.
30840      IF Nbins>0 THEN
30850          REDIM Histo(Nbins,3)
30860          ENTER @Lvdas USING "#,W",Histo(*)
30870      END IF
30880      ! Enter the number of samples for the previous and current histogram.
30890      ENTER @Lvdas USING "#,W",Nnew,Nold
30900      Hplot:      ! Scale part of the CRT for the histogram plotting.
30910      VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
30920      WINDOW Kw*Wndw(G,1),Kw*Wndw(G,2),Wndw(G,3),Wndw(G,4)
30930      Xpixel=Kw*(Wndw(Channel,2)-Wndw(Channel,1))/(Vwprt(Channel,2)-Vwprt(Channel,1))
30940      N1=N(Channel)
30950      N2=N(Channel)-Nold+Nnew
30960      N(Channel)=N(Channel)-Nold+Nnew
30970      FOR I=1 TO Nbins
30980          Old=MIN(Histo(I,3),Wndw(Channel,4))
30990          New=MIN(Histo(I,2),Wndw(Channel,4))
31000          AREA PEN SGN(New-Old)
31010          X1=Histo(I,1)*Ww+Min
31020          X2=Ww
31030          Y1=Old
31040          Y2=New-Old
31050          IF X1<Kw*Wndw(G,1) THEN X1=Kw*Wndw(G,1)
31060          IF X1>Kw*Wndw(G,2)-X2 THEN X1=Kw*Wndw(G,2)-X2
31070          MOVE X1,Y1
31080          RECTANGLE X2-Xpixel,Y2,FILL
31090      NEXT I
31100      NEXT Channel
31110      UNTIL KBD$<>" OR NOT Repeat

```

```

31120 SUBEXIT
31130 Hdone: Done=1
31140 RETURN
31150 SUBEND
31160 Histo: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
31170 Pt_histo: SUB Pt_histo(Symbols(*),Run,File,Pos,INTEGER Nsam)
31180 ! Description:
31190 ! This subprogram plots post time histograms within five of the nine plots on the CRT screen.
31200 ! The histogram data is acquired from the LVDAS over a specified acquisition time.
31210 ! Variables Defined in Main Program:
31220 ! Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
31230 ! U1(*),V1(*),W1(*),A1(*),B1(*)
31240 ! Local Variables:
31250 ! Histo(*) Array of histogram bin heights indexed by bin number.
31260 ! Data(*) Array of instantaneous U,V,W velocity or A,B voltage data.
31270 ! Nsam Number of samples acquired.
31280 ! Xmin Minimum value for histogram. Left side of histogram scale.
31290 ! Xmax Maximum value for histogram. right side of histogram scale.
31300 ! Xwin Window width of each vertical histogram bar.
31310 ! K Used as an index to the above arrays.
31320 ! L Used as an index to the Histo(*). Specifies one of 100 bins.
31330 ! Xpixel Horizontal length of one picture on the CRT in scale units.
31340 ! Channel Selects one of the 5 channels of U1(*),V1(*),W1(*),A1(*),B1(*) data.
31350 ! G Used as an index to the graphics arrays. Specifies one of nine plots.
31360 OPTION BASE 1
31370 COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
31380 COM /Graph/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),
Legend$(*)
31390 INTEGER Histo(0:100)
31400 REAL Data(1000)
31410 REDIM Data(Nsam)
31420 FOR Channel=1 TO 5
31430 ! Fill the data array with U1(*),V1(*),W1(*),A1(*), or B1(*) depending on Channel.
31440 G=Channel
31450 IF Channel=1 THEN MAT Data= U1
31460 IF Channel=2 THEN MAT Data= V1
31470 IF Channel=3 THEN MAT Data= W1
31480 IF Channel=4 THEN MAT Data= A1
31490 IF Channel=5 THEN MAT Data= B1
31500 ! Draw a new empty histogram plot.
31510 CALL Set_up(Channel,Symbols(*))
31520 Hsort: Xmin=Wndw(Channel,1)
31530 Xmax=Wndw(Channel,2)
31540 Xwin=(Xmax-Xmin)/100
31550 ! Sort the data into a histogram.
31560 MAT Data= Data-(Xmin)
31570 MAT Data= Data/((Xmax-Xmin)/100)
31580 MAT Histo= (0)
31590 FOR K=1 TO Nsam
31600 L=MAX(MIN(Data(K),100),0)
31610 Histo(L)=Histo(L)+1
31620 NEXT K
31630 Hplot: ! Scale part of the CRT for histogram plotting.
31640 VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
31650 WINDOW 0,100,Wndw(G,3),Wndw(G,4)
31660 Xpixel=(100-0)/(Vwprt(Channel,2)-Vwprt(Channel,1))
31670 ! Draw the histogram.
31680 FOR K=0 TO 100
31690 IF Histo(K) THEN
31700 MOVE K-.5,0
31710 AREA PEN SGN(1) ! Positive pens will plot while negative histograms erase.
31720 RECTANGLE 1-Xpixel,Histo(K),FILL
31730 END IF
31740 NEXT K
31750 NEXT Channel
31760 SUBEXIT
31770 SUBEND
DONE

```

REPORT DOCUMENTATION PAGE			Form Approved OMB No. 0704-0188	
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503.				
1. AGENCY USE ONLY (Leave blank)		2. REPORT DATE January 1992		3. REPORT TYPE AND DATES COVERED Final Contractor Report
4. TITLE AND SUBTITLE Measurement of Vortex Flow Fields			5. FUNDING NUMBERS C NAS1-18667 WU 505-62-30-01	
6. AUTHOR(S) T. Kevin McDevitt, Todd A. Ambur, Gary M. Orngard, F. Kevin Owen				
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Complere, Inc. P.O. Box 1697 Palo Alto, CA 94302			8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) National Aeronautics and Space Administration Langley Research Center Hampton, VA 23665-5225			10. SPONSORING / MONITORING AGENCY REPORT NUMBER NASA CR-189543	
11. SUPPLEMENTARY NOTES Technical Monitor - Bobby L. Berrier, MS 280, Langley Research Center, Hampton, VA 23665 (804) 864-3001 FTS 928-3001				
12a. DISTRIBUTION / AVAILABILITY STATEMENT Unclassified - Unlimited Subject Category 02			12b. DISTRIBUTION CODE	
13. ABSTRACT (Maximum 200 words) <p>A three-dimensional laser fluorescence anemometer has been designed, built, and demonstrated for use in the Langley 16-by 24-Inch Water Tunnel. Innovative optical design flexibility combined with compact and portable data acquisition and control systems have been incorporated into the instrument. This will allow its use by NASA in other test facilities. A versatile fiber optic system facilitates normal and off-axis laser beam alignment, removes mirror losses and improves laser safety. This added optical flexibility will also enable simple adaptation for use in the adjacent jet facility. New proprietary concepts in transmitting color separation, light collection, and novel prism separation of the scattered light have also been designed and built into the system. Off-axis beam traverse and alignment complexity led to the requirement for a specialized, programmable traverse controller and the inclusion of an additional traverse for the off-axis arm. To meet this challenge, an "in-house" prototype unit was designed and built and traverse control software developed specifically for the water tunnel traverse applications. A specialized data acquisition interface was also required. This was designed and built for the Laser Fluorescence Anemometer system.</p>				
14. SUBJECT TERMS Laser Velocimetry Water Tunnel Measurements			15. NUMBER OF PAGES 91	
			16. PRICE CODE A05	
17. SECURITY CLASSIFICATION OF REPORT Unclassified	18. SECURITY CLASSIFICATION OF THIS PAGE Unclassified	19. SECURITY CLASSIFICATION OF ABSTRACT Unclassified	20. LIMITATION OF ABSTRACT UL	